

PRODUCT CATALOG



AMSINO
Advancing Healthcare Worldwide®



AMSure® Enteral Feeding Systems..... 1

PUGGLE® Enteral Feeding Pump.....	2
Sentinel™ Enteral Feeding Pump.....	3
ENFit® and Non-ENFit Sets.....	3
ENFit and Non-ENFit Syringes.....	3
ENFit and Non-ENFit Irrigation Kits and Trays.....	5

AMSafe® Medication Delivery Systems 6

Gravity Administration Sets.....	7
10 Drops/mL.....	7
15 Drops/mL.....	9
60 Drops/mL.....	13
AMSafe-3® Selectable Drip Chamber	
10, 15 & 60 drops/mL.....	14
Secondary Piggyback Sets.....	15
Y-type Blood Sets.....	15
IV Extension Sets.....	16
Standard Extension Sets.....	16
Microbore Extension Sets.....	17
Smallbore Extension Sets.....	17
Flow Controller Extension Sets.....	17
Pressure Infusion Extension Sets.....	18
Needle-Free Connectors.....	18
Bag Spike.....	18
Needle-Free Vial Adapter.....	18
IV Pump Administration Sets.....	19
Sure-Lok® Needle-Free IV Pump Administration Sets.....	19
Pre-Pierced Y-Injection IV Pump Administration Sets.....	20

AMSure Urological Care 21

Foley Catheter Trays & Collection.....	22
Pre-Connected Foley Catheter Urine Drainage Bag Trays.....	22
Add-a-Foley Catheter Trays.....	23
Pre-Connected Foley Catheter Urine Meter Kits.....	24
Add-a-Foley Urine Meter Trays.....	25
Foley Insertion Trays.....	25
Urinary Drainage Bags.....	26
Urine Meter.....	27
ACCU-MEASURE® Pediatric Urine Meter.....	27
Foley Catheters.....	28
Hydrophilic Latex Foley Catheters.....	28
Silicone-Coated Latex Foley Catheters.....	29
100% Silicone Foley Catheters.....	30
Intermittent Urinary Catheters.....	31
Urethral Insertion Trays and Kits.....	31
Closed System Urethral Drainage Kits.....	32
Self-Catheterization Kits with R-Polished™ Eye Catheters.....	33
PVC Intermittent Urethral Catheter with R-Polished Eyes.....	34
PVC Intermittent Urethral Catheter.....	35
Red Rubber Latex Urethral Catheter.....	35
Irrigation Sets.....	36
Cystoscopy & TUR Irrigation Sets.....	36
Ancillary & Patient Care Products.....	37
Urinary Leg Bags.....	37
Extension Tube for Leg Bags.....	37

Pediatric Urine Collector.....	37
Irrigation Trays and Kits.....	38
Syringes.....	39
Specimen Containers.....	40

Sterile Water & Saline Products..... 41

Sterile Water for Inhalation, USP.....	42
Sterile Water Prefilled Humidifiers, USP.....	42
Nebulizer Kits.....	69
Saline Products.....	43
PreFilled Saline Flush Syringes.....	43
Sodium Chloride Inhalation Solution, (0.9%) USP.....	44
Device and Wound Irrigation.....	44

Suction & Waste Fluid Management Systems..... 45

EZE-VAC® Suction Canister System.....	46
EZE-VAC 2.0 Suction Canister System.....	47
RECEPTAL Liners.....	48
RECEPTAL Twist Liners.....	49
RECEPTAL Outer Shells.....	50
ACCU-MEASURE® Critical Measurement Devices.....	51
RECEPTAL Solidifiers.....	51
AMSure® Suction Tubing with Connectors.....	52
Suction Yankauers and Kits.....	53
SteriCadi™ Suction Device Holster System.....	54
Suction Catheters and Trays.....	55
Suction Catheter Kits.....	56
Accessories.....	57
Reusable Hardware.....	60
Regulator Mounts.....	62
iRECEPTAL® Surgical Waste Fluid Management System.....	63

AMSure Respiratory Therapy Systems..... 65

Respiratory Masks and Nasal Cannulas.....	66
Nebulizer Kits.....	69
Accessories.....	70
Aerosol Drainage Bag.....	70
Heat and Moisture Exchangers.....	71
Tracheostomy Clean and Care Trays.....	71
Sterile Water.....	72
Suction Catheters, Trays and Kits.....	73
Suction Tubing with Connectors.....	75
Suction Yankauers and Kits.....	76

AMSure Surgical & Nursing Systems..... 77

Suction Tubing with Connectors.....	78
Suction Yankauer and Kits.....	79
Suction Device Holster and Kits.....	79
Specimen Containers.....	80
Vaginal Speculums.....	80
Disposable Umbilical Cord Clamp.....	80
Ear/Ulcer Syringes.....	81
Cleansing Enema Sets.....	81
Connectors.....	82

Index of REFS 83

AMSure®

ENTERAL FEEDING SYSTEMS



Optimizing Nutritional Delivery

AMSure Enteral Feeding Systems enhance delivery from nutrition source to patient feeding tube. The comprehensive enteral product portfolio is available with standard and ENFit® connections. The systems are designed to provide best patient outcomes and includes:

PUGGLE® Enteral Feeding Pump	2
Sentinel™ Enteral Feeding Pump.....	3
ENFit & Non-ENFit Enteral Feeding Sets.....	3
ENFit & Non-ENFit Enteral Syringes	4
ENFit & Non-ENFit Irrigation Kits & Trays	5

PUGGLE® Enteral Feeding Pump

The versatile pump is designed for pediatric and adult ambulatory applications. Simplified program settings make feeding on-the-go a snap.

Puggle is only available the US.

- Large color LED screen for easy at-a-glance view of clock and settings
- Compact lightweight design fits into backpacks and transports easily
- Quiet motor with continuous peristaltic rotary pump action
- 20 hour battery life that fully charges in only 2.5 hours



Dimensions: 5.4 x 5.0 x 2.5 inches (13.7 x 12.7 x 6.4 cm)
Weight: 1.43 lbs (649 g)
Operating Range: 1 to 400 mL/h in 1 mL increments
Delivered Volume: 1 – 9999 mL in 1 mL increments
Accuracy: +/- 5%



REF	DESCRIPTION	QTY/UNIT
PG001	PUGGLE Enteral Feeding Pump with Pole Clamp	1/Case
PE90ST	PUGGLE Enteral Pump Spike Set, Anti-free Flow Valve, ENFit® and Transition Connector	30/Case
PE90B5T	PUGGLE Enteral Pump Set with 500 mL Bag, Anti-free Flow Valve, ENFit and Transition Connector	30/Case
PE90B12T	PUGGLE Enteral Pump Set with 1200 mL Bag, Anti-free Flow Valve, ENFit and Transition Connector	30/Case

Sentinel™ Enteral Feeding Pump



The Sentinel enteral feeding pump is easy to use and program which minimizes training time and helps to improve programming efficacy. The AMSure® enteral feeding pump sets are recommended to be used with the Sentinel pump.

- Large display and buttons promote ease of programming
- Intuitive design simplifies tubing set-up
- Unique drop sensor minimizes nuisance alarms
- Tubing securements help support uninterrupted feeding
- Polycarbonate pump housing increases durability
- Designed for feeding or fluid applications

REF	DESCRIPTION	QTY/UNIT
S-1000-SI	Sentinel Enteral Pump	1/Each

ENFit® and Non-ENFit Enteral Feeding Sets



ENF1200K

AMSure Enteral Feeding Sets are available with and without ENFit connectors. The feeding sets are intended for administration of enteral nutrition, fluids or medication. AMSure enteral feeding pump sets are recommended to use with the Sentinel™ enteral pump and are compatible with the Compat™ and EntraFlo™ pumps.

- ENFit connector sets are compliant with the ISO 80369-3 to help reduce risk of misconnections
- Patient label provides convenient patient documentation
- ENFit sets include pre-attached stepped transition connector
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	ENFit CONNECTOR	FEEDING BAG SIZE			FEEDING SET TYPE		PUMP COMPATIBILITY		QTY/UNIT
		Spike Only- No Bag	500 mL	1200 mL	Gravity Set	Pump Set	Sentinel®	Compat™/ EntraFlo™	
E-1004		•				•	•	•	30/Case
E-1200				•		•	•	•	30/Case
E-1203				•	•				30/Case
ENF1004K ©	•	•				•	•	•	30/Case
ENF1200K ©	•			•		•	•	•	30/Case
ENF1203K ©	•			•	•				30/Case
ENF500K ©	•		•			•	•	•	30/Case

EntraFlo™ is a registered trademark of Medline Industries, Inc.
 Compat™ is a registered trademark of Société des Produits Nestlé, S.A Corporation.
 Sentinel™ is a registered trademark of Alcor Scientific Inc. ENFit® is a registered trademark of GEDSA.

ENFit® and Non-ENFit Enteral Syringes

AMSure® Enteral Syringes are available with or without a transition adapter ENFit tip. The ENFit adapter syringe tip is compliant with the ENFit ISO 80369-3 design standard which supports the goal to improve patient safety by reducing the risk of misconnections. All syringes include a tip protector and are DEHP-free, BPA-free and not manufactured with natural rubber latex.

Bulb Syringe

- Flexible soft bulb promotes one-handed administration



Thumb Controlled Piston Syringe

- Designed for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage



Flat Top Piston Syringe

- Visible dosing marks designed for ease of use
- Precision crafted O-ring plunger ensures a positive seal
- Visible dosing marks designed for ease of use



REF	ENFIT TIP	TIP ADAPTER	60 ML SYRINGE TYPE			PACKAGE TYPE			BULK	STERILE	QTY/UNIT
			Bulb	Thumb Control Ring	Flat Top	Poly Pouch	Form Fill Seal	Resealable Pole Bag			
AS010			•						•		50/Case
AS011 ©			•				•			•	50/Case
AS115 ©		•			•		•			•	50/Case
AS116 ©		•			•						30/Case
ENS015 ©	•			•			•			•	50/Case
ENS016	•			•							30/Case
ENS115 ©	•				•		•			•	50/Case
ENS116	•				•						30/Case
ENS115NS	•				•						30/Case
ENS115T*	•				•						30/Case

*Includes Transition Adapter

ENFit® and Non-ENFit Enteral Feeding Irrigation Kits and Trays

AMSure® Enteral Feeding Irrigation Kits and Trays are designed for simplicity and efficiency of care.

Kits are non-sterile and include a patient label, 500 mL graduated container, either a poly pouch or re-sealable bag and a 60 mL syringe with or without an ENFit tip.

Trays are sterile and include 500 mL graduated container, 1000 mL tray, alcohol prep pad and moisture-proof under pad, a 60 mL syringe with or without an ENFit tip.

- ENFit® connector sets are compliant with ISO 80369-3 to help reduce risk of misconnections
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector cap



Enteral Feeding Irrigation Kits

REF	ENFIT® SYRINGE TIP	500 mL GRADUATED CONTAINER	60 mL BULB SYRINGE	60 mL THUMB CONTROL PISTON SYRINGE	POLY POUCH	RE-SEALABLE POLE BAG	PATIENT LABEL	QTY/UNIT
AS121		•	•				•	30/Case
AS126		•		•			•	30/Case
AS127		•		•		•	•	30/Case
ENK127 ©	•	•		•		•	•	30/Case

Enteral Feeding Irrigation Trays

REF	ENFIT® SYRINGE TIP	1000 mL OUTER TRAY	500 mL GRADUATED CONTAINER	60 mL BULB SYRINGE	60 mL THUMB CONTROL PISTON SYRINGE	STERILE	ALCOHOL PAD	LARGE UNDERPAD	QTY/UNIT
AS130 ©		•	•	•		•	•	•	20/Case
AS136 ©		•	•		•	•	•	•	20/Case
ENK136 ©	•	•	•		•	•	•	•	20/Case

AMSafe®

MEDICATION DELIVERY SYSTEMS



Safety First

The AMSafe IV product portfolio provides solutions to help improve patient and caregiver safety and infection prevention. Our IV product selection includes:

Gravity Administration IV Sets.....	7
10 Drops/mL.....	7
15 Drops/mL.....	9
60 Drops/mL.....	13
AMSafe-3® Selectable Drip Chamber, 10, 15 & 60 drops/mL.....	14
Secondary Piggyback Sets	15
Y-type Blood Sets	15
IV Extension Sets	16
Standard Extension Sets.....	16
Microbore Extension Sets.....	17
Smallbore Extension Sets.....	17
Flow Controller Extension Sets.....	17
Pressure Infusion Extension Sets	18
Needle-Free Connectors	18
Bag Spike.....	18
Needle-Free Vial Adapter	18
IV Pump Administration Sets	19
Sure-Lok® Needle-Free IV Pump Administration Sets.....	19
Pre-Pierced Y-Injection IV Pump Administration Sets.....	20







Gravity Administration IV Sets

Amsino's comprehensive line of gravity administration IV sets are designed to improve safety and infection prevention for both caregivers and patients.

- Available in needle-free, pre-pierced and injection Y-sites in a variety of configurations to meet your clinical needs
- Available with or without pre-attached extension sets
- 10, 15 or 60 drops/mL sets available
- Includes a rotating male luer lock
- Sterile, non-pyrogenic, DEHP-free, not manufactured with natural rubber latex

10 Drops/mL

General Primary Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
107201	10	72 in. (183 cm)	Injection Y-site	14 mL	50/Case
108301 ©	10	83 in. (211 cm)	Injection Y-site	15.6 mL	50/Case
					
108303	10	83 in. (211 cm)	Sure-Lok® Needle-free Y-site	16 mL	50/Case
					
108304	10	89 in. (226 cm)	Non-vented, Roller Clamp (2), Sure-Lok Needle-free Y-site with Pre-attached Ext Set with Pinch Clamp	17.5 mL	50/Case
					
109602 ©	10	96 in. (244 cm)	(2) Injection Y-sites	17.5 mL	50/Case
					
108306 ©	10	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site	16 mL	50/Case
					
MRSA10-84	10	84 in. (213 cm)	Vented/Non-vented, Pre-pierced Y-site, MicroClave® Clear Needle-free Y-site	10.5 mL	50/Case
					
AC3110 ©	10	105 in. (266 cm)	(2) Sure-Lok Needle-free Y-sites	20 mL	50/Case
					
MRA10E-95	10	95 in. (241 cm)	Vented/Non-vented, Pre-pierced Y-site, Slide Clamp, Pre-attached Ext Set, Slide Clamp, MicroClave Clear Needle-free Y-site	18.5 mL	50/Case
					

10 Drops/mL

General Primary Sets (continued)


REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
108305 ©	10	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok® Needle-free Y-site with Pre-attached Ext Set, Pinch Clamp	17.5 mL	50/Case




MRA10E-104	10	104 in. (264 cm)	Vented/Non-vented, MicroClave® Clear Needle-free Y-sites, Pre-pierced Y-site, with Pre-attached Ext Set, Slide Clamp, Pre-pierced Y-site, MicroClave Clear Needle-free Y-site	13.4 mL	50/Case
-------------------	----	------------------	---	---------	---------



MRSA10E-107	10	107 in. (272 cm)	Vented/Non-vented, (2) MicroClave Clear Needle-free Y-sites, Pre-attached Ext Set	13.4 mL	50/Case
--------------------	----	------------------	---	---------	---------




108308 ©	10	83 in. (211 cm)	(2) Sure-Lok Needle-free Y-sites, Pre-pierced Y-site	16.6 mL	50/Case
-----------------	----	-----------------	--	---------	---------




10 Drops/mL

Flow Controller Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AFS102 ©	10	100 in. (254 cm)	Vented/Non-vented, Flow Controller, Slide Clamp, Sure-Lok Needle-free Y-site	18.5 mL	50/Case



AFS208 ©	10	100 in. (254 cm)	Vented/Non-vented, Flow Controller, Slide Clamp, (2) AMSafe® Needle-free Y-sites	19.2 mL	50/Case
-----------------	----	------------------	--	---------	---------



AFS208B	10	96 in. (244 cm)	Vented/Non-vented, Back Check Valve, Flow Controller, (2) AMSafe Needle-free Y-sites, Slide Clamp	19 mL	50/Case
----------------	----	-----------------	---	-------	---------








Gravity Administration IV Sets

10 Drops/mL



General Primary Sets with Back Check Valve

- Provides the flexibility to infuse more than one bag of medication or solution
- Featuring a back check valve to prevent fluid exchange or backflow into the primary infusion line

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
1010001B ©	10	100 in. (254 cm)	Vented/Non-vented, Back Check Valve, Sure-Lok® Needle-free Y-site	19.3 mL	50/Case
					
108302 ©	10	83 in. (211 cm)	Vented/Non-vented, Back Check Valve, Pre-pierced Y-site, Sure-Lok Needle-free Y-site, Rotating Male Luer	16.4 mL	50/Case
					
MRA102-84 ©	10	84 in. (213 cm)	Vented/Non-vented, Back Check Valve, Pre-pierced Y-site, (2) MicroClave® Clear Needle-free Y-sites	10.7 mL	50/Case
					
MRSA103-94	10	94 in. (239 cm)	Vented/Non-vented, Back Check Valve, Pre-pierced Y-site, (3) MicroClave Clear Needle-free Y-sites	17.5 mL	50/Case
					
MRSA10E-110	10	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (2) Pre-pierced Y-sites, (2) MicroClave Clear with Pre-attached Ext Set	21.1 mL	50/Case
					

15 Drops/mL

General Primary Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AA3101 ©	15	78 in. (198 cm)	Vented/Non-vented, Slide Clamp, Injection Y-site	15 mL	50/Case
					
158310 ©	15	83 in. (211 cm)	Vented/Non-vented, Pre-pierced Y-site	17 mL	50/Case
					

15 Drops/mL General Primary Sets








REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
158309 ©	15	83 in. (211 cm)	Vented/Non-vented, Sure-Lok® Needle-free Y-site	17 mL	50/Case
					
MRSA15-100	15	100 in. (254 cm)	Vented/Non-vented, MicroClave® Clear Needle-free Y-site	12.4 mL	50/Case
					
158306 ©	15	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site	16 mL	50/Case
					
MRSA15-91	15	91 in. (231 cm)	Vented/Non-vented, Pre-pierced Y-site, MicroClave Clear Needle-free Y-site	11.4 mL	50/Case
					
159602 ©	15	96 in. (244 cm)	(2) Injection Y-sites	17.7 mL	50/Case
					
MRA15-100	15	100 in. (254 cm)	Vented/Non-vented, Pre-pierced Y-site, MicroClave Clear Needle-free Y-site	13 mL	50/Case
					
1510501 ©	15	105 in. (266 cm)	Vented/Non-vented, (2) Sure-Lok Needle-free Y-sites	20 mL	50/Case
					
158305 ©	15	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site with Preattached Ext Set, Pinch Clamp	17.5 mL	50/Case
					
MRSA15E-88	15	88 in. (224 cm)	Vented/Non-vented, Pinch Clamp, Pre-pierced Y-site with Pre-attached Ext Set, Pinch Clamp, (2) MicroClave Clear Needle-free Y-sites	11.1 mL	50/Case
					
158307 ©	15	89 in. (226 cm)	Sure-Lok Needle-free Y-site, Pre-pierced Y-site with Preattached Ext Set, Sure-Lok Needle-free Y-site, Slide Clamp, Pre-pierced Y-site	18.7 mL	50/Case
					
MRSA15E-120	15	120 in. (305 cm)	Vented/Non-vented, Pre-pierced Y-site, MicroClave Clear Needle-free Y-site, Rotating Male Luer Lock, Drop-In Ext Set MicroClave Clear Y-site	15.8 mL	50/Case
					

Gravity Administration IV Sets

15 Drops/mL







Piggyback Sets with Back Check Valve

- Provides the flexibility to infuse more than one bag of medication or solution
- Featuring a back check valve to prevent fluid exchange or backflow into the primary infusion line

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
MRSA15E-96	15	96 in. (244 cm)	Vented/Non-vented, Back Check Valve, (2) MicroClave® Clear Needle-free Y-sites, Pre-attached Ext Set, Slide Clamp	12.1 mL	50/Case
					
MRA15B-105	15	105 in. (267 cm)	Vented/Non-vented, Back Check Valve, (2) MicroClave Clear Needle-free Y-sites	13.2 mL	50/Case
					
AS3100 ©	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (2) Injection Y-sites	19.7 mL	50/Case
					
AS3101 ©	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, Slide Clamp, (2) Sure-Lok™ Needle-free Y-sites, Roller Clamp	20.8 mL	50/Case
					
AS3102 ©	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, Slide Clamp, (2) Pre-pierced Y-sites, Roller Clamp	20.4 mL	50/Case
					
AS3110A ©	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (3) Sure-Lok Needle-free Y-sites, Roller Clamp, Slide Clamp	22.3 mL	50/Case
					
AS3110B ©	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (3) Pre-pierced Y-sites, Slide Clamp	22 mL	50/Case
					





60 Drops/mL General Primary Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
607201 ©	60	72 in. (183 cm)	Injection Y-site	14 mL	50/Case
					
608304	60	89 in. (226 cm)	Non-vented, Roller Clamp (2), Sure-Lok® Needle-free Y-site with Pre-attached Ext Set with Pinch Clamp	17.5 mL	50/Case
					
608305 ©	60	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site with Pre-attached Ext Set, Pinch Clamp	17.5 mL	50/Case
					
608306 ©	60	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site	16 mL	50/Case
					
MRSA60-84	60	84 in. (213 cm)	Vented/Non-vented, Pre-pierced Y-site, (2) MicroClave® Clear Needle-free Y-sites	10.7 mL	50/Case
					
MRSA60E-88	60	88 in. (224 cm)	Vented/Non-vented, Pre-pierced Y-site with Pre-attached Ext Set, (2) MicroClave Clear Needle-free Y-sites	11.1 mL	50/Case
					
MRSA60-100	60	100 in. (254 cm)	Vented/Non-vented, MicroClave Clear Needle-free Y-site, Pre-pierced Y-site	12.5 mL	50/Case
					


60 Drops/mL Piggyback Sets with Back Check Valve

- Provides the flexibility to infuse more than one bag of medication or solution
- Featuring a back check valve to prevent fluid exchange or backflow into the primary infusion line


REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
MRA60B-92	60	92 in. (234 cm)	Vented/Non-vented, Back Check Valve, (2) MicroClave Clear Needle-free Y-sites, Slide Clamp	11.4 mL	50/Case
					
MRSA602-100	60	100 in. (254 cm)	Vented/Non-vented, Back Check Valve, (2) MicroClave Clear Needle-free Y-sites	12.5 mL	50/Case
					

60 Drops/mL Burette Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
MRS60B-106 ©	60	106 in. (269 cm)	150 mL Burette Set (w/MicroClave®, Shut-Off), (2) MicroClave Clear Needle-free Y-sites	12.7 mL	25/Case



MRS60B-117 ©	60	117 in. (297 cm)	24-hour 150 mL Burette Set (w/15 Micron Filter, Shut-off) with (3) MicroClave Clear Needle-free Y-sites	13.5 mL	25/Case
---------------------	----	------------------	--	---------	---------



AMSafe-3® Selectable Drip Chamber

10, 15 and 60 Drops/mL

AMSafe-3 is a 3-in-1 multi-drip chamber IV set that allows providers to select a 10, 15 or 60 drops/mL setting. Designed for efficiency where inventory reduction and cost containment are essential.

- Drop setting can be adjusted during medication delivery
- No need for clamps above and below chamber
- Includes male rotating luer lock
- Sterile, DEHP-free, not manufactured with natural rubber latex



REF	DROPS/mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
D38301 ©	10, 15 or 60	83 in. (210 cm)	Pre-pierced Y-site, Sure-Lok® Needle-free Y-site	27.7 mL	50/Case



D38302 ©	10, 15 or 60	91 in. (231 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site with Pre-attached Ext Set with Pre-pierced Y-site, Slide Clamp, Sure-Lok Needle-free Y-site	30.2 mL	50/Case
-----------------	--------------	-----------------	---	---------	---------




D38305 ©	10, 15 or 60	91 in. (231 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site, Back Check Valve with Pre-attached Ext Set with Pre-pierced Y-site, Slide Clamp, Sure-Lok Needle-free Y-site	30.5 mL	50/Case
-----------------	--------------	-----------------	--	---------	---------




Secondary Piggyback Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AS3515 ©	15	40 in. (102 cm)	Vented/Non-vented, Roller Clamp, Hanger	10 mL	50/Case
AS3500 ©	20	40 in. (102 cm)	Vented/Non-vented, Roller Clamp, Hanger	9.2 mL	50/Case




Y-Type Blood Sets


REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
MRSB-94 ©	10	94 in. (239 cm)	170 Micron Filter, Roller Clamps, Hand Pump, Slide Clamp, Roller Clamp, MicroClave® Y-site	57.2 mL	25/Case



MRSB-102	20	102 in. (259 cm)	200 Micron Filter, Roller Clamps, Hand Pump, MicroClave Y-site	20.2 mL	25/Case
-----------------	----	------------------	--	---------	---------



MRSB-105	10	105 in. (267 cm)	170 Micron Filter, Roller Clamps, Hand Pump, Pre-pierced Y-site, Slide Clamp, Roller Clamp, MicroClave Y-site, Pre-attached Ext Set	59.3 mL	25/Case
-----------------	----	------------------	---	---------	---------



IV Extension Sets

Amsino's extension sets are available in various lengths and dimensions to accommodate clinical applications in healthcare settings. All extension sets are DEHP-free, not manufactured with natural rubber latex, sterile and include a rotating male luer lock unless otherwise noted.







AE3108

Standard Extension Sets




REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
000606	6 in. (15 cm)	Female Luer Lock, Pre-pierced Y-site, Sure-Lok® Needle-free Y-site, Slide Clamp	2.5 mL	50/Case
				
MR4T07 ©	6 in. (15 cm)	MicroClave® Clear Needle-free Connector, T-connector	0.86 mL	50/Case
				
AE3107 ©	7 in. (18 cm)	Removeable Sure-Lok Needle-free Connector, Slide Clamp	1.76 mL	50/Case
				
MR4107 ©	7 in. (18 cm)	Removable MicroClave Clear Needle-free Connector, Slide Clamp	0.69 mL	50/Case
				
MR4207 ©	7 in. (18 cm)	Bonded MicroClave Clear Needle-free Connector, Slide Clamp	0.69 mL	50/Case
				
AE3108 ©	8 in. (20 cm)	Removeable Sure-Lok Needle-free Connector, Pinch Clamp	1.76 mL	100/Case
				
MR4R10 ©	10 in. (25 cm)	Removable MicroClave Clear Needle-free Connector, Pinch Clamp	0.98 mL	50/Case
				
003201 ©	27 in. (69 cm)	Female Luer Lock, (2) Slide Clamps, (2) Pre-pierced Y-sites	4.9 mL	50/Case
				
003202	27 in. (69 cm)	Female Luer Lock, (2) Slide Clamps, (2) Sure-Lok Needle-free Y-sites	5.3 mL	50/Case
				
MR4Y32 ©	32 in. (81 cm)	(2) MicroClave Clear Needle-free Connectors, (2) Slide Clamps	5.02 mL	50/Case
				

IV Extension Sets

Microbore Extension Sets



REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AE0106 ©	6 in. (15 cm)	Removeable Sure-Lok® Needle-free Connector, Slide Clamp	0.63 mL	50/Case
				
AE1106 ©	6 in. (15 cm)	Female Luer Lock, Slide Clamp, Fixed Male Luer Lock	0.5 mL	50/Case
				
AE0107 ©	7 in. (18 cm)	Female Luer Lock, Slide Clamp, Sure-Lok Needle-free Y-site	0.9 mL	50/Case
				
AE1112 ©	12 in. (31 cm)	Female Luer Lock, Slide Clamp	0.7 mL	50/Case
AE1138 ©	38 in. (97 cm)	Female Luer Lock, Slide Clamp	1.0 mL	50/Case
				

Smallbore Extension Sets

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
MR4307 ©	7 in. (18 cm)	Bonded MicroClave® Clear Needle-free Connector, Slide Clamp	0.24 mL	50/Case
				
MR4B07 ©	7 in. (18 cm)	Bifuse Ext Set with (2) MicroClave Clear Needle-free Connector, (2) Slide Clamps	0.35 mL	50/Case
				
MR4TF04 ©	4 in. (10 cm)	Trifuse Ext Set with (3) MicroClave Clear Needle-free Connectors, (3) Slide Clamps	0.4 mL	50/Case
				

Flow Controller Extension Sets

Amsino's IV extension sets with flow control devices are designed to control IV fluid and medications using a 'regulator dial' to control the rate of infusion. Includes a rotating male luer lock.

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AFS107 ©	20 in. (51 cm)	Flow Controller, Injection Y-site, Slide Clamp	3.5 mL	50/Case
				
MRSFC-19 ©	19 in. (48 cm)	Flow Controller, MicroClave Clear Needle-free Connector, Clamp	3.2 mL	50/Case
				

Pressure Infusion Extension Sets

Amsino's single-patient-use pressure infusion IV extension set is made to withstand up to 400 psig for rapid administration of medication or fluids in an emergency. Includes a rotating male luer lock.

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
MR4P07 ©	7 in. (18 cm)	Pressure Infusion (400psig) Ext Set with MicroClave® Clear Needle-free Connector, Purple Slide Clamp	0.43 mL	50/Case




Needle-Free Connectors

Our AMSafe® needle free valves are manufactured to provide both clinicians and patients safe procedural outcomes.

- Swabbable
- MRI-safe
- Not limited to 24-hour use
- Sterile, DEHP-free, not manufactured with natural rubber latex

REF	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
MR4001	MicroClave Clear Needle-free Pressure Infusion (400psig) Connector	0.04 mL	100/Case
AY0200 ©	Sure-Lok® Needle-free Connector	0.33 mL	100/Case
AY0206 ©	NeuFlo® Neutral Pressure Valve	0.08 mL	100/Case



Bag Spike

REF	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
MBSS001	MicroClave Needle-free Connector IV Bag Spike	0.34 mL	50/Case



Needle-Free Vial Adapter

REF	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
MB4V01	MicroClave Clear Needle-free Connector Single Vial Adapter	0.16 mL	50/Case



IV Pump Administration Sets

Sure-Lok® Needle-Free IV Pump Administration Sets


- FDA cleared IV pump sets for use with Baxter Spectrum and Flo-Gard® 6201 and 6301 series pumps
- Sure-Lok Needle-free connector and Pre-Pierced Y-injection site options
- 0.22 micron filter option available
- TUBETECH® technology provides +/- 5% accurate fluid delivery over 24 hour timeframe
- DEHP-free, non-pyrogenic, not manufactured with natural rubber latex

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
FP092	10	92 in. (234 cm)	Roller Clamp, Rotating Male Luer Lock with Retractable Collar	14 mL	50/Case
					
FP092Y1	10	92 in. (234 cm)	1 Needle-free Valve, Roller Clamp, 1 Slide Clamp, Rotating Male Luer Lock with Retractable Collar	14 mL	50/Case
					
FP100BY2	10	100 in. (254 cm)	Back check Valve, 2 Needle-free Valves, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	16 mL	50/Case
					
FP109BY3	10	109 in. (277 cm)	Back check Valve, 3 Needle-free Valves, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	18 mL	50/Case
					
FP169BY3	60	109 in. (277 cm)	Needle-free Valve, Roller Clamp, 2 Slide Clamps, Back check Valve, Rotating Male Luer Lock with Retractable Collar	18 mL	50/Case
					
FP106FY2	10	106 in. (269 cm)	Back check Valve, 2 Needle-free Valves, Roller Clamp, 2 Slide Clamps, 0.22 Micron Filter, Rotating Male Luer lock with Retractable Collar	22 mL	50/Case
					
FP100FY1	10	100 in. (254 cm)	1 Needle-free Valve, Roller Clamp, 2 Slide Clamps, 0.22 Micron Filter, Rotating Male Luer lock with Retractable Collar	22 mL	50/Case
					


Pre-Pierced Y-Injection IV Pump Administration Sets

- 510(k) cleared IV pump sets for use with Baxter Spectrum and Flo-Gard® 6201 and 6301 series pumps


REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
FP100BI2	10	100 in. (254 cm)	Back check Valve, 2 Pre-pierced Injection Y-sites, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	16 mL	50/Case



FP109BI3	10	109 in. (277 cm)	Back check Valve, 3 Pre-pierced Injection Y-Sites, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	18 mL	50/Case
-----------------	----	------------------	--	-------	---------



FP106FI2	10	106 in. (269 cm)	Back check Valve, 2 Pre-pierced Injection Y-sites, Roller Clamp, 2 Slide Clamps, 0.22 Micron Filter, Rotating Male Luer lock with Retractable Collar	22 mL	50/Case
-----------------	----	------------------	--	-------	---------



AMSafe® Prefilled Normal Saline Flush Syringe

- 0.9% Sodium Chloride Injection, USP
- 3 mL, 5 mL, 10 mL and 20 mL
- Sterile Field Ready (SFR) and Sterile Fluid Path packaging configurations



For more details and ordering information, please refer to page 55.



UROLOGICAL CARE

Solutions for Infection Prevention

AMSure Urological Care product portfolio provides solutions to help improve patient safety and infection prevention. Our product selection includes:

- Foley Catheter Trays & Collection 22**
 - Pre-Connected Foley Catheter Urine Drainage Bag Trays..... 22
 - Add-a-Foley Catheter Trays 23
 - Pre-Connected Foley Catheter Urine Meter Kits 24
 - Add-a-Foley Urine Meter Trays 25
 - Foley Insertion Trays..... 25
 - Urinary Drainage Bags 26
 - Urine Meter 27
 - ACCU-MEASURE® Pediatric Urine Meter 27
- Foley Catheters 28**
 - Hydrophilic Latex Foley Catheters 28
 - Silicone-Coated Latex Foley Catheters 29
 - 100% Silicone Foley Catheters 30
- Intermittent Urinary Catheters 31**
 - Urethral Insertion Trays and Kits 31
 - Closed System Urethral Drainage Kits 32
 - Self-Catheterization Kits with R-Polished™ Eye Catheters 33
 - PVC Intermittent Urethral Catheter with R-Polished™ Eyes..... 34
 - PVC Intermittent Urethral Catheter 35
 - Red Rubber Latex Urethral Catheter 35
- Irrigation Sets 36**
 - Cystoscopy & TUR Irrigation Sets 36
- Ancillary & Patient Care Products..... 37**
 - Urinary Leg Bags 37
 - Extension Tube for Leg Bags 37
 - Pediatric Urine Collector 37
 - Irrigation Trays and Kits 38
 - Syringes 39
 - Specimen Containers 40

Foley Catheter Trays & Collection



AMSure® Pre-Connected Foley Catheter Urine Drainage Bag Trays

AMSure® Pre-Connected Drainage Bag Foley Trays improve simplicity and efficiency of care. The catheter is pre-attached to a low-profile drainage bag which helps to minimize potential contamination. (HCPCS codes: See product chart below)

- Tray includes all necessary components for aseptic catheter insertion
- Foley catheter connected with tamper evident seal to prevent disconnection and maintain sterile fluid path
- Tray improves clinical efficacy through cost-effective inventory management
- Gloves and povidone iodine solution help facilitate aseptic technique and reduce the risk of infection
- Needle-free urinary sampling port reduces the risk of needle stick injuries
- Trays are available with 100% silicone, silicone-coated latex, or hydrophilic-coated latex Foley catheters to suit patient needs
- Sterile

Trays include: 2000 mL low-profile urinary drainage bag, pre-connected 16 in. Foley catheter with 5 cc balloon, waterproof and fenestrated drapes, Nitrile powder-free gloves, (5) Rayon balls, 10% Povidone Iodine solution, forceps, lubricating jelly, syringe prefilled with 10 cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap and Tyvek header bag

CATHETER TYPE	REF	FOLEY SIZE	BALLOON SIZE	QTY/UNIT
100% Silicone (HCPCS Code A4315)	AS89314S	14 Fr.	5 cc	10/Case
	AS89316S	16 Fr.	5 cc	10/Case
	AS89318S	18 Fr.	5 cc	10/Case
Silicone Coated Latex (HCPCS Code A4314)	AS89314	14 Fr.	5 cc	10/Case
	AS89316	16 Fr.	5 cc	10/Case
	AS89318	18 Fr.	5 cc	10/Case
Hydrophilic Coated Latex (HCPCS Code A4314)	AS89314H	14 Fr.	5 cc	10/Case
	AS89316H	16 Fr.	5 cc	10/Case
	AS89318H	18 Fr.	5 cc	10/Case

Foley Catheter Trays & Collection



AS89300

AMSure® Add-a-Foley Catheter Trays

AMSure Add-a-Foley Catheter Trays improve simplicity and efficiency of care by including the urine collection device and necessary supplies to safely insert a Foley catheter. (HCPCS code A4354)

- **Needle-free sampling port on urinary bag reduces the risk of needle stick injuries**
- **Gloves and povidone iodine swab sticks help maintain sterility and minimized risk of infection**
- **Sterile, not manufactured with natural rubber latex**

Trays include: 2000 mL low-profile urinary drainage bag, waterproof and fenestrated drapes, Nitrile powder-free gloves, (5) Rayon balls, 10% povidone iodine solution, forceps, lubricating jelly, 2 oz. specimen container with patient label, central supply room wrap and Tyvek header bag

REF	DESCRIPTION	PRE-FILLED SYRINGE	QTY/UNIT
AS89300	Add-a-Foley Tray 2000 mL Urine Drainage Bag	10 cc	10/Case
AS89330 ©	Add-a-Foley Tray 2000 mL Urine Drainage Bag	30 cc	10/Case

AMSure® Pre-Connected Foley Catheter Urine Meter Kits

The AMSure pre-connected urine meter trays offer the benefit of a closed system along with the accuracy of a meter to measure exact urine output. These truly closed systems provide all the necessary supplies for an aseptic catheter insertion to promote the US Centers for Disease Control (CDC) recommendations to help prevent catheter associated urinary tract infections (CAUTIs). (HCPCS codes: See product chart below)

Trays include: 400 mL Urine Meter Preconnected with 5 cc Foley Catheter, waterproof and fenestrated drapes, (2) Nitrile gloves, 10% Povidone Iodine solution, (5) Rayon balls, forceps, lubricating jelly, syringe prefilled with 10 cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap, Tyvek header bag



AS89814S

CATHETER TYPE	REF	FOLEY SIZE	BALLOON SIZE	QTY/UNIT
100% Silicone (HCPCS Code A4315)	AS89814S ©	14 Fr.	5 cc	10/Case
	AS89816S ©	16 Fr.	5 cc	10/Case
	AS89818S ©	18 Fr.	5 cc	10/Case
Silicone Coated Latex (HCPCS Code A4314)	AS89814 ©	14 Fr.	5 cc	10/Case
	AS89816 ©	16 Fr.	5 cc	10/Case
	AS89818 ©	18 Fr.	5 cc	10/Case
Hydrophilic Coated Latex (HCPCS Code A4314)	AS89814H ©	14 Fr.	5 cc	10/Case
	AS89816H ©	16 Fr.	5 cc	10/Case
	AS89818H ©	18 Fr.	5 cc	10/Case

Foley Catheter Trays & Collection

AMSure® Add-a-Foley Urine Meter Trays

The AMSure Add-a-Foley urine meter kits provide all the supplies necessary for an aseptic insertion. They also offer the caregiver the flexibility to choose the best Foley catheter for the patient and the situation, along with providing a urine meter to obtain the most accurate urine output measurements. (HCPCS Code A4354)



AS89800

- **Sterile, not manufactured with natural rubber latex**

Trays include: 400 ml Urine Meter Preconnected with 5 cc Foley Catheter, waterproof and fenestrated drapes, (2) Nitrile gloves, 10% Povidone Iodine solution, (5) Rayon balls, forceps, lubricating jelly, syringe prefilled with 10 cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap, Tyvek header bag

REF	DESCRIPTION	PRE-FILLED SYRINGE	QTY/UNIT
AS89800	Add-a-Foley Tray 400 mL Urine Meter with 2600 mL Urine Drainage Bag	10 cc	10/Case
AS89830	Add-a-Foley Tray 400 mL Urine Meter with 2600 mL Urine Drainage Bag	30 cc	10/Case

AMSure® Foley Insertion Trays

AMSure Foley Insertion Trays improve simplicity and efficiency of care by providing the basic supplies required to safely insert a Foley catheter. (HCPCS code A4310)



AS880

- **Deep tray basin designed to prevent contamination**
- **Gloves and povidone iodine swab sticks help maintain sterility and minimized risk of infection**
- **Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex**

Trays include: 1000 mL collection basin, (3) Povidone-iodine swab sticks or benzalkonium chloride prep pad, vinyl powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), lubrication jelly, fenestrated drape

SWAB STICK	REF	PRE-FILLED SYRINGE	QTY/UNIT
PVP Swab Sticks	AS880 ©	10 cc	10/Case
	AS890 ©	30 cc	10/Case
BZK Swab Sticks	AS880K ©	10 cc	10/Case
	AS890K ©	30 cc	10/Case

AMSure® Urinary Drainage Bags

AMSure urine drainage bags are available in a variety of different configurations to meet the clinicians' needs to best provide patient safety and comfort while helping to minimize infections. All drainage bags include a sheet clamp for the positioning and securing of the drainage tube. The drainage tube has an ID of 11/32" (0.87 cm) for unrestricted flow and are 4 feet (1.2 M) in length, terminating with a standard terraced catheter fitting for a secure and leak-proof connection. Each bag includes a storage pocket for the emptying tube. (HCPCS code A4358)

- 2,000 cc capacity graduated in 200 cc increments
- Anti-reflux chambers help improve patient safety by restricting urine backflow
- Leak-proof seams ensure fluid containment
- Needle-free sampling port reduces the risk of needle stick injuries
- Air vents help to maintain consistent urine flow, delivering safer patient outcomes
- Sterile fluid pathway, BPA-free, not manufactured with natural rubber latex



AS322 Teardrop



AS326 Oval



AS32600 Low-Profile

REF	BAG SHAPE	ANTI-REFLUX TYPE	NEEDLE-FREE SAMPLING PORT	DRAIN TYPE	HANGER TYPE	QTY/UNIT
AS312 ©	Diamond	Drip Device	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	20/Case
AS322 ©	Teardrop	Drip Device	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	20/Case
AS332	Teardrop	Small Chamber	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	20/Case
AS302 ©	Pear	Flutter Valve	Pre-Pierced Blunt Cannula	T-Tap	Single-Hook	20/Case
AS32200 ©	Low-Profile	Large Chamber	Pre-Pierced Blunt Cannula	T-Tap	Swivel Single-Hook & Rope	20/Case
AS32600* ©	Low-Profile	Large Chamber	Male Luer-Lock or Luer-Slip	T-Tap	Swivel Single-Hook & Rope	20/Case
AS326* ©	Oval	Drip Device	Pre-Pierced Blunt Cannula	Metal Clamp with Doghouse Storage	Double-Hook & Rope	20/Case

*100% sterile packed in peel pouch

Foley Catheter Trays & Collection

AMSure® Urine Meter

AMSure Urine Meter Drain Bags are designed to provide all caregivers the ability to easily and accurately measure urine output. The meter device is pre-attached to a urine drainage bag, for convenient draining of the meter while still obtaining incremental urine measurements. (HCPCS code A4358)



- **Anti-Reflux Chamber helps improve patient safety by reducing urine backflow**
- **Twist valve drains at the bottom of the urine meter for easy drainage**
- **Leak-proof seams ensure fluid containment**
- **Needle-free sampling port reduces the risk of needle stick injuries**
- **Air vents help to maintain consistent urine flow delivering safer patient outcomes**
- **Sterile, not manufactured with natural rubber latex**

REF	DESCRIPTION	QTY/UNIT
AS328 ©	400 mL Urine Meter with 2600 mL Urine Drainage Bag	10/Case

ACCU-MEASURE® Pediatric Urine Meter

2000 mL Bag with 150 mL Chamber

ACCU-MEASURE Pediatric Urine Meter is specifically designed for pediatric patients. With 1 mL increments, the 150 mL meter chamber's white background provides accurate and precise output measurements. The needle-free sample port reduces the risk of needle stick injuries, helping keep nurses safe while making their job easier. (HCPCS code A4358)



- **Double outlet valves – Twist valve for taking sample, T-valve for emptying the bag**
- **Needle-free sampling port with luer lock access**
- **Leak-proof seams ensure fluid containment**
- **Meter bottle capacity 150 mL**
- **Adjustable bed sheet clamp for securing drainage bag tubing**
- **Sterile, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	QTY/UNIT
AS327100 ©	ACCU-MEASURE® Pediatric Urine Meter 2000 mL Bag with 150 mL Chamber	30/Case



Foley Catheters

AMSure® Hydrophilic Latex Foley Catheters

(HCPCS code A4338)

- Hydrophilic coating provides for lubricious surface when contacting moisture to help reduce irritation of the urethra
- Smooth eyes for improved patient comfort
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Sterile, DEHP-free, BPA-free



BALLOON SIZE	REF	FOLEY SIZE	QTY/UNIT
5 cc	AS41014H ©	14 Fr. x 16 in. (40.6 cm)	10/Box
	AS41016H ©	16 Fr. x 16 in. (40.6 cm)	10/Box
	AS41018H ©	18 Fr. x 16 in. (40.6 cm)	10/Box
30 cc	AS42014H	14 Fr. x 16 in. (40.6 cm)	10/Box
	AS42016H ©	16 Fr. x 16 in. (40.6 cm)	10/Box
	AS42018H ©	18 Fr. x 16 in. (40.6 cm)	10/Box

Foley Catheters

AMSure® Silicone-Coated Latex Foley Catheters

(HCPCS code A4338)

- Silicone coating helps reduce urethral irritation
- Smooth beveled symmetrical eyes minimize tissue damage
- Large oval eyes aid in maximum drainage
- Firm catheter tip designed to improve ease of insertion
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Individually packaged, sterile, DEHP-free, BPA-free



BALLOON SIZE	REF	FOLEY SIZE	QTY/UNIT
3 cc	AS40008 ©	8 Fr. × 10.6 in. (27cm)	10/Box
	AS40010 ©	10 Fr. × 10.6 in. (27cm)	10/Box
5 cc	AS41012 ©	12 Fr. × 16 in. (40.6cm)	10/Box
	AS41014 ©	14 Fr. × 16 in. (40.6cm)	10/Box
	AS41016 ©	16 Fr. × 16 in. (40.6cm)	10/Box
	AS41018 ©	18 Fr. × 16 in. (40.6cm)	10/Box
	AS41020 ©	20 Fr. × 16 in. (40.6cm)	10/Box
	AS41022 ©	22 Fr. × 16 in. (40.6cm)	10/Box
	AS41024 ©	24 Fr. × 16 in. (40.6cm)	10/Box
	AS41026 ©	26 Fr. × 16 in. (40.6cm)	10/Box
	AS41028 ©	28 Fr. × 16 in. (40.6cm)	10/Box
	AS41030 ©	30 Fr. × 16 in. (40.6cm)	10/Box
30 cc	AS42012 ©	12 Fr. × 16 in. (40.6cm)	10/Box
	AS42014 ©	14 Fr. × 16 in. (40.6cm)	10/Box
	AS42016 ©	16 Fr. × 16 in. (40.6cm)	10/Box
	AS42018 ©	18 Fr. × 16 in. (40.6cm)	10/Box
	AS42020 ©	20 Fr. × 16 in. (40.6cm)	10/Box
	AS42022 ©	22 Fr. × 16 in. (40.6cm)	10/Box
	AS42024 ©	24 Fr. × 16 in. (40.6cm)	10/Box
	AS42026 ©	26 Fr. × 16 in. (40.6cm)	10/Box
	AS42028 ©	28 Fr. × 16 in. (40.6cm)	10/Box
	AS42030 ©	30 Fr. × 16 in. (40.6cm)	10/Box

AMSure® 100% Silicone Foley Catheters

(HCPCS code A4344)

- Eliminates latex allergy concerns
- Smooth eyes improve patient comfort
- Reinforced tip provides a stylet for optimized insertion efficacy
- 100% silicone material has non-adherent characteristics and reduces encrustation and calcifications
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Individually packaged, sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



BALLOON SIZE	REF	FOLEY SIZE	QTY/UNIT
5 cc	AS41012S ©	12 Fr. × 16 in. (40.6cm)	10/Box
	AS41014S ©	14 Fr. × 16 in. (40.6cm)	10/Box
	AS41016S ©	16 Fr. × 16 in. (40.6cm)	10/Box
	AS41018S ©	18 Fr. × 16 in. (40.6cm)	10/Box
	AS41020S ©	20 Fr. × 16 in. (40.6cm)	10/Box
	AS41022S ©	22 Fr. × 16 in. (40.6cm)	10/Box
	AS41024S ©	24 Fr. × 16 in. (40.6cm)	10/Box
	AS41026S ©	26 Fr. × 16 in. (40.6cm)	10/Box
30 cc	AS41028S ©	28 Fr. × 16 in. (40.6cm)	10/Box
	AS42012S ©	12 Fr. × 16 in. (40.6cm)	10/Box
	AS42014S ©	14 Fr. × 16 in. (40.6cm)	10/Box
	AS42016S ©	16 Fr. × 16 in. (40.6cm)	10/Box
	AS42018S ©	18 Fr. × 16 in. (40.6cm)	10/Box
	AS42020S ©	20 Fr. × 16 in. (40.6cm)	10/Box

Intermittent Urinary Catheters



AS870

AMSure® Urethral Insertion Trays and Kits

AMSure® Urethral Trays improve simplicity and efficiency of care by including all components needed to perform short-term aseptic catheterization. (HCPCS code A4353)

- **Gloves and povidone iodine swab sticks help maintain aseptic technique and reduce the risk of infection**

PVC Catheter Tray

- **Staggered eyes help maintain tip strength during insertion and promote efficient drainage**
- **Short-term catheter use helps decrease infections**
- **Sterile, BPA-free, not manufactured with natural rubber latex**

Red Rubber Catheter Trays

- **Smooth tip and quality rubber latex promote patient comfort during insertion**
- **Clearly marked French size on funnel promotes ease of selection**
- **Sterile, DEHP-free, BPA-free**

Trays include: 1000 mL outer basin tray, 16 in. urethral catheter, waterproof and fenestrated drapes, vinyl powder-free gloves, (3) 10% povidone iodine swab sticks, lubricating jelly, 4 oz. specimen container

REF	URETHRAL CATHETER TYPE	QTY/UNIT
AS870 ©	14 Fr PVC*	20/Case
AS87114 ©	14 Fr Red Rubber Latex	20/Case
AS87118 ©	18 Fr Red Rubber Latex	20/Case

*Not manufactured with natural rubber latex.

AMSure® Closed System Urethral Drainage Kits

AMSure® Pre-connected Drainage Bag Urethral Catheter Kits improve simplicity and efficiency of care. The catheter is pre-attached to the drainage bag which helps to minimize contamination. (HCPCS code A4353)

- **Gloves and povidone iodine swab sticks help maintain sterility and reduce the risk of infection**
- **Sterile**

Red Rubber Catheter Kits

- **Smooth tip and quality rubber latex promote patient comfort during insertion**
- **Clearly marked French size on funnel promotes ease of selection**
- **Sterile, DEHP-free, BPA-free**

PVC Catheter Kits

- **Staggered eyes help maintain tip strength during insertion and promote efficient drainage**
- **Short-term catheter use helps decrease infections**
- **Sterile, BPA-free, not manufactured with natural rubber latex**



AS89214R

Kit includes: 1200 mL urinary drainage bag, 16 in. (40.64 cm) urethral catheter, CSR wrap, waterproof and fenestrated drapes, (2) Nitrile powder-free gloves, (3) Povidone-iodine swab sticks, lubricating jelly, 2 oz. specimen container with ID label

SWAB STICK	REF	FOLEY SIZE	FOLEY MATERIAL	QTY/UNIT
PVP Swab Sticks	AS89214R ©	14 Fr. × 16 in. (40.6cm)	Red Rubber (Contains Latex)	20/Case
	AS89214 ©	14 Fr. × 16 in. (40.6cm)	PVC Latex- Free	20/Case
BZK Swab Sticks	AS89214RK ©	14 Fr. × 16 in. (40.6cm)	Red Rubber (Contains Latex)	20/Case
	AS89216RK ©	16 Fr. × 16 in. (40.6cm)	Red Rubber (Contains Latex)	20/Case

Intermittent Urinary Catheters



AMSure® Self-Catheterization Kits with R))POLISHED™ Eye Catheters

AMSure Self-Catheterization kits with R-Polish eyelets technology allow the patient, even with limited dexterity, to perform incontinence care with confidence. R-Polish technology uses radio frequency to create ultra-smooth eyelets for optimal patient comfort.

The pre-lubricated, closed system with introducer tip helps to reduce potential for infections. The kits contain all the necessary components needed to perform aseptic self-catheterization. (HCPCS code A4353)

- **Advanced R-Polished eyelets enhance patient comfort during catheter insertion and removal**
- **Introducer tip reduces the potential for urinary tract infection (UTI)**
- **Collection bag with built-in thumb holes allows for ease of use**
- **Sterile, BPA-free, not manufactured with natural rubber latex**

Kit includes: Pre-lubricated 16 in. (40.64 cm) long PVC urethral catheter, 1100 mL collection bag with insertion tip and protector cap, (3) Povidone-iodine swab sticks and benzalkonium chloride prep pad, vinyl powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), waste disposal bag, gauze pad 4 in. x 4 in. (10.16 cm x 10.16 cm)

SWAB STICK	REF	CATHETER SIZE	QTY/UNIT
PVP Swab Sticks	AS85008 ©	8 Fr. x 16 in. (40.6cm)	50/Case
	AS85010 ©	10 Fr. x 16 in. (40.6cm)	50/Case
	AS85012 ©	12 Fr. x 16 in. (40.6cm)	50/Case
	AS85014 ©	14 Fr. x 16 in. (40.6cm)	50/Case
	AS85016 ©	16 Fr. x 16 in. (40.6cm)	50/Case
BZK Swab Sticks	AS85008K ©	8 Fr. x 16 in. (40.6cm)	50/Case
	AS85010K ©	10 Fr. x 16 in. (40.6cm)	50/Case
	AS85012K ©	12 Fr. x 16 in. (40.6cm)	50/Case
	AS85014K ©	14 Fr. x 16 in. (40.6cm)	50/Case
	AS85016K ©	16 Fr. x 16 in. (40.6cm)	50/Case

AMSure® PVC Intermittent Urethral Catheter with R-POLISHED™ Eyelets

AMSure PVC Intermittent Catheters with R-Polish technology use radio frequency to create ultra-smooth eyelets for optimal patient comfort. They are designed to provide a smooth surface in a catheter that is firm enough for insertion, yet soft enough to be slightly pliable. These urethral catheters are designed for intermittent or short-term use and are available in various French sizes, as well as female, pediatric and male lengths. A Coudé tip is also available to assist in more difficult catheterizations where strictures or other anomalies may be present in the urethra. (HCPCS code A4351)

- **Advanced R-Polished eyelets enhance patient comfort during catheter insertion and removal**
- **Staggered eyes help maintain tip strength during insertion and promote efficient drainage**
- **Designed to minimize trauma and reduce risk of infection**
- **PVC material eliminates latex allergy concerns**
- **Clear shaft provides easy viewing as you measure**
- **Short-term catheter use helps decrease infections**
- **Individually packaged, sterile**
- **BPA-free, not manufactured with natural rubber latex**



FEATURES	REF	DESCRIPTION	QTY/UNIT
Female	AS860608	8 Fr, 6 in. (15.24 cm)	50/Case
	AS860610	10 Fr, 6 in. (15.24 cm)	50/Case
	AS860612	12 Fr, 6 in. (15.24 cm)	50/Case
	AS860614	14 Fr, 6 in. (15.24 cm)	50/Case
	AS860616	16 Fr, 6 in. (15.24 cm)	50/Case
Pediatric	AS861006	6 Fr, 10 in. (25.40 cm)	50/Case
	AS861008	8 Fr, 10 in. (25.40 cm)	50/Case
	AS861010	10 Fr, 10 in. (25.40 cm)	50/Case
Male	AS861608	8 Fr, 16 in. (40.64 cm)	50/Case
	AS861610	10 Fr, 16 in. (40.64 cm)	50/Case
	AS861612	12 Fr, 16 in. (40.64 cm)	50/Case
	AS861614	14 Fr, 16 in. (40.64 cm)	50/Case
	AS861616	16 Fr, 16 in. (40.64 cm)	50/Case
	AS861618	18 Fr, 16 in. (40.64 cm)	50/Case
Coudé Tip	AS861610C	10 Fr, 16" (40.65 cm)	30/Case
	AS861612C	12 Fr, 16" (40.65 cm)	30/Case
	AS861614C	14 Fr, 16" (40.65 cm)	30/Case
	AS861616C	16 Fr, 16" (40.65 cm)	30/Case
	AS861618C	18 Fr, 16" (40.65 cm)	30/Case

Intermittent Urinary Catheters

AMSure Red Rubber Urethral Catheter

AMSure Red Rubber Urethral Catheters are pliable and designed for intermittent use. The urethral catheters are available in various French sizes as well as female and male lengths. (HCPCS code A4351)

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Clearly marked French size on funnel promotes ease of selection
- Individually packaged, sterile, DEHP-free




REF	DESCRIPTION	QTY/UNIT
AS44008 ©	8 Fr, 11 in. (27.94 cm)	100/Case
AS44010 ©	10 Fr, 11 in. (27.94 cm)	100/Case
AS44012 ©	12 Fr, 16 in. (40.64 cm)	100/Case
AS44014 ©	14 Fr, 16 in. (40.64 cm)	100/Case
AS44016 ©	16 Fr, 16 in. (40.64 cm)	100/Case
AS44018 ©	18 Fr, 16 in. (40.64 cm)	100/Case
AS44020 ©	20 Fr, 16 in. (40.64 cm)	100/Case
AS44022 ©	22 Fr, 16 in. (40.64 cm)	100/Case
AS44024 ©	24 Fr, 16 in. (40.64 cm)	100/Case
AS44028 ©	28 Fr, 16 in. (40.64 cm)	100/Case

Irrigation Sets

AMSure® Cystoscopy & TUR Irrigation Sets

The AMSure Cystoscopy and TUR Irrigation sets are used for continuous bladder irrigation after prostate or bladder surgeries. These are available with one, two, or four leads to meet a variety of needs. (HCPCS code A4351)

- Robert Clamp(s)
- Drip Chamber(s)
- Sterile, not manufactured with natural rubber latex



REF	TYPE	QTY/UNIT
ABL400 ©	Single Lead	20/Case



REF	TYPE	QTY/UNIT
ABL402 ©	Two Leads	20/Case



REF	TYPE	QTY/UNIT
ABL404 ©	Four Leads	20/Case

Urinary Leg Bags

(HCPCS code A4358)

- Adjustable, soft, breathable pre-attached elastic leg straps for patient comfort
- Anti-reflux flutter valve improves patient safety by reducing urine backflow
- Leak-proof seams ensure fluid containment
- Sterile fluid pathway
- BPA-free, not manufactured with natural rubber latex



AS307N

SIZE	REF	DRAIN PORT STYLE	18 in. (45.72 cm) EXTENSION TUBE WITH ADAPTER & PROTECTIVE CAP	QTY/UNIT
600 mL	AS306N ©	Push-Pull		48/Case
	AS308N ©	Twist-Turn		48/Case
	LBE306N	Push-Pull	•	48/Case
	LBE308N	Twist-Turn	•	48/Case
900 mL	AS307N ©	Push-Pull		48/Case
	AS309N ©	Twist-Turn		48/Case
	LBE307N	Push-Pull	•	48/Case
	LBE309N	Twist-Turn	•	48/Case

Extension Tube for Leg Bags

(HCPCS code A4331)

- Kink-resistant tubing safeguards against fluid buildup
- Flexible tubing permits leg bag to be positioned on the thigh or calf
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
AS311	18 in. (45.72 cm) Extension Tube	50/Case

Pediatric Urine Collector

(HCPCS code A4328)

- Adhesive surface folds to form a spill-proof specimen transport container
- Hypoallergenic feature makes it safe to use with most pediatric patients
- Visible graduation marks designed for ease of use
- Single-use disposable
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
AS409 ©	200 mL Pediatric Urine Bag	50/Box

Irrigation Trays and Kits

AMSure® Irrigation Kits and Trays are designed for simplicity and efficiency of care. The trays are available sterile or non-sterile to suit your clinical need. (HCPCS code A4320)

Bulb Syringe

- Flexible soft bulb promotes one-handed administration
- Designed to stand on the bulb end to help minimize contamination
- Patient label provides convenient patient documentation
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector



AS121

Thumb Control Ring Syringe

- Thumb control ring top allows for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage
- Patient label provides convenient patient documentation
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip adapter and protector cap



AS126

Non-sterile irrigation kits include: 500 mL graduated container, patient ID label and tip protector

Sterile irrigation trays include: 1000 mL outer tray, 500 mL graduated container, waterproof underpad, alcohol prep pad, patient ID label and tip protector

Non-sterile Irrigation Trays

REF	THUMB CONTROL RING SYRINGE	BULB SYRINGE	RESEALABLE POLE BAG	QTY/UNIT
AS121		•		30/Case
AS126	•			30/Case
AS127	•		•	30/Case

Sterile Irrigation Trays

REF	THUMB CONTROL RING SYRINGE	BULB SYRINGE	QTY/UNIT
AS130 ©		•	20/Case
AS136 ©	•		20/Case



AS136

Ancillary & Patient Care Products

Syringes

(HCPCS code A4322)

Bulb Syringe

- Flexible soft bulb promotes one-handed administration
- Designed to stand on the bulb end to help minimize contamination
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector

Piston Thumb Control Ring Syringe

- Piston thumb control ring top allows for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe luer tip adapter and protector cap

Flat Top Syringe

- Visible dosing marks designed for ease of use
- Precision crafted O-ring plunger ensures a positive seal
- Visible dosing marks designed for ease of use
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Available in sterile or non-sterile packaging
- Includes syringe luer tip adapter and protector cap



AS010

AS016

AS115

TYPE	REF	PACKAGING	STERILITY	TIP PROTECTOR CAP	SMALL TUBE ADAPTER	QTY/UNIT
Bulb	AS011 ©	Poly Pouch	Sterile (EO)	•		50/Case
	AS010	Bulk	Non-Sterile			50/Case
Piston Thumb Control Ring	AS015 ©	Poly Pouch	Sterile (EO)	•	•	50/Case
	AS016	Resealable Hanger Bag	Non-Sterile	•	•	30/Case
Piston Flat Top	AS115 ©	Poly Pouch	Sterile (EO)	•	•	50/Case
	AS116 ©	Resealable Pole Bag	Non-Sterile	•	•	30/Case

AMSure® Specimen Containers

AMSure 4 oz., leak-resistant specimen containers are designed for the collection and safe transport of patient samples. (HCPCS code P9615)

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	PACKAGING	STERILITY	LABELS	QTY/UNIT
AS340 ©	Poly Pouch	Interior	Patient Information	100/Case
AS341 ©	Bulk	Non-Sterile	Tamper Evident	100/Case
AS343	Bulk	Non-Sterile	None	500/Case
AS345 ©	Form-Fill & Seal	Fully Sterile	None	100/Case



AS340



AMSure®

STERILE WATER AND SALINE PRODUCTS



Inhalation, Irrigation, and Injection Solutions

AMSure Respiratory Therapy Systems are designed to deliver optimal FiO₂ levels as well as keep upper airways unobstructed and clear to help improve quality of care. The product portfolio includes:

Sterile Water for Inhalation, USP	42
Sterile Water Prefilled Humidifiers, USP	42
Nebulizer Kits	69
Saline Products	43
Device and Wound Irrigation	44

AMSure® Sterile Water

All products are DEHP-Free, PVC-Free, not manufactured with natural rubber latex, preservative-free, and Non-Pyrogenic.



AS1065

AS1000

Sterile Water for Inhalation (Bottles)

- 500 mL and 750 mL, and 1000 mL bottles Sterile Water, USP
- Used in nebulizers and humidifiers

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
AS1565	AS1565E ©	500 mL Sterile Water for Inhalation, USP	12/Case
AS1765	AS1765E ©	750 mL Sterile Water for Inhalation, USP	12/Case
AS1065	AS1065E ©	1000 mL Sterile Water for Inhalation, USP	12/Case
AS1566	—	500 mL Sterile Water for Inhalation with non-sterile cap	12/Case
AS1066	—	1000 mL Sterile Water for Inhalation with non-sterile cap	12/Case
AS1000	—	1000 mL Sterile Water for Inhalation, Hanging Spikeable Bottle	12/Case

Sterile Water for Inhalation (Vials)

- Available in two sizes: 3 mL and 5 mL, USP
- Easy open twist-off cap is simple to use and allows dispensing by the drop or milliliter
- Used for suctioning and respiratory inhalation



ASR0035

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
ASR0035	ASR0035E ©	3 mL Sterile Water Inhalation, USP	1000/Case
ASR0055	ASR0055E ©	5 mL Sterile Water Inhalation, USP	1000/Case

Sterile Water Prefilled Humidifiers

- 350 mL and 550 mL bottles Sterile Water, USP for inhalation therapy and humidification
- For respiratory therapy, adapters included with each bottle
- Reservoir bottles have built-in oxygen tubing connectors
- Micro-diffuser produces smaller bubbles and greater surface agitation for quiet operation
- Audible alarm alerts clinician of flow restriction or occlusion of oxygen tubing



AS0352

AS0552

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
AS0352	AS0352E ©	350 mL Sterile Water for Inhalation Prefilled Humidifier, USP	20/Case
AS0552	AS0552E ©	550 mL Sterile Water for Inhalation Prefilled Humidifier, USP	12/Case

Saline Products

AMSafe® Prefilled Saline Flush Syringe

- 3 mL 5 mL, 10 mL and 20 mL syringes, 0.9% Sodium Chloride Injection, USP
- Sterile fluid path (dust cover) and sterile field configurations
- Complies with PICC and INS standards
- Compatible with standard ISO Luer connectors
- FDA cleared
- Sterile, non-pyrogenic, DEHP-free, PVC-free, not manufactured with natural rubber latex



12 mL and 20 mL Syringes

Standard Dust Cover Flush Syringes (Sterile Fluid Pathway)

REF	DESCRIPTION	QTY/UNIT
IVF1203 ©	3 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1205 ©	5 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1210 ©	10 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF2020 ©	20 mL USP Sodium Chloride Injection Fill in 20 mL Syringe	120/Case (6 Boxes – 20 per Box)

Sterile Field Ready Flush Syringes

REF	DESCRIPTION	QTY/UNIT
IVF1203SS ©	3 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1205SS ©	5 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1210SS ©	10 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF2020SS ©	20 mL USP Sodium Chloride Injection Fill in 20 mL Syringe	120/Case (6 Boxes – 20 per Box)



10 mL Syringes

Standard Dust Cover Flush Syringes (Sterile Fluid Pathway)

REF	DESCRIPTION	QTY/UNIT
IVF1003 ©	3 mL USP Sodium chloride Injection Fill in 10 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1005 ©	5 mL USP Sodium chloride Injection Fill in 10 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1010 ©	10 mL USP Sodium chloride Injection Fill in 10 mL Syringe	240/Case (8 Boxes – 30 per Box)
IVF1010TM	10 mL USP Sodium chloride Injection Fill in 10 mL Syringe	240/Case (8 Boxes – 30 per Box)

AMSure® Sterile Sodium Chloride Inhalation Solution (vials)

- 3 mL, 5 mL and 15 mL Sterile Sodium Chloride Inhalation Solution (0.9%), USP
- Single-use vials with graduations for accurate dispensing
- Easy open twist-off cap is simple to use and allows dispensing by the drop or milliliter
- Used for suctioning and respiratory inhalation

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
ASR0039	ASR0039E ©	3 mL Sterile Sodium Chloride Inhalation Solution (0.9%) USP	1000/Case
ASR0059	ASR0059E ©	5 mL Sterile Sodium Chloride Inhalation Solution (0.9%) USP	1000/Case
ASR0159	ASR0159E ©	15 mL Sterile Sodium Chloride Inhalation Solution (0.9%) USP	144/Case



ASR0039



ASR0059



ASR0159

AMSure Device and Wound Irrigation

- 100 mL sterile water or saline bottle with screw cap and inner seal
- Easy open twist screw top container
- Used for device and wound irrigation

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
AS399 ©	AS399E ©	100 mL Sterile Water for Device and Wound Irrigation	48/Case
AS398 ©	AS398E ©	100 mL Sterile Saline 0.9% for Device and Wound Irrigation	48/Case



AS399



AS398



AMSure Sal-Sense™ Single Use Saline Topical Solution 0.9% Sodium Chloride

- Single patient use avoids cross contamination
- Save time from standard practice of drawing solution
- Twist design to open effortlessly
- Easy to irrigate area in any position



REF	DESCRIPTION	QTY/UNIT
AS0159	15 mL Single Use Saline Topical Solution	144/Case–3 Boxes (48 per Box)
AS0359	35 mL Single Use Saline Topical Solution	120/Case–6 Boxes (20 per Box)
AS0559	50 mL Single Use Saline Topical Solution	120/Case–6 Boxes (20 per Box)

RECEPTAL®

SUCTION & WASTE FLUID MANAGEMENT SYSTEMS



Waste Fluid Management Solutions


Amsino is a leading global provider of suction fluid management systems. The breadth of Amsino's RECEPTAL suction fluid management systems enhance patient quality of care through clinical efficacy, ease-of-use, product selection and key safety features. Our product selection includes:

EZE-VAC® Suction Canister System	46
EZE-VAC 2.0 Suction Canister System	47
RECEPTAL Liners	48
RECEPTAL Twist Liners	49
RECEPTAL Outer Shells	50
ACCU-MEASURE® Critical Measurement Devices	51
RECEPTAL Solidifiers	51
AMSure® Suction Tubing with Connectors	52
Suction Yankauers and Kits	53
SteriCadi™ Suction Device Holster System	54
Suction Catheters and Trays	55
Suction Catheter Kits	56
Accessories	57
Reusable Hardware	60
Regulator Mounts	62
iRECEPTAL® Surgical Waste Fluid Management System	63

EZE-VAC® Suction Canister System

EZE-VAC suction canisters are a single-use disposable, designed for convenience by removing the need for a secondary waste liner. The canisters include VAC-GARD, a non-mechanical shut-off valve that supports best practice by eliminating 99.9% of pseudomonas and preventing backflow of aerosols and fluids into the vacuum system.

- Easy-to-read volume marks
- Reflux protector on patient port reduces contamination exposure risk
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
43203-01	1200 mL Canister, Lid and 2 Elbows	32/Case
43206-01	2000 mL Canister, Lid and 2 Elbows	24/Case
43213-01	3300 mL Canister, Lid and 2 Elbows	12/Case

EZE-VAC® 2.0 Suction Canister System

The new EZE-VAC 2.0 suction canister system is very economical and compatible with most hardware in the market. This canister system also includes VAC-GARD® non-mechanical shut-off valve and an anti-foaming emulsion that prevents froth formation and aerosolization, reducing the risk of contaminant exposure while maximizing fluid capacity.

- Easy-to-read volume marks
- VAC-GARD eliminates 99.9% of pseudomonas and protects the vacuum system from fluids
- Compatible with most hardware in the market
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



ITEM CODE	DESCRIPTION	QTY/UNIT
43210-01	800 mL EZE-VAC 2.0 Canister and Lid	32/Case
43202-01	1200 mL EZE-VAC 2.0 Canister and Lid	32/Case
43202-820821	1200 mL EZE-VAC 2.0 Canister Procedure Kit with Non-sterile Suction Tube 3/16 in. ID × 6 ft. (5 mm × 1.8 m) + 3/16 in ID × 18 in. (5 mm × 46 cm)	30/Case

RECEPTAL Liners

RECEPTAL Disposable Liners are designed with permanently attached lids and a multifunctional lid tube which help decrease the risk of fluid exposure. RECEPTAL Liners are to be used with the same size RECEPTAL Outer Shells.

- Anti-foaming emulsion helps prevent froth formation and aerosolization, reducing the risk of contaminant exposure
- 100% leak tested seamless liners minimize risk of contamination
- Non-sterile, BPA-free, not manufactured with natural rubber latex



SIZE (mL)	REF	DESCRIPTION	QTY/UNIT
1000	BULK PACKAGED		
	43010-05	Liner with VAC-GARD® and Pour Spout	50/Case
	43010S ©	Liner with VAC-GARD and 30g Solidifier	50/Case
1500	BULK PACKAGED		
	43025-01	Liner with Pour Spout	50/Case
	43044-01 ©	Liner with VAC-GARD and Pour Spout	50/Case
	43015S ©	Liner with VAC-GARD and 45g Solidifier	50/Case
	INDIVIDUALLY WRAPPED		
	0G903-01 ©	Liner with Pour Spout	50/Case
	0G904-01 ©	Liner with VAC-GARD	50/Case
0L654-01 ©	Liner with VAC-GARD and Pour Spout	50/Case	
2000	BULK PACKAGED		
	43027-01	Liner with Pour Spout,	50/Case
	43041-01	Liner with Ortho Port	50/Case
	43042-01	Liner with VAC-GARD	50/Case
	43044-05 ©	Liner with VAC-GARD and Pour Spout	50/Case
	43020S ©	Liner with VAC-GARD and 60g Solidifier	50/Case
	INDIVIDUALLY WRAPPED		
	0G901-01 ©	Liner	50/Case
	0G902-01 ©	Liner with VAC-GARD	50/Case
	0P508-01 ©	Liner with VAC-GARD and Pour Spout	50/Case
	0G905-01 ©	Liner with Ortho Port	50/Case
0P507-01 ©	Liner with Pour Spout	50/Case	
3000	BULK PACKAGED		
	0L656-01*	Liner with Dual Ports (Ortho + Standard) and Pour Spout	50/Case
	0L655-01*	Liner with VAC-GARD, Dual Ports (Ortho + Standard) and Pour Spout	50/Case
	43030S* ©	Liner with VAC-GARD, Dual Ports (Ortho + Standard), Pour Spout and 90g Solidifier	50/Case

* Dual port allows for dual suctioning during procedures.
Standard Port = 5/16" (7.9mm) OD, 1/4" (6.4mm) ID; Ortho Port = 3/8" (9.6mm) OD, 5/16" (7.8mm) ID

RECEPTAL® Twist Liners

The RECEPTAL Twist liner is a patented all-in-one waste fluid collection and solidifier solution available in 1.5 and 2 L sizes, in both VAC-GARD® and non-VAC-GARD options.

After a procedure, simply twist the integrated and preloaded Solidifier Dispenser Cap at the top of the lid clockwise 180° to easily and safely release solidifier into the liner. The Twist eliminates the need to open ports, creating a closed system where waste fluid is safe, sealed, solid and ready for disposal.

- **Premeasured, preloaded and contact-free Solidifier Dispenser Cap**
- **Integrated solidifier eliminates the need for open ports**
- **All-in-one solidifier and liner solution reduces costs associated with extra inventory**
- **100% leak tested seamless liners minimize risk of contamination**
- **Non-sterile, BPA-free, not manufactured with natural rubber latex**



43121S



43120S shown above
in canister 43445-01
(canister not included)

Premeasured, preloaded and contact-free solidifier bottle

Anti-foaming emulsion helps prevent froth formation and aerosolization, reducing the risk of contaminant exposure

Designed with permanently attached lids and a multi-functional lid tube which helps decrease the risk of fluid exposure

VAC-GARD® option, a non-mechanical shut-off valve that supports best practice by eliminating 99.9% of pseudomonas and preventing backflow of aerosols and fluids into the vacuum system

REF	DESCRIPTION	QTY/UNIT
43110S	1000 mL Liner with VAC-GARD and Pre-attached Solidifier Dispenser Cap	40/Case
43115S ©	1500 mL Liner with VAC-GARD and Pre-attached Solidifier Dispenser Cap	40/Case
43116S ©	1500 mL Liner with Pre-attached Solidifier Dispenser Cap	40/Case
43120S ©	2000 mL Liner with VAC-GARD and Pre-attached Solidifier Dispenser Cap	40/Case
43121S ©	2000 mL Liner with Pre-attached Solidifier Dispenser Cap	40/Case

RECEPTAL® Outer Shells

RECEPTAL Outer Shells are designed with durability and efficiency in mind. The reusable canister is to be used with the same size RECEPTAL Disposable Liner.

- Easy to read volume marks
- Large patient label strip for convenient documentation
- Non-sterile, DEHP-free, PVC-free, not manufactured with natural rubber latex




Color-coded by size for easy selection.

REF	DESCRIPTION	QTY/UNIT
43449-05 ©	1000 mL Outer Shell	10/Case
43423-01 ©	1500 mL Outer Shell	10/Case
43445-01 ©	2000 mL Outer Shell	10/Case
43503-01 ©	3000 mL Outer Shell	8/Case



The canister ring seal firmly secures the liner lid to the canister, maximizing vacuum power to provide the best clinical outcomes.



Large, easy to read volume marks in 25 mL increments help with patient documentation.



Canister "T" designates vacuum source connection to simplify canister set up.

ACCU-MEASURE® Critical Measurement Devices

RECEPTAL
ACCU-MEASURE
Device



Canister
attachment
bracket

ACCU-MEASURE critical measurement canisters enable the user to make contact-free critical measurements of smaller quantities of aspirate (i.e. during pediatric procedures).

The RECEPTAL® ACCU-MEASURE is a disposable, tapered 250 mL critical measurement canister calibrated to scale. Graduations in 5 mL increments from 0 to 50 mL, 10 mL increments from 50 mL to 150 mL and 20 mL increments from 150 mL to 250 mL.

- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
43306-01	250 mL Critical Measurement Canister	Tandem Tube	12/Case
43305-10	Canister Attachment Bracket for 250 mL Canister	Compatible with 1000 mL, 1500 mL and 2000 mL RECEPTAL Canisters	10/Case

RECEPTAL® Solidifier Bottles and Dissolvable Pouches

RECEPTAL Solidifier Bottles and Dissolvable Pouches effectively solidify potentially serious infectious waste within the suction canister. It penetrates rapidly throughout the canister/liner to solidify blood and body fluids and gels. It retains solidity to help prevent spills and leaks for safe and easy bagging and transport to disposal.

- Multiple sizes to match any need
- Dissolvable pouch option eliminates plastic bottle waste
- No shaking required
- Chlorine-free for non-toxic incineration
- Not made with natural rubber latex



4320SP

RECEPTAL Solidifier Pouches

REF	DESCRIPTION	QTY/UNIT
4305SP	500 Solidifier Pouch, solidifies 500 mL	96/Case
4312SP	1200 Solidifier Pouch, solidifies 1200 mL	96/Case
4315SP	1500 Solidifier Pouch, solidifies 1500 mL	96/Case
4320SP	2000 Solidifier Pouch, solidifies 2000 mL	96/Case
4330SP	3000 Solidifier Pouch, solidifies 3000 mL	96/Case



4315SB

4320SB

RECEPTAL Solidifier Bottles

REF	DESCRIPTION	QTY/UNIT
4305SB	500 Solidifier Bottle, solidifies 500 mL	96/Case
4312SB	1200 Solidifier Bottle, solidifies 1200 mL	96/Case
4315SB	1500 Solidifier Bottle, solidifies 1500 mL	96/Case
4320SB	2000 Solidifier Bottle, solidifies 2000 mL	96/Case
4330SB	3000 Solidifier Bottle, solidifies 3000 mL	96/Case
43160SB	16000 Solidifier Bottle, solidifies up to 16,000 mL	20/Case

AMSure® Suction Tubing with Connectors

AMSure Suction Tubing, with molded female connectors and male adapters, provide secure connections to support uninterrupted procedures. Flexible, low-coil memory allows for easy draping.

- Sterile and non-sterile options, BPA-free, not manufactured with natural rubber latex




AS825

DIAMETER (ID)	LENGTH	STERILE	NON-STERILE
3/16 in. (5 mm) Green Connectors	18 in. (46 cm)	AS820 (50/Case)	AS820B (100/Case)
	6 ft. (1.8 m)	AS821 (50/Case)	AS821B (50/Case)
	10 ft. (3.0 m)	AS822 (50/Case)	AS822B (50/Case)
	12 ft. (3.7 m)	AS823 (20/Case)	AS823B (50/Case)
	20 ft. (6.1 m)	AS824 (20/Case)	AS824B (50/Case)
1/4 in. (6 mm) Light Blue Connectors	6 ft. (1.8 m)	AS825 (50/Case)	AS825B (50/Case)
	10 ft. (3.0 m)	AS826 (50/Case)	AS826B (50/Case)
	10 ft. (6 mm × 3.0 m) with Pinch Clamps	AS826C (50/Case)	—
	12 ft. (3.7 m)	AS827 (20/Case)	AS827B (50/Case)
	20 ft. (6.1 m)	AS828 (20/Case)	AS828B (50/Case)
	100 ft. (30.0 m) with Connectors every 6 ft. (1.8 m)	—	AS829B-01 (1Roll/Case)
9/32 in. (7 mm) Dark Blue Connectors	6 ft. (1.8 m)	AS825M (50/Case)	AS825BM (50/Case)
	10 ft. (3.0 m)	AS826M (50/Case)	AS826BM (50/Case)
	12 ft. (3.7 m)	AS827M (20/Case)	AS827BM (50/Case)
	20 ft. (6.1 m)	AS828M (20/Case)	AS828BM (50/Case)

Suction Yankauers and Kits

AMSure® Suction Yankauers are designed with a slip-resistant, lightweight handle that helps provide optimal control and maneuverability.

- **Translucent pathway allows for fluid visualization**
- **Smooth tips are designed with side eyes to minimize the chance of tissue grab and damage**
- **Sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex**



Yankauers

DESCRIPTION			STERILE	NON-STERILE
Rigid	Regular Beveled Tip	Non-Vented	AS830 © 50/Case	AS830B 500/Case
Rigid	Regular Beveled Tip	Vented	AS831 © 50/Case	AS831B 500/Case
Rigid	Bulbous Tip	Non-Vented	AS832 © 50/Case	AS832B 500/Case
Rigid	Bulbous Tip	Vented	AS833 © 50/Case	AS833B 500/Case
Flexible	Regular Beveled Tip	Non-Vented	AS834 © 50/Case	AS834B 500/Case
Flexible	Regular Beveled Tip	Vented	AS835 © 50/Case	AS835B 500/Case
Flexible	Bulbous Tip	Non-Vented	AS836 © 50/Case	AS836B 500/Case
Flexible	Bulbous Tip	Vented	AS837 © 50/Case	AS837B 500/Case

Yankauer & Tubing Kits

REF	DESCRIPTION	QTY/UNIT
AS825P ©	Rigid Yankauer, Bulbous Tip, Non-Vented, pre-attached to 1/4 in. (6 mm) × 6 ft. (1.8 m) Suction Tubing, Sterile	50/Case
AS827P* ©	Rigid Yankauer, Non-Vented with Bulbous Tip, pre-attached to 1/4 in. (6 mm) × 12 ft. (3.7 m) Suction Tubing, Sterile	20/Case


*Made to order

SteriCadi™ Suction Device Holster System

SteriCadi is a sterile disposable holder to store suction devices. The intuitive design is easy to set up and use to help support better infection control outcomes.


- Suction devices can be isolated to help reduce cross-contamination and improve safety measures
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

ITEM CODE	DESCRIPTION	QTY/UNIT
SIH0001	SteriCadi Suction Device Holder, Sterile	30/Case



Flex-Field Angle is designed to house a variety of suction devices to help support storage standardization

Suspended tip stays dry



Holder clamp attaches easily to bed rail

Tubing clasp provides a secure attachment

Suction Catheters and Trays

Suction Catheters

AMSure® Suction Catheters improve patient care by safely removing respiratory secretions from the upper airway. The suction catheters are available in both straight and coiled packs in various French sizes.

- **Opposing staggered eyelets reduce the potential for blocked suction**
- **Whistle tip port enhances controlled suction**
- **Sterile, BPA-free, not manufactured with natural rubber latex**

AS361



REF	DESCRIPTION	FEATURES	QTY/UNIT
AS361 ©	6 Fr, 17 ½ in. (44.5 cm)	Straight Pack	50/Case
AS361C ©	6 Fr, 17 ½ in. (44.5 cm)	Coiled Pack	50/Case
AS362 ©	8 Fr, 17 ½ in. (44.5 cm)	Straight Pack	50/Case
AS362C ©	8 Fr, 17 ½ in. (44.5 cm)	Coiled Pack	50/Case
AS363 ©	10 Fr, 23 in. (58.4 cm)	Straight Pack	50/Case
AS363C ©	10 Fr, 23 in. (58.4 cm)	Coiled Pack	50/Case
AS364 ©	12 Fr, 23 in. (58.4 cm)	Straight Pack	50/Case
AS364C ©	12 Fr, 23 in. (58.4 cm)	Coiled Pack	50/Case
AS365 ©	14 Fr, 23 in. (58.4 cm)	Straight Pack	50/Case
AS365C ©	14 Fr, 23 in. (58.4 cm)	Coiled Pack	50/Case
AS366 ©	16 Fr, 23 in. (58.4 cm)	Straight Pack	50/Case
AS366C ©	16 Fr, 23 in. (58.4 cm)	Coiled Pack	50/Case
AS367 ©	18 Fr, 23 in. (58.4 cm)	Straight Pack	50/Case
AS367C ©	18 Fr, 23 in. (58.4 cm)	Coiled Pack	50/Case

Mini Suction Catheter Trays

The AMSure Mini Suction Catheter Trays' space-saving design offers everything needed for effective procedures. The mini tray serves as a basin to hold sterile water or saline.

- **Opposing staggered eyelets reduce the potential for blocked suction**
- **DeLee Tip's smooth and rounded tip minimizes trauma for pediatric patients**
- **Whistle tip enhances controlled suction**
- **Sterile, BPA-free, not manufactured with natural rubber latex**



SCT10

REF	DESCRIPTION	CATHETER SIZE	FEATURES	QTY/UNIT
SCT08	Tray, Vinyl Gloves and 17 ½ in. Suction Catheter	8 Fr	DeLee Tip	100/Case
SCT10	Tray, Vinyl Gloves and 23 in. Suction Catheter	10 Fr	DeLee Tip	100/Case
SCT12	Tray, Vinyl Gloves and 23 in. Suction Catheter	12 Fr	Whistle Tip	100/Case
SCT14	Tray, Vinyl Gloves and 23 in. Suction Catheter	14 Fr	Whistle Tip	100/Case
SCT16	Tray, Vinyl Gloves and 23 in. Suction Catheter	16 Fr	Whistle Tip	100/Case

Suction Catheter Kits

AMSure® Suction Catheter Kits provide hassle-free and convenient extraction of mucus or saliva from the upper respiratory airway.

- Uniquely designed two-eyed whistle tip offers safe and effective suctioning
- Sterile, BPA-free, not manufactured with natural rubber latex



AS371

REF	DESCRIPTION	CATHETER SIZE	FEATURES	QTY/UNIT
AS371 ©	17 ½ in. (44.5 cm) Whistle Tip Catheter	6 Fr	1 Vinyl Glove	50/Case
AS381 ©	17 ½ in. (44.5 cm) Whistle Tip Catheter	6 Fr	2 Vinyl Gloves	50/Case
AS372 ©	17 ½ in. (44.5 cm) Whistle Tip Catheter	8 Fr	1 Vinyl Glove	50/Case
AS382 ©	17 ½ in. (44.5 cm) Whistle Tip Catheter	8 Fr	2 Vinyl Gloves	50/Case
AS373 ©	23 in. (58.4 cm) Whistle Tip Catheter	10 Fr	1 Vinyl Glove	50/Case
AS383 ©	23 in. (58.4 cm) Whistle Tip Catheter	10 Fr	2 Vinyl Gloves	50/Case
AS374 ©	23 in. (58.4 cm) Whistle Tip Catheter	12 Fr	1 Vinyl Glove	50/Case
AS384 ©	23 in. (58.4 cm) Whistle Tip Catheter	12 Fr	2 Vinyl Gloves	50/Case
AS375 ©	23 in. (58.4 cm) Whistle Tip Catheter	14 Fr	1 Vinyl Glove	50/Case
AS385 ©	23 in. (58.4 cm) Whistle Tip Catheter	14 Fr	2 Vinyl Gloves	50/Case
AS376 ©	23 in. (58.4 cm) Whistle Tip Catheter	16 Fr	1 Vinyl Glove	50/Case
AS386 ©	23 in. (58.4 cm) Whistle Tip Catheter	16 Fr	2 Vinyl Gloves	50/Case
AS387 ©	23 in. (58.4 cm) Whistle Tip Catheter	18 Fr	2 Vinyl Gloves	50/Case

Accessories



43303-01

RECEPTAL® Specimen Traps

RECEPTAL Specimen Traps safely collect and retain specimens for analysis. The traps are compatible with RECEPTAL and EZE-VAC® lids with a pour spout.

- Pour spout cover provides a tight seal to ensure proper suction
- Specimen traps with stiffener offer simplified insertion into the collection port
- Non-sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES		QTY/UNIT
		Suction Port Size	Stiffener	
43303-01 ©	4 in. (10.16 cm), Pour Spout Cover and Tethered Cap	3/8 in. (0.95 cm)		24/Case
43303-02 ©	4 in. (10.16 cm), Pour Spout Cover and Tethered Cap	3/8 in. (0.95 cm)	•	24/Case
43303-05 ©	8 in. (20.32 cm), Pour Spout Cover and Tethered Cap	1/2 in. (1.27 cm)	•	24/Case



External VAC-GARD® Shut-Off Valve

VAC-GARD external non-mechanical shut-off valve helps support best practice by eliminating 99.9% of pseudomonas and preventing backflow of aerosols and fluids into the vacuum system.

- Non-sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
43309-01	External Shut-Off Valve, ¼ in. (0.64 cm) Male x Female Elbow	10/Case



RECEPTAL® Accu-Measure® Device Bracket

This bracket is used to secure 250 mL canister (43306-01) to a RECEPTAL Reusable Canister. Compatible with 1000 mL, 1500 mL and 2000 mL RECEPTAL Canisters.

REF	DESCRIPTION	QTY/UNIT
43305-10	Canister Attachment Bracket for 250 mL Canister	10/Case

Suction Tubing

For permanent Outer Shell-to-vacuum source, or Outer Shell-to-Outer Shell connections.

REF	DESCRIPTION	QTY/UNIT
43455-01	100 ft. (30.48 m) PVC Non-Sterile Suction Tubing 9/32 in. (7 mm) I.D.	1/Case
43457-25	24 in. (0.61 m) PVC Non-Sterile Suction Tubing 9/32 in. (7 mm) I.D.	25/Case



EZE-VAC® Tandem Tubing

EZE-VAC Tandem disposable tubing allows two or more EZE-VAC canisters to work together to meet procedural requirements.

- Non-sterile, BPA-free, not manufactured with natural rubber latex


REF	DESCRIPTION	QTY/UNIT
43317-25	18 in. (0.45 m) Tandem Tubing	25/Case



RECEPTAL Outer Shell Connection Tubing

RECEPTAL Outer Shell connection disposable tubing allows to easily connect the vacuum "T" of the RECEPTAL canisters between them for easy tandem connection.


REF	DESCRIPTION	QTY/UNIT
CTM50B ©	1/4 in. x 20 in. (6.35 mm x 500 mm), Female Connector	100/Case
CTM20B ©	1/4 in. x 8 in. (6.35 mm x 200 mm), Female Connector	300/Case



Tubing Manifold

For connection of 2, 3 or 4 canisters.

REF	DESCRIPTION	QTY/UNIT
TMF001	Tubing Manifold	10/Case



Accessories

RECEPTAL® Elbow Adapters and Taper

Use with Patient Port.



REF	DESCRIPTION	QTY/UNIT
43318-01 ©	1/4 in. (0.64 cm) Male x Female Elbow (use with 3/8 in. (0.95 cm) Orthopedic Port)	100/Case
43316-01 ©	1/4 in. (0.64 cm) Male x Female Elbow	100/Case
43494-25	1/4 in. (0.64 cm) Male x Female Elbow	25/Case
43473-25	Double 1/4 in. (0.64 cm) Male Taper (use with 3/8 in. (0.95 cm) Orthopedic Port)	25/Case

RECEPTAL Outer Shell Parts

Attach to Canister for vacuum suction.



REF	DESCRIPTION	QTY/UNIT
44300-05	T, Double Male 1/4 in. (0.64 cm) with O-Ring and Screw	5/Case
43448-25	Taper, Male 1/4 in. (0.64 cm) x Female (use with Canister Vacuum "T")	25/Case

5-in-1 Straight Connector

AMSure® 5-in-1 Connector is a durable connector used to connect same size tubing or use as a reducer and fits five (3/16 in., 1/4 in., 5/16 in., 3/8 in., 9/16 in.) sizes of tubing.

- **DEHP-free, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	QTY/UNIT
AS807 ©	5-in-1 Connector, Sterile	100/Case
AS807B	5-in-1 Connector, Non-sterile	500/Case

6-in-1 Universal Y Connector

Lightweight, durable polypropylene plastic connectors adapt for almost any use allowing for multiple pieces of tubing to be joined together. Tapered legs fit six sizes of tubing, from 6 mm to 15 mm (1/4 in. to 9/16 in.).

- **DEHP-free, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	QTY/UNIT
UYC001S	6-in-1 Universal Y Connector, Sterile	100/Case
UYC001	6-in-1 Universal Y Connector, Non-sterile	500/Case



Reusable Hardware


Roll Stands

REF	DESCRIPTION	QTY/UNIT
43474-01	Large head, stainless steel Adjustable Roll Stand, conductive casters, holds up to 12 RECEPTAL® or 8 EZE-VAC® Canisters, use with 1500 mL, 2000 mL and 3000 mL canisters	1/Case
43475-01	Small head, stainless steel Adjustable Roll Stand, conductive casters, holds multiple RECEPTAL and EZE-VAC Canisters	1/Case
43476-P2	Built-in Manifold with ON-OFF Valve Adjustable Roll Stand with 4 polycarbonate canister support, use with 1500 mL and 2000 mL canisters	1/Case
43476-P3	Built-in Manifold with ON-OFF Valve Adjustable Roll Stand with 4 polycarbonate canister support, use with 3000 mL canisters	1/Case



Single and In-Line Floor Stands

REF	DESCRIPTION	QTY/UNIT
43431-01	Single Canister Floor Stand with tip-resistant base, use with 1500 mL and 2000 mL RECEPTAL Canisters	1/Case
43432-01	Dual Floor Stand, stainless steel with rollers and conductive casters, use with 1500 mL and 2000 mL RECEPTAL Canisters	1/Case
43447-01	In-Line Quad Floor Stand, stainless steel with rollers and conductive casters, use with 1500 mL and 2000 mL RECEPTAL Canisters	1/Case



Reusable Hardware

RECEPTAL® Liner System Tandem Canister Supports



REF	DESCRIPTION	QTY/UNIT
43429-01*	For 1500 mL and 2000 mL canisters	1/Case
43439-01	For 3000 mL canisters	1/Case

RECEPTAL Liner System Single Canister Supports



REF	DESCRIPTION	QTY/UNIT
43425-05*	For 1000 mL canister, metal	5/Case
43424-10	For 1000 mL canister support, with mounting plate, metal	5/Case
43426-05	For 1500 mL and 2000 mL canisters, metal	5/Case
43438-10	For 3000 mL canisters, with mounting plate for vacuum regulator, metal	10/Case
43428-05*	For 1500 mL and 2000 mL canisters, with mounting plate for vacuum regulator, metal	5/Case
43436-05	For 3000 mL canisters, metal	5/Case
43426P	For 1000 ml, 1500 mL, and 2000 mL canisters, plastic	20/Case
43436P	For 3000 mL canisters, plastic	10/Case
43424-V	For 1000/1500 mL canister support ring with valve	10/Case
43428-V	For 2000 mL canister support ring with valve	10/Case
43438-V	For 3000 mL canister support ring with valve	10/Case

EZE-VAC® Rigid System Single Canister Supports



REF	DESCRIPTION	QTY/UNIT
43467-01	For EZE-VAC 2.0 800 mL and 1200 mL canisters, metal, small: ID=4.4 in. (11.7 cm), wall to outer edge = 8.5 in. (21.6 cm)	12/Case
43469-12*	For EZE-VAC 1200 mL, 2000 mL and 3300 mL canisters, metal, large: ID=5.6 in. (14.2 cm), wall to outer edge = 9.4 in. (23.9 cm)	12/Case
43469-V	For EZE-VAC 1200 mL, 2000 mL, and 3300 mL with valve	10/Case

*Shown on Wall Bracket (REF 43444-05)

Canister Brackets

REF	DESCRIPTION	QTY/UNIT
43444-05	Wall channel bracket for mounting canister, stainless steel	5/Case
43444P	Wall channel bracket for mounting canister, polycarbonate	50/Case
44303-10	Bracket Clamp Adapter for adapting Wall Bracket (43444-05) to fit stands, carts or vertical poles, stainless steel	10/Case
44304-05	Canister Support Bracket to attach Canister directly to a vacuum regulator, use with 1000 mL, 1200 mL and 1500 mL canisters	5/Case
RCB001	Metal Rail Bracket with Universal Holder Clamp	30/Case



Regulator Mounts

REF	DESCRIPTION	QTY/UNIT
43408-05	Use to attach regulator alone or to other supports. 1/8 in. (0.32 cm) male pipe thread × 1/4 in. (0.64 cm) male taper with 1-1/8 in. (2.86 cm) threaded shank and 2 locknuts	5/Case
43422-05	Use when regulator is wall mounted to vacuum and a simple taper is needed to attach tubing to regulator. 1/8 in. (0.32 cm) male pipe thread × 1/4 in. (0.64 cm) male taper	5/Case
43401-05	Use to attach regulator alone or to other supports. 1/4 in. (0.64 cm) taper with 1-1/4 in. (3.18 cm) shank	5/Case
43400-05	Use when regulator is wall mounted to vacuum and a simple taper is needed to attach tubing to regulator. 1/4 in. (0.64 cm), male taper	5/Case
OFV001	Metal on-off valve, 1/8 in. NPT	25/Case



Note: Amsino does not sell Regulators

iRECEPTAL® Surgical Waste Fluid Management System

The iReceptal Digital Surgical Suction System is an all-in-one, closed suction, container, and disposal system. This advanced, easy-to-use surgical suction solution is completely self-contained to provide for independent collection and disposal of waste fluids, along with a docking station for automatic and closed emptying of the rover. iReceptal is currently available in EMEA and China markets. iReceptal is not currently cleared by the US FDA.

- Constantly fully closed system, from collection to disposal
- Staff is never exposed to potentially infectious fluid waste
- Helps optimize surgical turnaround times
- Manifold provides multiple suction lines at the same time with backflow protection
- Separates fluid from tissue and bone fragments
- Waste fluid discharges directly into facility's drain system
- High volume capacity of 20L
- Volume gradation indications
- Easy to maneuver with large wheels
- Lowers disposal costs by eliminating full suction canisters in the waste stream
- Automatically cleans when docked in under 5 minutes

The system is comprised of:

- Mobile Rover
- Docking Station
- Disposable Suction Manifold
- Single-Use Smoke Filter
- Cleaning Detergent




Mobile Rover



Docking Station

REF	DESCRIPTION	QTY/UNIT
iRR100-01 ©	iRECEPTAL 1.0 Rover	1/Each
iRR150-01 ©	iRECEPTAL 1.5 Rover with 7cm (2.7") LCD Screen	1/Each
iRD102-01 ©	iRECEPTAL Docking Station (220VAC)	1/Each
iRD101-01 ©	iRECEPTAL Docking Station (110VAC)	1/Each
iRM100-01 ©	4-Port Manifold (disposable)	1/Each
iRC001-01 ©	Detergent for Docking Station	1/Each
iRC002-01 ©	Detergent for Rover	1/Each



Manifold Storage

- 4 ports for suction tubing with backflow control



Intelligent Control Panel

- Fluid read-out (Version 1.5 only)

Completely Closed System

- Avoid direct touch of waste fluid

Electromagnetic Couplers

- Allow for easy docking




20L Volume

- Suitable for large amounts of fluid

Direct to Drain Systems

- Hands-free disposal and self-cleaning





Enhancing Oxygen Delivery

AMSure Respiratory Therapy Systems are designed to deliver optimal FiO2 levels as well as keep upper airways unobstructed and clear to help improve quality of care. The product portfolio includes:

Respiratory Masks and Nasal Cannulas.....	66
Nebulizer Kits	69
Accessories.....	70
Aerosol Drainage Bag	70
Heat and Moisture Exchangers	71
Tracheostomy Clean and Care Trays	71
Sterile Water	72
Suction Catheters, Trays and Kits	73
AMSure Suction Tubing with connectors	75
Suction Yankauers and Kits	76

Respiratory Masks and Nasal Cannulas



AMSure® ROLLINS7™

AMSure Rollins7 is a replacement mask for virtually every type of oxygen mask on the market with its ability to provide flow ranges from 5-15 liters per minute. Rollins7 is designed with the ability to administer aerosol and nebulizer medication without interrupting oxygen therapy helping improve patient safety, patient comfort and clinical efficiency.

- Medication delivery without sacrificing FiO2. Efficient FiO2 escalation and de-escalation management
- Oxygen tubing provides versatility by connecting to flow meters with or without a Christmas tree nipple nut adapter
- Access for hydration, oral suction or oral medicine
- Non-sterile, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
OM001 ©	Mask, Adapter, 7 ft. Star-shaped Lumen Oxygen Tubing with Label	50/Case

Oxygen Mask

AMSure Oxygen Mask provides medium-flow oxygen delivery. The oxygen mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
AS74010 ©	Adult Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	50/Case

Respiratory Masks and Nasal Cannulas

High Flow Nasal Oxygen Cannula

Amsino's AMSure® high flow oxygen nasal cannula is a patient interface for delivery of humidified respiratory gases.

- Deliver flows of up to 60 lpm
- Rotatable silicone clip for optimizing patient comfort and fit
- Flexible EVA tube with detachable and slide-in pliable nasal prongs for connection to device on either side of patient
- Universal connector to be applied with most high flow nasal cannula devices
- Single patient use, not manufactured with natural rubber latex



AS75280

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS75280 ©	High Flow Nasal Oxygen Cannula	Small, Adult	10/Case
AS75285	High Flow Nasal Oxygen Cannula	Medium, Adult	10/Case
AS75282 ©	High Flow Nasal Oxygen Cannula	Large, Adult	10/Case

Nasal Oxygen Cannula

AMSure Nasal Oxygen Cannulas provide low-to-medium flow oxygen delivery and are designed with soft and pliable prongs which are ideal for patients who require oxygen therapy.

- Over the ear design allows for complete freedom of movement
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex




AS75080

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS75080 ©	Standard Nasal Oxygen Cannula	Flared Prongs, Adult	50/Case
AS75085 ©	Standard Nasal Oxygen Cannula	Non-Flared Prongs, Adult	50/Case

Cushion Nasal Oxygen Cannula

AMSure Cushion Nasal Oxygen Cannulas with soft foam ear covers provide low-to-medium flow oxygen delivery. Soft, pliable prongs are ideal for patients who require long-term oxygen therapy.

- Over the ear design allows for complete freedom of movement
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex




REF	DESCRIPTION	FEATURES	QTY/UNIT
AS75181 ©	Soft Nasal Oxygen Cannula	Curved Prongs with Cushion, Adult	25/Case

Soft Nasal Oxygen Cannula

AMSure® Soft Nasal Oxygen Cannulas provide low-to-medium flow oxygen delivery and are designed with soft and pliable prongs which are ideal for patients who require long-term oxygen therapy.

- Over the ear design allows for complete freedom of movement
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex



AS75180

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS75180 ©	Soft Nasal Oxygen Cannula	Curved Prongs, Adult	50/Case
AS75190 ©	Soft Nasal Oxygen Cannula	Curved Prongs, Pediatric	50/Case
AS75199 ©	Soft Nasal Oxygen Cannula	Curved Prongs, Infant	50/Case

Non-Rebreather Oxygen Mask

AMSure Non-rebreather Oxygen Masks provide high-flow oxygen delivery. The non-rebreather mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- Soft and comfortable elastic strap reduces facial pressure while providing a tight seal
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex



AS75010

REF	DESCRIPTION	QTY/UNIT
AS75010 ©	Adult Non-rebreather Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	50/Case
AS75020 ©	Pediatric Non-rebreather Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	50/Case

Aerosol Mask

AMSure Aerosol Mask for pediatric patients fits any standard 22 mm tubing or nebulizer. The aerosol mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- Soft and comfortable elastic strap reduces facial pressure
- Non-sterile, BPA-free, not manufactured with natural rubber latex



AS75050

REF	DESCRIPTION	QTY/UNIT
AS75040 ©	Adult Aerosol Mask	50/Case
AS75050 ©	Pediatric Aerosol Mask	50/Case

Respiratory Masks and Nasal Cannulas



Tracheostomy Mask

AMSure® Tracheostomy Mask features a soft cuff to improve patient comfort. The 22 mm connector fits standard corrugated tubing for humidified oxygen and nebulizer treatments.

- Soft cuff reduces pressure and irritation to patient skin
- Adapter swivels and rotates 360 degrees, ensuring optimum position for attachments
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS75060 ©	Adult Tracheostomy Mask	50/Case



Nebulizer Kits

AMSure Nebulizer Kits are designed for ease of use and optimal particle size delivery. The nebulizer kits feature a 6 mL medication cup, 7 ft. tubing and either a T-mouthpiece with corrugated tubing or curved mouthpiece.

- Anti-spill medication container lid helps to prevent leakage
- Visible dosing marks designed for ease of use
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex

T-Mouthpiece

- Disposable nebulizers deliver effective treatment to a wide range of patients

Curved Mouthpiece

- Disposable and reusable nebulizers delivers effective treatment to a wide range of patients

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS78010 ©	Nebulizer, 5 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	T-mouthpiece	50/Case
AS78015 ©	Nebulizer, 3 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	Curved mouthpiece	50/Case
AS78020 ©	Nebulizer, 3 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	Aerosol mask	50/Case
AS78025 ©	Nebulizer, 3 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	Curved mouthpiece, Reusable	50/Case

Accessories

Oxygen Tubing

7 ft. Star-shaped Lumen Tubing with Standard 22 mm Connector.

- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	FEATURES	QTY/UNIT
AS76007 ©	7 ft. Oxygen Tubing	Star-shaped Lumen	50/Case

Respiratory Adapters

AMSure® Respiratory Adapters are available in 15 mm I.D. and 22 mm I.D.

- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS730100	Universal Cuff Adapter	22 mm ID Connector	50/Case
AS730200	Multi Adapter	15 mm ID / 22 mm OD Connector	50/Case



AS730100

AS730200

Aerosol Drainage Bag

AMSure Aerosol Drainage System collects moisture accumulating from continuous in-line aerosol therapy. The system includes a 750 mL drainage bag, Y adapter, self-locking plastic chain and accepts 22 mm I.D. corrugated tubing.

- Pre-connected self-locking plastic chain secures the system to a stable object
- Convenient drain port allows for easy one-handed drainage
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS75035	750 mL Aerosol Drainage Bag	50/Case



Heat and Moisture Exchangers



HME001

Amsino is offering Passive Heat and Moisture Exchanger (HME) for Tracheostomy Patients. HME is required when a patient can no longer inspire and expire through the nose and mouth. The nose and mouth provide natural warmth, moisture and filtration for the air we breathe. Having a tracheostomy tube, however, by-passes these mechanisms so humidification must be provided to keep secretions thin and to avoid mucus plugs. In order to keep the environment at an optimal humidity level, Heat Moisture Exchangers are used.

REF	DESCRIPTION	QTY/UNIT
HME001 ©	Heat and Moisture Exchangers, Non-Sterile	50/Case
HME002 ©	Heat and Moisture Exchangers, Sterile	50/Case

Tracheostomy Clean and Care Trays



AS861

AMSure® Tracheostomy Clean and Care Trays feature the necessary components to perform cleaning in and around the stoma. The trays provide clinicians the versatility to perform care based on the patient's individual needs.

- **Sterile, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS860	Cleaning Brush, Vinyl Powder-free Gloves, (2) Cotton Tip Applicators with Plastic Shaft, Twill Tape, Tracheostomy Dressing and Waterproof Drape	Tray, (2) 4" x 4" Gauze, (2) Pipe Cleaners, Pop-up Solution Cup	20/Case
AS861	Cleaning Brush, Vinyl Powder-free Gloves, (2) Cotton Tip Applicators with Plastic Shaft, Twill Tape, Tracheostomy Dressing and Waterproof Drape	2 Compartment Tray, Removable Inner Tray, (4) 4" x 4" Gauze, (4) Pipe Cleaners	20/Case

AMSure® Sterile Water

All products are DEHP-Free, PVC-Free, not manufactured with natural rubber latex, preservative-free, and Non-Pyrogenic.



AS1065

AS1000

Sterile Water for Inhalation (Bottles)

- 500 mL and 750 mL, and 1000 mL bottles Sterile Water, USP
- Used in nebulizers and humidifiers

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
AS1565	AS1565E ©	500 mL Sterile Water for Inhalation, USP	12/Case
AS1765	AS1765E ©	750 mL Sterile Water for Inhalation, USP	12/Case
AS1065	AS1065E ©	1000 mL Sterile Water for Inhalation, USP	12/Case
AS1566	—	500 mL Sterile Water for Inhalation with non-sterile cap	12/Case
AS1066	—	1000 mL Sterile Water for Inhalation with non-sterile cap	12/Case
AS1000	—	1000 mL Sterile Water for Inhalation, Hanging Spikeable Bottle	12/Case

Sterile Water for Inhalation (Vials)

- Available in two sizes: 3 mL and 5 mL, USP
- Easy open twist-off cap is simple to use and allows dispensing by the drop or milliliter
- Used for suctioning and respiratory inhalation



ASR0035

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
ASR0035	ASR0035E ©	3 mL Sterile Water Inhalation, USP	1000/Case
ASR0055	ASR0055E ©	5 mL Sterile Water Inhalation, USP	1000/Case

Sterile Water Prefilled Humidifiers

- 350 mL and 550 mL bottles Sterile Water, USP for inhalation therapy and humidification
- For respiratory therapy, adapters included with each bottle
- Reservoir bottles have built-in oxygen tubing connectors
- Micro-diffuser produces smaller bubbles and greater surface agitation for quiet operation
- Audible alarm alerts clinician of flow restriction or occlusion of oxygen tubing



AS0352

AS0552

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
AS0352	AS0352E ©	350 mL Sterile Water for Inhalation Prefilled Humidifier, USP	20/Case
AS0552	AS0552E ©	550 mL Sterile Water for Inhalation Prefilled Humidifier, USP	12/Case

Suction Catheters, Trays and Kits

Suction Catheters

AMSure® Suction Catheters improve patient care by safely removing respiratory secretions from the upper airway. The suction catheters are available in both straight and coiled packs in various French sizes.

- **Opposing staggered eyelets reduce the potential for blocked suction**
- **Whistle tip port enhances controlled suction**
- **Sterile, BPA-free, not manufactured with natural rubber latex**



REF	DESCRIPTION	FEATURES	QTY/UNIT
AS361 ©	6 Fr, 17 ½ in.	Straight Pack	50/Case
AS361C ©	6 Fr, 17 ½ in.	Coiled Pack	50/Case
AS362 ©	8 Fr, 17 ½ in.	Straight Pack	50/Case
AS362C ©	8 Fr, 17 ½ in.	Coiled Pack	50/Case
AS363 ©	10 Fr, 23 in.	Straight Pack	50/Case
AS363C ©	10 Fr, 23 in.	Coiled Pack	50/Case
AS364 ©	12 Fr, 23 in.	Straight Pack	50/Case
AS364C ©	12 Fr, 23 in.	Coiled Pack	50/Case
AS365 ©	14 Fr, 23 in.	Straight Pack	50/Case
AS365C ©	14 Fr, 23 in.	Coiled Pack	50/Case
AS366 ©	16 Fr, 23 in.	Straight Pack	50/Case
AS366C ©	16 Fr, 23 in.	Coiled Pack	50/Case
AS367 ©	18 Fr, 23 in.	Straight Pack	50/Case
AS367C ©	18 Fr, 23 in.	Coiled Pack	50/Case

Suction Catheter Kits

AMSure Suction Catheters improve patient care by safely removing respiratory secretions from the upper airway. The suction catheters are available in both straight and coiled packs in various French sizes.

- **Uniquely designed two-eyed whistle tip offers safe and effective suctioning**
- **Sterile, BPA-free, not manufactured with natural rubber latex**



REF	DESCRIPTION	FEATURES	QTY/UNIT
AS371 ©	6 Fr, 17 ½ in.	Whistle Tip	50/Case
AS381 ©	6 Fr, 17 ½ in.	Whistle Tip	50/Case
AS372 ©	8 Fr, 17 ½ in.	Whistle Tip	50/Case
AS382 ©	8 Fr, 17 ½ in.	Whistle Tips	50/Case
AS373 ©	10 Fr, 23 in.	Whistle Tip	50/Case
AS383 ©	10 Fr, 23 in.	Whistle Tip	50/Case
AS374 ©	12 Fr, 23 in.	Whistle Tip	50/Case
AS384 ©	12 Fr, 23 in.	Whistle Tip	50/Case
AS375 ©	14 Fr, 23 in.	Whistle Tip	50/Case
AS385 ©	14 Fr, 23 in.	Whistle Tip	50/Case
AS376 ©	16 Fr, 23 in.	Whistle Tip	50/Case
AS386 ©	16 Fr, 23 in.	Whistle Tip	50/Case
AS387 ©	18 Fr, 23 in.	Whistle Tip	50/Case

AS371

Mini Suction Catheter Trays

The AMSure® Mini Suction Catheter Trays' space-saving design offers everything needed for effective procedures. The mini tray serves as a basin to hold sterile water or saline.

- Opposing staggered eyelets reduce the potential for blocked suction
- DeLee Tip's smooth and rounded tip minimizes trauma for pediatric patients
- Whistle tip enhances controlled suction
- Sterile, BPA-free, not manufactured with natural rubber latex




SCT10

REF	DESCRIPTION	FEATURES	QTY/UNIT
SCT08	8 Fr, 17 ½ in. Suction Catheter, Tray and Vinyl Gloves	DeLee Tip	100/Case
SCT10	10 Fr, 23 in. Suction Catheter, Tray and Vinyl Gloves	DeLee Tip	100/Case
SCT12	12 Fr, 23 in. Suction Catheter, Tray and Vinyl Gloves	Whistle Tip	100/Case
SCT14	14 Fr, 23 in. Suction Catheter, Tray and Vinyl Gloves	Whistle Tip	100/Case
SCT16	16 Fr, 23 in. Suction Catheter, Tray and Vinyl Gloves	Whistle Tip	100/Case

Closed Suction Catheter Trays

24-hour Adult Catheter from 10 FR to 14 FR Endotracheal Catheter Kits and Tracheostomy Catheter Kits with MDI Port Adapter.

- Closed system to prevent cross infections in facility
- Clean easily with needle free port
- Avoid interruptions of oxygen supply due to suction



VS365TP124

VS365TC124

REF	DESCRIPTION	FLEX CONNECTOR	QTY/UNIT
Endotracheal Kits with MDI Port Adapter			
VS363RC124 ©	10 Fr, 21.3 in. Endotracheal Catheter, Double Swivel Elbow	15 mm x 22 mm	20/Case
VS364RC124 ©	12 Fr, 21.3 in. Endotracheal Catheter, Double Swivel Elbow	15 mm x 22 mm	20/Case
VS365RC124 ©	14 Fr, 21.3 in. Endotracheal Catheter, Double Swivel Elbow	15 mm x 22 mm	20/Case
VS363TP124 ©	10 Fr, 21.3 in. Endotracheal Catheter, T-Piece	22 mm x 22 mm	20/Case
VS364TP124 ©	12 Fr, 21.3 in. Endotracheal Catheter, T-Piece	22 mm x 22 mm	20/Case
VS365TP124 ©	14 Fr, 21.3 in. Endotracheal Catheter, T-Piece	22 mm x 22 mm	20/Case
Tracheostomy Kits with MDI Port Adapter			
VS363TC124 ©	10 Fr, 12 in. Tracheostomy Catheter, Double Swivel Elbow	15 mm x 22 mm	20/Case
VS364TC124 ©	12 Fr, 12 in. Tracheostomy Catheter, Double Swivel Elbow	15 mm x 22 mm	20/Case
VS365TC124 ©	14 Fr, 12 in. Tracheostomy Catheter, Double Swivel Elbow	15 mm x 22 mm	20/Case
VS363TTP124 ©	10 Fr, 12 in. Tracheostomy Catheter, T-Piece	22 mm x 22 mm	20/Case
VS364TTP124 ©	12 Fr, 12 in. Tracheostomy Catheter, T-Piece	22 mm x 22 mm	20/Case
VS365TTP124 ©	14 Fr, 12 in. Tracheostomy Catheter, T-Piece	22 mm x 22 mm	20/Case

Suction Tubing with Connectors

AMSure® Suction Tubing, with molded female connectors and male adapters, provide secure connections to support uninterrupted procedures. Flexible, low-coil memory allows for easy draping.

- Sterile and non-sterile options, BPA-free, not manufactured with natural rubber latex



AS825

DIAMETER (ID)	LENGTH	STERILE	NON-STERILE
3/16 in. (5 mm) Green Connectors	18 in. (46 cm)	AS820 (50/Case)	AS820B (100/Case)
	6 ft. (1.8 m)	AS821 (50/Case)	AS821B (50/Case)
	10 ft. (3.0 m)	AS822 (50/Case)	AS822B (50/Case)
	12 ft. (3.7 m)	AS823 (20/Case)	AS823B (50/Case)
	20 ft. (6.1 m)	AS824 (20/Case)	AS824B (50/Case)
1/4 in. (6 mm) Light Blue Connectors	6 ft. (1.8 m)	AS825 (50/Case)	AS825B (50/Case)
	10 ft. (3.0 m)	AS826 (50/Case)	AS826B (50/Case)
	10 ft. (6 mm × 3.0 m) with Pinch Clamps	AS826C (50/Case)	—
	12 ft. (3.7 m)	AS827 (20/Case)	AS827B (50/Case)
	20 ft. (6.1 m)	AS828 (20/Case)	AS828B (50/Case)
	100 ft. (30.0 m) with Connectors every 6 ft. (1.8 m)	—	AS829B-01 (1 Roll/Case)
9/32 in. (7 mm) Dark Blue Connectors	6 ft. (1.8 m)	AS825M (50/Case)	AS825BM (50/Case)
	10 ft. (3.0 m)	AS826M (50/Case)	AS826BM (50/Case)
	12 ft. (3.7 m)	AS827M (20/Case)	AS827BM (50/Case)
	20 ft. (6.1 m)	AS828M (20/Case)	AS828BM (50/Case)

Suction Yankauers and Kits

AMSure® Suction Yankauers are designed with a slip-resistant, lightweight handle that helps provide optimal control and maneuverability.

- **Translucent pathway allows for fluid visualization**
- **Smooth tips are designed with side eyes to minimize the chance of tissue grab and damage**
- **Sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex**



AS835

AS832

Yankauers

DESCRIPTION			STERILE	NON-STERILE
Rigid	Regular Beveled Tip	Non-Vented	AS830 © 50/Case	AS830B 500/Case
Rigid	Regular Beveled Tip	Vented	AS831 © 50/Case	AS831B 500/Case
Rigid	Bulbous Tip	Non-Vented	AS832 © 50/Case	AS832B 500/Case
Rigid	Bulbous Tip	Vented	AS833 © 50/Case	AS833B 500/Case
Flexible	Regular Beveled Tip	Non-Vented	AS834 © 50/Case	AS834B 500/Case
Flexible	Regular Beveled Tip	Vented	AS835 © 50/Case	AS835B 500/Case
Flexible	Bulbous Tip	Non-Vented	AS836 © 50/Case	AS836B 500/Case
Flexible	Bulbous Tip	Vented	AS837 © 50/Case	AS837B 500/Case

Pre-Connected Yankauer & Tubing

REF	DESCRIPTION	QTY/UNIT
AS825P ©	Rigid Yankauer, Bulbous Tip, Non-Vented, pre-attached to 1/4 in. (6 mm) × 6 ft. (1.8 m) Suction Tubing, Sterile	50/Case
AS827P* ©	Rigid Yankauer, Non-Vented with Bulbous Tip, pre-attached to 1/4 in. (6 mm) × 12 ft. (3.7 m) Suction Tubing, Sterile	20/Case

*Made to order

AMSure®

SURGICAL AND NURSING SOLUTIONS



Everyday Patient Care Solutions

AMSure Surgical and Nursing Systems offer a variety of frequently used products for patient care. These essentials include:

Suction Tubing	78
Suction Yankauers and Kits	79
Suction Device Holster System	79
Specimen Containers.....	80
Vaginal Speculums.....	80
Disposable Umbilical Cord Clamp.....	80
Ear/Ulcer Syringes	81
Cleansing Enema Sets.....	81
Connectors.....	82

Suction Tubing with Connectors

AMSure® Suction Tubing, with molded female connectors and male adapters, provide secure connections to support uninterrupted procedures. Flexible, low-coil memory allows for easy draping.

- Sterile and non-sterile options, BPA-free, not manufactured with natural rubber latex



AS825

DIAMETER (ID)	LENGTH	STERILE	NON-STERILE
3/16 in. (5 mm) Green Connectors	18 in. (46 cm)	AS820 © (50/Case)	AS820B © (100/Case)
	6 ft. (1.8 m)	AS821 © (50/Case)	AS821B © (50/Case)
	10 ft. (3.0 m)	AS822 © (50/Case)	AS822B © (50/Case)
	12 ft. (3.7 m)	AS823 © (20/Case)	AS823B © (50/Case)
	20 ft. (6.1 m)	AS824 © (20/Case)	AS824B © (50/Case)
1/4 in. (6 mm) Light Blue Connectors	6 ft. (1.8 m)	AS825 © (50/Case)	AS825B © (50/Case)
	10 ft. (3.0 m)	AS826 © (50/Case)	AS826B © (50/Case)
	10 ft. (6 mm × 3.0 m) with Pinch Clamps	AS826C © (50/Case)	—
	12 ft. (3.7 m)	AS827 © (20/Case)	AS827B © (50/Case)
	20 ft. (6.1 m)	AS828 © (20/Case)	AS828B © (50/Case)
	100 ft. (30.0 m) with Connectors every 6 ft. (1.8 m)	—	AS829B-01 © (1Roll/Case)
9/32 in. (7 mm) Dark Blue Connectors	6 ft. (1.8 m)	AS825M © (50/Case)	AS825BM © (50/Case)
	10 ft. (3.0 m)	AS826M © (50/Case)	AS826BM © (50/Case)
	12 ft. (3.7 m)	AS827M © (20/Case)	AS827BM © (50/Case)
	20 ft. (6.1 m)	AS828M © (20/Case)	AS828BM © (50/Case)

Suction Yankauers and Kits

AMSure® Suction Yankauers are designed with a slip-resistant, lightweight handle that helps provide optimal control and maneuverability.

- **Translucent pathway allows for fluid visualization**
- **Smooth tips are designed with side eyes to minimize the chance of tissue grab and damage**
- **Sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex**



AS835

AS832

Yankauers

DESCRIPTION			STERILE	NON-STERILE
Rigid	Regular Beveled Tip	Non-Vented	AS830 © 50/Case	AS830B 500/Case
Rigid	Regular Beveled Tip	Vented	AS831 © 50/Case	AS831B 500/Case
Rigid	Bulbous Tip	Non-Vented	AS832 © 50/Case	AS832B 500/Case
Rigid	Bulbous Tip	Vented	AS833 © 50/Case	AS833B 500/Case
Flexible	Regular Beveled Tip	Non-Vented	AS834 © 50/Case	AS834B 500/Case
Flexible	Regular Beveled Tip	Vented	AS835 © 50/Case	AS835B 500/Case
Flexible	Bulbous Tip	Non-Vented	AS836 © 50/Case	AS836B 500/Case
Flexible	Bulbous Tip	Vented	AS837 © 50/Case	AS837B 500/Case

Pre-Connected Yankauer & Tubing

REF	DESCRIPTION	QTY/UNIT
AS825P ©	Rigid Yankauer, Bulbous Tip, Non-Vented, pre-attached to 1/4 in. (6 mm) × 6 ft. (1.8 m) Suction Tubing, Sterile	50/Case
AS827P* ©	Rigid Yankauer, Non-Vented with Bulbous Tip, pre-attached to 1/4 in. (6 mm) × 12 ft. (3.7 m) Suction Tubing, Sterile	20/Case

*Made to order

Suction Device Holster System

SteriCadi™ is a sterile disposable holder to store suction devices. The intuitive design is easy to set up and use to help support better infection control outcomes.

- **Suction devices can be isolated to help reduce cross-contamination and improve safety measures**
- **Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	QTY/UNIT
SIH0001 ©	SteriCadi Suction Device Holder with Clip	30/Case



Specimen Containers

AMSure® 4 oz., leak-resistant specimen containers are designed for the collection and safe transport of patient samples.

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	FORM FILL SEAL	POLY POUCH	BULK	STERILE INTERIOR	TAMPER EVIDENT LABEL	PATIENT LABELS	QTY/UNIT
AS340 ©		•		•		•	100/Case
AS341 ©			•	•	•		100/Case
AS343			•				500/Case
AS345 ©	•			•			100/Case



AS340

Vaginal Speculums

AMSure Vaginal Speculum helps to increase patient comfort during insertion and may provide greater visualization. The speculums are available in small, medium and large sizes.

- Round-edge design helps to increase patient comfort
- Single-use disposable
- Non-sterile, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS032S ©	Small Vaginal Speculum	100/Case
AS032M ©	Medium Vaginal Speculum	100/Case
AS032L ©	Large Vaginal Speculum	100/Case



AS032M

Disposable Umbilical Cord Clamp

AMSure Umbilical Cord Clamp with safety lock helps provide a secure connection.

- Safety lock is designed to decrease risk of accidental disconnection
- Grooved teeth help provide an effective grip
- Sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
UCC100	Umbilical Cord Clamp, sterile	50/Case





Ear/Ulcer Syringes

AMSure® Ear/Ulcer syringes are ideal for irrigating smaller wounds and alleviating congestion. The syringes are available in several sizes to meet clinical requirements.

- Gentle tip design reduces irritation
- Ribbed bulb provides steady grip control
- Designed to stand on the bulb end to help minimize contamination
- Sterile and non-sterile options, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS00501	1 oz. syringe	Non-sterile, Bulk	300/Case
AS00502	2 oz. syringe	Non-sterile, Bulk	200/Case
AS00503	3 oz. syringe	Non-sterile, Bulk	150/Case
AS00501S ©	1 oz. syringe	Sterile, Single Pack	50/Case
AS00502S ©	2 oz. syringe	Sterile, Single Pack	50/Case
AS00503S ©	3 oz. syringe	Sterile, Single Pack	50/Case



AS330

Cleansing Enema Sets

AMSure Cleansing Enema Bag and Bucket Sets provide the necessary components to perform a comfortable and safe procedure. The set includes 60 in. (152.4 cm) tubing with pre-lubricated tip, soap packet, waterproof drape and either a 1500 mL bag or 1500 mL bucket.

Bag Set

- One-way bag valve prevents solution spills
- Pre-lubricated tip allows for easy and comfortable insertion
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

Bucket Set

- Graduated markings designed for precision and ease of use
- Built-in handles provide convenient transport
- Pre-lubricated tip allows for easy and comfortable insertion
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



AS333

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS330	60 in. (152.4 cm) Tubing with Pre-Lubricated Tip, Slide Clamp, Soap Packet and Waterproof Drape	1500 mL Bag	50/Case
AS333	60 in. (152.4 cm) Tubing with Pre-Lubricated Tip, Slide Clamp, Soap Packet and Waterproof Drape	1500 mL Bucket	50/Case

Connectors

5-in-1 Straight Connector

AMSure® 5-in-1 Connector is a durable connector used to connect same size tubing or use as a reducer and fits five (3/16 in., 1/4 in., 5/16 in., 3/8 in., 9/16 in.) sizes of tubing.

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS807 [Ⓒ]	5-in-1 Connector	Sterile	100/Case
AS807B	5-in-1 Connector	Non-sterile, Bulk	500/Case



AS807

6-in-1 Universal Y Connector

Lightweight, durable polypropylene plastic connectors adapt for almost any use allowing for multiple pieces of tubing to be joined together. Tapered legs fit six sizes of tubing, from 6 mm to 15 mm (1/4 in. to 9/16 in.).

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
UYC001S	6-in-1 Universal Y Connector, Sterile	100/Case
UYC001	6-in-1 Universal Y Connector, Non-sterile	500/Case



Index

REF	PAGE	REF	PAGE	REF	PAGE	REF	PAGE
000606.....	16	43316-01	59	AS010	4, 39	AS372	56, 73
003201.....	16	43317-25	58	AS011	4, 39	AS373	56, 73
003202.....	16	43318-01	55	AS015	39	AS374	56, 73
0G901-01	48	43400-05	62	AS0159	44	AS375	56, 73
0G902-01	48	43401-05	62	AS016	39	AS376	56, 73
0G903-01	48	43408-05	62	AS032L	80	AS381	56, 73
0G904-01	48	43422-05	62	AS032M	80	AS382	56, 73
0G905-01	48	43423-01	50	AS032S	80	AS383	56, 73
0L654-01	48	43424-10	61	AS0352	42, 72	AS384	56, 73
0L655-01	48	43424-V	61	AS0352E	42, 72	AS385	56, 73
0L656-01	48	43425-05	61	AS0359	44	AS386	56, 73
0P507-01	48	43426-05	61	AS0552	42, 72	AS387	56, 73
0P508-01	48	43426P	61	AS0552E	42, 72	AS398	44
1010001B	9	43428-05	61	AS0559	44	AS398E	44
107201.....	7	43428-V	61	AS1000	42, 72	AS399	44
108301.....	7	43429-01	61	AS1065	42, 72	AS399E	44
108302.....	9	43431-01	60	AS1065E	42, 72	AS409	37
108303.....	7	43432-01	60	AS1066	42, 72	AS40008	29
108304.....	7	43436-05	61	AS115	4, 39	AS40010	29
108305.....	8	43436P	61	AS116	4, 39	AS41012	29
108306.....	7	43438-10	61	AS121	5, 38	AS41012S	30
108308.....	8	43438-V	61	AS126	5, 38	AS41014	29
109602.....	7	43439-01	61	AS127	5, 38	AS41014H	28
1510501.....	10	43444-05	62	AS130	5, 38	AS41014S	30
158305.....	10	43444P	62	AS136	5, 38	AS41016	29
158306.....	10	43445-01	50	AS1565	42, 72	AS41016H	28
158307.....	10	43447-01	60	AS1565E	42, 72	AS41016S	30
158309.....	10	43448-25	59	AS1566	42, 72	AS41018	29
158310.....	9	43449-05	50	AS1765	42, 72	AS41018H	28
159602.....	10	43455-01	58	AS1765E	42, 72	AS41018S	30
43010-05	48	43457-25	58	AS302	26	AS41020	29
43010S.....	48	43467-01	61	AS306N	37	AS41020S	30
43015S.....	48	43469-12	61	AS307N	37	AS41022	29
43020S.....	48	43469-V	61	AS308N	37	AS41022S	30
43025-01	48	43473-25	59	AS309N	37	AS41024	29
43027-01	48	43474-01	60	AS3100	11	AS41024S	30
43030S.....	48	43475-01	60	AS3101	11	AS41026	29
43041-01	48	43476-P2	60	AS3102	11	AS41026S	30
43042-01	48	43476-P3	60	AS3110A	11	AS41028	29
43044-01	48	43494-25	59	AS3110B	11	AS41028S	30
43044-05	48	43503-01	50	AS312	26	AS41030	29
4305SB.....	51	44300-05	59	AS322	26	AS42012	29
4305SP.....	51	44303-10	62	AS32200	26	AS42012S	30
43110S.....	49	44304-05	62	AS326	26	AS42014	29
43115S.....	49	607201.....	13	AS328	27	AS42014H	28
43116S.....	49	608304.....	13	AS32600	26	AS42014S	30
43120S.....	49	608305.....	13	AS327100	27	AS42016H	28
43121S.....	49	608306.....	13	AS330	81	AS42016	29
4312SB.....	51	AA3101	9	AS333	81	AS42016S	30
4312SP.....	51	ABL400	36	AS340	40, 80	AS42018	29
4315SB.....	51	ABL402	36	AS341	40, 80	AS42018H	28
4315SP.....	51	ABL404	36	AS343	40, 80	AS42018S	30
43160SB.....	51	AC3110	7	AS345	40, 80	AS42020	29
43202-01	47	AE0106	17	AS3500	15	AS42020S	30
43202-820821	47	AE0107	17	AS3515	15	AS42022	29
43203-01	46	AE1106	17	AS361	55, 73	AS42024	29
43206-01	46	AE1112	17	AS361C	55, 73	AS42026	29
4320SB.....	51	AE1138	17	AS362	55, 73	AS42028	29
4320SP.....	51	AE3107	16	AS362C	55, 73	AS42030	29
43210-01	47	AE3108	16	AS363	55, 73	AS44008	35
43213-01	46	AFS102	8	AS363C	55, 73	AS44010	35
43303-01	57	AFS107	17	AS364	55, 73	AS44012	35
43303-02	57	AFS208	8	AS364C	55, 73	AS44014	35
43303-05	57	AFS208B	8	AS365	55, 73	AS44016	35
43305-10	51	AS00501	81	AS365C	55, 73	AS44018	35
43305-10	57	AS00501S	81	AS366	55, 73	AS44020	35
43306-01	51	AS00502	81	AS366C	55, 73	AS44022	35
43309-01	57	AS00502S	81	AS367	55, 73	AS44024	35
4330SB.....	51	AS00503	81	AS367C	55, 73	AS44028	35
4330SP.....	51	AS00503S	81	AS371	56, 73	AS730100	70

REF	PAGE	REF	PAGE	REF	PAGE	REF	PAGE
AS730200	70	AS85008	33	ASR0059	44	MR4T07	16
AS74010	66	AS85008K	33	ASR0059E	44	MR4TF04	17
AS75010	68	AS85010	33	ASR0159	44	MR4Y32	16
AS75020	68	AS85010K	33	ASR0159E	44	MRA102-84	9
AS75035	70	AS85012	33	AY0200	18	MRA10E-104	8
AS75040	68	AS85012K	33	AY0206	18	MRA10E-95	7
AS75050	68	AS85014	33	CTM20B	58	MRA15-100	10
AS75060	69	AS85014K	33	CTM50B	58	MRA15B-105	11
AS75080	67	AS85016	33	D38301	14	MRA60B-92	13
AS75085	67	AS85016K	33	D38302	14	MRS60B-106	14
AS75180	68	AS860	71	D38305	14	MRS60B-117	14
AS75181	67	AS860608	34	E-0500M	3	MRSA10-84	7
AS75190	68	AS860610	34	E-1004	3	MRSA103-94	9
AS75199	68	AS860612	34	E-1200	3	MRSA103-94	9
AS75280	67	AS860614	34	E-1203	3	MRSA10E-107	8
AS75282	67	AS860616	34	ENF1004K	3	MRSA15-100	10
AS75285	67	AS861	71	ENF1200K	3	MRSA15-91	10
AS76007	70	AS861006	34	ENF1203K	3	MRSA15E-120	10
AS78010	69	AS861008	34	ENF500K	3	MRSA15E-88	10
AS78015	69	AS861010	34	ENK127	5	MRSA15E-96	11
AS78020	69	AS861608	34	ENK136	5	MRSA60-100	13
AS78025	69	AS861610	34	ENS015	4	MRSA60-84	13
AS807	59, 82	AS861610C	34	ENS016	4	MRSA602-100	13
AS807B	59, 82	AS861612	34	ENS115	4	MRSA60E-88	13
AS820	52, 75, 78	AS861612C	34	ENS115NS	4	MRSB-102	15
AS820B	52, 75, 78	AS861614	34	ENS116	4	MRSB-105	15
AS821	52, 75, 78	AS861614C	34	FP092	19	MRSB-94	15
AS821B	52, 75, 78	AS861616	34	FP092Y1	19	MRSFC-19	17
AS822	52, 75, 78	AS861616C	34	FP100BI2	20	OFV001	62
AS822B	52, 75, 78	AS861618	34	FP100BY2	19	OM001	66
AS823	52, 75, 78	AS861618C	34	FP100FY1	19	PE90B12T	2
AS823B	52, 75, 78	AS870	31	FP106FI2	20	PE90B5T	2
AS824	52, 75, 78	AS87114	31	FP106FY2	19	PE90ST	2
AS824B	52, 75, 78	AS87118	31	FP109BI3	20	PG001	2
AS825	52, 75, 78	AS880	25	FP109BY3	19	RCB001	62
AS825B	52, 75, 78	AS880K	25	FP169BY3	19	S-1000-SI	3
AS825BM	52, 75, 78	AS890	25	HME001	71	SCT08	55, 74
AS825M	52, 75, 78	AS890K	25	HME002	71	SCT10	55, 74
AS825P	53, 76, 79	AS89214	32	iRC001-01	63	SCT12	55, 74
AS826	52, 75, 78	AS89214R	32	iRC002-01	63	SCT14	55, 74
AS826B	52, 75, 78	AS89214RK	32	iRD101-01	63	SCT16	55, 74
AS826BM	52, 75, 78	AS89216RK	32	iRD102-01	63	SIH0001	54, 79
AS826C	52, 75, 78	AS89300	23	iRM100-01	63	TMF001	58
AS826M	52, 75, 78	AS89314	22	iRR100-01	63	UCC100	80
AS827	52, 75, 78	AS89314H	22	iRR150-01	63	UYC001	59, 82
AS827B	52, 75, 78	AS89314S	22	IVF1003	43	UYC001S	59, 82
AS827BM	52, 75, 78	AS89316	22	IVF1005	43	VS363RC124	74
AS827M	52, 75, 78	AS89316H	22	IVF1010	43	VS363TC124	74
AS827P	53, 76, 79	AS89316S	22	IVF1010TM	43	VS363TP124	74
AS828	52, 75, 78	AS89318	22	IVF1203	43	VS363TTP124	74
AS828B	52, 75, 78	AS89318H	22	IVF1203SS	43	VS364RC124	74
AS828BM	52, 75, 78	AS89318S	22	IVF1205	43	VS364TC124	74
AS828M	52, 75, 78	AS89330	23	IVF1205SS	43	VS364TP124	74
AS829B-01	52, 75, 78	AS89800	25	IVF1210	43	VS364TTP124	74
AS830	53, 76, 79	AS89814	24	IVF1210SS	43	VS365RC124	74
AS830B	53, 76, 79	AS89814H	24	IVF2020	43	VS365TC124	74
AS831	53, 76, 79	AS89814S	24	IVF2020SS	43	VS365TP124	74
AS831B	53, 76, 79	AS89816	24	LBE306N	37	VS365TTP124	74
AS832	53, 76, 79	AS89816H	24	LBE307N	37		
AS832B	53, 76, 79	AS89816S	24	LBE308N	37		
AS833	53, 76, 79	AS89818	24	LBE309N	37		
AS833B	53, 76, 79	AS89818H	24	MB4V01	18		
AS834	53, 76, 79	AS89818S	24	MB5S001	18		
AS834B	53, 76, 79	AS89830	25	MR4001	18		
AS835	53, 76, 79	ASR0035	42, 72	MR4107	16		
AS835B	53, 76, 79	ASR0035E	42, 72	MR4207	16		
AS836	53, 76, 79	ASR0039	44	MR4307	17		
AS836B	53, 76, 79	ASR0039E	44	MR4B07	17		
AS837	53, 76, 79	ASR0055	42, 72	MR4P07	18		
AS837B	53, 76, 79	ASR0055E	42, 72	MR4R10	16		



Advancing Healthcare Worldwide[®]

Amsino Medical Group Co Ltd.

708 Corporate Center Drive
Pomona, CA 91768
USA

P: +1 800.632.6746

F: +1 909.626.3888

info@amsino.com

www.amsino.com

All products listed are legally marketed and sold within the US and, where required, cleared by the US FDA. Products listed may not be available or licensed in all geographies. Amsino's manufacturing facilities are certified to ISO 13485:2016 and MDSAP. Please contact our customer service department or our regional office for product availability in your country or region.

Items listed with © are CE marked and legally available within the European Union (EU) and other regions where CE marking may be required.

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMSINO®, ACCU-MEASURE®, AMSafe®, AMSafe-3®, AMShield®, AMSure®, AMSure HG®, EZE-VAC®, iRECEPTAL®, One Mask™, Puggie®, R-Polished®, Receptal®, Safgard™, Sal-Sense™, SteriCadi®, Sure-Lok™, THE-POLE-SYRINGE® and VAC-GARD® are registered trademarks of Amsino Medical Group Co. Ltd. © Copyright 2023, Amsino Medical Group Co. Ltd. All rights reserved.