



SURGICAL PRODUCT  
CATALOGUE 2022

AUSTRALIA &  
NEW ZEALAND



# QUICK LINKS

## SURGICAL GOWNS

<b>AERO CHROME*</b>	page 3
<b>AERO BLUE*</b>	page 4
<b>Surgical Gowns</b>	page 5

## SURGICAL DRAPES & PACKS

<b>Cardiac Catheterisation / Radiology Drapes</b>	page 6
<b>Cardiovascular Drapes</b>	page 8
<b>ENT Drapes &amp; Packs</b>	page 9
<b>Equipment Covers</b>	page 12
<b>General Drapes &amp; Packs:</b>	
Accessories	page 13
Basic Packs	page 14
Basic Drapes	page 15
Breast / Chest Drapes	page 17
Fenestrated Drapes	page 18
Laparoscopy / Laparotomy Drapes	page 19
Minor Procedure Drapes	page 21
Thyroid Drapes	page 22
Universal Packs	page 23
Universal Drapes	page 26
<b>Neurosurgical Drape</b>	page 28
<b>Obstetrics and Gynaecology Drapes &amp; Packs:</b>	
C-Section Drapes	page 29
Laparoscopy Packs	page 30
Laparoscopy Drapes	page 31
Lithotomy Packs	page 32
Lithotomy Drapes	page 33
Under Buttocks Drapes	page 34

<b>Ophthalmic Drapes &amp; Packs</b>	page 35
<b>Orthopaedic Drapes &amp; Packs:</b>	
Ortho Packs	page 36
Extremity Drapes	page 38
Hand Drapes	page 40
Hip Drapes	page 41
Knee Arthroscopy Drapes	page 42
Shoulder Drapes	page 44
U & Bar Drapes	page 46
<b>Urology Drapes &amp; Packs</b>	page 50

## STANDARD PACKS

<b>ENT Pack</b>	page 51
<b>Universal Pack</b>	page 51
<b>Cysto-Gynae Pack</b>	page 52
<b>Caesarean Pack</b>	page 52
<b>Lower Extremity Pack</b>	page 53

<b>CUSTOM PROCEDURE PACKS (CPP)</b>	page 54
-------------------------------------	---------

<b>INDEX</b>	page 55
--------------	---------

# SURGICAL GOWNS



## AERO CHROME\* BREATHABLE PERFORMANCE SURGICAL GOWN

- Soft Knitted V-neck Collar
- Clearly Identifiable Silver Colour
- Soft, Light and Breathable Fabric
- Raglan Sleeves
- Adjustable Hook and Loop Neckline
- Book Fold
- 1 x Overwrap, 90cm x 90cm
- 1 x Blue Huck Towel Included for all Sterile Codes, 51cm x 64cm
- AAMI Level 4 / EN 13795 High Performance

Code	Description	Packaging
<b>44680</b>	AERO CHROME* Breathable Performance Surgical Gown, Sterile, Small	Items/Case = 30
44661	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, Small	Items/Case = 36
<b>44681</b>	AERO CHROME* Breathable Performance Surgical Gown, Sterile, Large	Items/Case = 30
44662	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, Large	Items/Case = 40
<b>44682</b>	AERO CHROME* Breathable Performance Surgical Gown, Sterile, X-Large	Items/Case = 28
44664	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, X-Large	Items/Case = 32
<b>44683</b>	AERO CHROME* Breathable Performance Surgical Gown, Sterile, XX-Large	Items/Case = 26
44666	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, XX-Large	Items/Case = 32
<b>44684</b>	AERO CHROME* Breathable Performance Surgical Gown, Sterile, XXX-Large	Items/Case = 24
<b>44685</b>	AERO CHROME* Breathable Performance Surgical Gown, Sterile, Large, X-Long	Items/Case = 28
44667	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, Large, X-Long	Items/Case = 36
<b>44686</b>	AERO CHROME* Breathable Performance Surgical Gown, Sterile, X-Large, X-Long	Items/Case = 30
44668	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, X-Large, X-Long	Items/Case = 32
<b>44687</b>	AERO CHROME* Breathable Performance Surgical Gown, Sterile, XX-Large, X-Long	Items/Case = 28
44669	AERO CHROME* Breathable Performance Surgical Gown, Non-Sterile, XX-Large, X-Long	Items/Case = 30

All codes in **bold** indicate a sterile product.

NOTE: Please see package inserts for directions for use for each of these products.

# SURGICAL GOWNS




## AERO BLUE\* PERFORMANCE SURGICAL GOWNS

- Patented Breathable Fabric
- Low Linting Fabric
- Raglan Sleeves
- Colour-Coded Neck Binding – Green
- Adjustable Hook and Loop Neckline
- Book Fold
- 1 x Overwrap, 90cm x 90cm
- 1 x Blue Huck Towel Included for all Sterile Codes, 51cm x 64cm
- AAMI Level 3 / EN 13795 High Performance

<b>Code</b>	<b>Description</b>	<b>Packaging</b>
<b>41870</b>	AERO BLUE* Performance Surgical Gown, Sterile, Small	Items/Case = 30
43941	AERO BLUE* Performance Surgical Gown, Non-Sterile, Small	Items/Case = 36
<b>41871</b>	AERO BLUE* Performance Surgical Gown, Sterile, Large	Items/Case = 30
41739	AERO BLUE* Performance Surgical Gown, Non-Sterile, Large	Items/Case = 40
<b>42620</b>	AERO BLUE* Performance Surgical Gown, Sterile, Large, X-Long	Items/Case = 28
41725	AERO BLUE* Performance Surgical Gown, Non-Sterile, Large, X-Long	Items/Case = 36
<b>41872</b>	AERO BLUE* Performance Surgical Gown, Sterile, X-Large	Items/Case = 28
41740	AERO BLUE* Performance Surgical Gown, Non-Sterile, X-Large	Items/Case = 32
<b>42621</b>	AERO BLUE* Performance Surgical Gown, Sterile, X-Large, X-Long	Items/Case = 30
41727	AERO BLUE* Performance Surgical Gown, Non-Sterile, X-Large, X-Long	Items/Case = 32
<b>41873</b>	AERO BLUE* Performance Surgical Gown, Sterile, XX-Large	Items/Case = 26
41741	AERO BLUE* Performance Surgical Gown, Non-Sterile, XX-Large	Items/Case = 32
<b>43940</b>	AERO BLUE* Performance Surgical Gown, Sterile, XX-Large, X-Long	Items/Case = 28
43623	AERO BLUE* Performance Surgical Gown, Non-Sterile, XX-Large, X-Long	Items/Case = 30
<b>41874</b>	AERO BLUE* Performance Surgical Gown, Sterile, XXX-Large	Items/Case = 24

All codes in **bold** indicate a sterile product.  
NOTE: Please see package inserts for directions for use for each of these products.

# SURGICAL GOWNS




## SURGICAL GOWNS

- EVOLUTION\* 4 Fabric
- Low Linting Fabric
- Set-In-Sleeves
- Colour-Coded Neck Binding – Yellow
- Adjustable Hook and Loop Neckline
- Book Fold
- 1 x Overwrap, 90cm x 90cm
- 1 x Blue Huck Towel Included for all Sterile Codes, 51cm x 64cm
- EN 13795 Standard Performance

<b>Code</b>	<b>Description</b>	<b>Packaging</b>
<b>90003</b>	Surgical Gown, Sterile, Small–Medium	Items/Case = 36
<b>90018</b>	Surgical Gown, Sterile, Large	Items/Case = 36
72978	Surgical Gown, Non-Sterile, Large	Items/Case = 640
<b>90019</b>	Surgical Gown, Sterile, Large, X-Long	Items/Case = 34
<b>90048</b>	Surgical Gown, Sterile, X-Large	Items/Case = 32
72988	Surgical Gown, Non-Sterile, X-Large	Items/Case = 640
<b>90049</b>	Surgical Gown, Sterile, X-Large, X-Long	Items/Case = 34
70392	Surgical Gown, Non-Sterile, X-Large, X-Long	Items/Case = 60
<b>90050</b>	Surgical Gown, Sterile, XX-Large	Items/Case = 28
<b>90138</b>	Surgical Gowns Double Pack	Items/Case = 18
	Contents:	
72978	1-Surgical Gown, Large	
72988	1-Surgical Gown, X-Large	

All codes in **bold** indicate a sterile product.  
 NOTE: Please see package inserts for directions for use for each of these products.

# CARDIAC CATHETERISATION / RADIOLOGY DRAPES




## PACEMAKER DRAPE, 221cm x 343cm

- Two Clear Side Windows 55cm x 342cm
- CONTROL\* Plus Fabric Reinforcement
- Conformable Breather Bar 35cm
- Two Adhesive Fenestrations 15cm x 20cm
- Bilateral Pockets 15cm x 81cm

Code	Description	Packaging
<b>89368</b>	Pacemaker Drape, Sterile	Items/Case = 11
79368	Pacemaker Drape, Non-Sterile	Items/Case = 12

## BRACHIAL ANGIOGRAPHY DRAPE, 111cm x 78cm


- CONTROL\* Plus Fabric Reinforcement
- 7cm Adhesive Circular Fenestration



Code	Description	Packaging
<b>89471</b>	Brachial Angiography Drape, Sterile	Items/Case = 50

## RADIAL FEMORAL ANGIOGRAPHY DRAPE, 240cm x 380cm

- Two 63cm x 350cm Windows
- Two 12cm x 7cm Oval Radial Fenestrations
- Two 12cm Round Femoral Fenestrations



Code	Description	Packaging
<b>89583</b>	HALYARD BASICS* Radial Femoral Angiography Drape, Sterile	Items/Case = 12

## RADIAL FEMORAL ANGIOGRAPHY DRAPE, 244cm x 381cm

- Two 7cm Adhesive Circular Fenestrations
- Two 4cm x 6cm Oval Fenestrations
- Two 61cm x 381cm Clear Windows
- CONTROL\* Fabric Reinforcement

Code	Description	Packaging
<b>44976</b>	Radial/Femoral Angiography Drape, Sterile	Items/Case = 6

All codes in **bold** indicate a sterile product.

† Reinforcements include:


**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


# CARDIAC CATHETERISATION / RADIOLOGY DRAPES



## EXTRA LONG FEMORAL ANGIOGRAPHY DRAPE WITH WINDOW, 208cm x 343cm

- Two 9cm Adhesive Circular Fenestrations
- One 68cm x 315cm Windows
- CONTROL\* Plus Fabric Reinforcement


Code	Description	Packaging
<b>89709</b>	Extra Long Femoral Angiography Drape with Window, Sterile	Items/Case = 14
79709	Extra Long Femoral Angiography Drape with Window, Non-Sterile	Items/Case = 16



## EXTRA LONG FEMORAL ANGIOGRAPHY DRAPE WITH WINDOWS, 221cm x 343cm

- Two 9cm Adhesive Circular Fenestrations
- Two 68cm x 343cm Windows
- CONTROL\* Plus Fabric Reinforcement

Code	Description	Packaging
<b>89711</b>	Extra Long Femoral Angiography Drape with Dual Windows, Sterile	Items/Case = 16
13296	Extra Long Femoral Angiography Drape with Dual Windows, Non-Sterile	Items/Case = 180



## EXTRA-LONG FEMORAL ANGIOGRAPHY DRAPE, 220cm x 342cm

- Extra-Long 220cm x 342cm Film Based Drape
- 81cm x 342cm Absorbent Fabric Reinforcement
- 78cm x 342cm Patient Comfort Layer
- Two 8cm Adhesive Circular Fenestrations
- Two 68cm x 342cm Clear Plastic Windows

Code	Description	Packaging
<b>89749</b>	Extra-Long Femoral Angiography Drape, Sterile	Items/Case = 16
79749	Extra-Long Femoral Angiography Drape, Non-Sterile	Items/Case = 144

All codes in **bold** indicate a sterile product.

† Reinforcements include:


**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# CARDIOVASCULAR DRAPES




187cm x 241cm

## PERINEAL DRAPE, 187cm x 241cm

- Split 44cm x 95cm
- CONTROL\* Plus Fabric Reinforcement
- Positioning Tape

Code	Description	Packaging
76356	Perineal Drape, Non-Sterile	Items/Case = 17




## CVARTS\* HEAVY-DUTY UNIVERSAL DRAPE SYSTEM

- Heavy-Duty Fabric
- CONTROL\* Plus Fabric Reinforcement

Code	Description	Packaging
<b>89454</b>	CVARTS* Heavy-Duty Universal Drape System, Sterile	Items/Case = 5

Contents:


- 1-CVARTS\* Top Drape, 254cm x 152cm
- 1-CVARTS\* Bottom Drape, 193cm x 152cm
- 2-CVARTS\* Side Drapes, 96cm x 213cm, Longitudinal Flaps with Hook and Loop



## CVARTS\* HEAVY-DUTY CV SPLIT DRAPE, 274cm x 203cm x 315cm

- Armboard Covers
- 15cm x 213cm Split
- Longitudinal Flaps with Hook and Loop
- Heavy-Duty Fabric
- CONTROL\* Plus Fabric Reinforcement

Code	Description	Packaging
<b>89452</b>	CVARTS* Heavy-Duty CV Split Drape, Sterile	Items/Case = 6
79452	CVARTS* Heavy-Duty CV Split Drape, Non-Sterile	Items/Case = 7



## CVARTS\* BILATERAL UNDERLEG / PERINEAL DRAPE, 213cm x 203cm

- Perineal Cover
- Umbilicus Cover
- CONTROL\* Plus Fabric Reinforcement

Code	Description	Packaging
<b>89453</b>	CVARTS* Bilateral Underleg / Perineal Drape, Sterile	Items/Case = 14
79453	CVARTS* Bilateral Underleg / Perineal Drape, Non-Sterile	Items/Case = 18

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency


**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.



# ENT DRAPES & PACKS




79351




79430

## EENT PACK I

Code	Description	Packaging
<b>88711</b>	EENT Pack I	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
<b>89618</b>	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
<b>95111</b>	71311	1-ULTRA Surgical Gown in Overwrap, with Towel, Large
<b>89601</b>	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm
		1-Absorbent Towel, 38cm x 56cm
		1-Suture Bag
<b>89431</b>	79430	1-Head Drape, 127cm x 91cm
<b>89352</b>	79351	1-EENT Split Drape with CONTROL* Reinforcement, 193cm x 315cm with 7cm x 71cm Split




79351



79430

## EENT PACK II

Code	Description	Packaging
<b>88716</b>	EENT Pack II	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
<b>89618</b>	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
<b>95111</b>	71311	1-ULTRA Surgical Gown In Overwrap, Fabric Reinforced, with Towel, Large
<b>95211</b>	74911	1-ULTRA Surgical Gown, Fabric Reinforced, with Towel, Large
		1-Microsurgery Mayo Stand Cover, 58cm x 135cm
		1-Suture Bag
<b>89431</b>	79430	1-Head Drape, 127cm x 91cm
<b>89352</b>	79351	1-EENT Split Drape with CONTROL* Reinforcement, 193cm x 315cm with 7cm x 71cm Split



## ENT PACK III

Code	Description	Packaging
<b>88223</b>	ENT Pack III	Items/Case = 6
Sterile Component Code	Non-Sterile Component Code	Contents
<b>89603</b>	79603	1-Back Table Cover, Zone-Reinforced, 152cm x 191cm
		1-Mayo Stand Cover 80cm x 145cm
		1-Suture Bag
<b>89431</b>	79430	1-Head Drape 91cm x 127cm:
		• Turban
		• Positioning Tape
<b>89352</b>	79351	1-ENT Split Drape 193cm x 315cm:
		• with 7 x 71cm Split
		• Positioning Tape
		• CONTROL* Fabric Reinforcement
		• Tube Holders

All codes in **bold** indicate a sterile product.

† Reinforcements include:


**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

# ENT DRAPES & PACKS



## EENT PACK IV


Code	Description	Packaging
<b>88719</b>	EENT Pack IV	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
<b>89618</b>	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm 1-Surgical Gown with Set-in Sleeves, In Overwrap, with Towel, Large
<b>89701</b>	79720	1-Surgical Gown, with Set-in Sleeves, Large
<b>89601</b>	79601	1-Absorbent Towel, 38cm x 56cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
<b>89438</b>	79438	1-Head Drape, 104cm x 81cm 1-EENT Split Drape with Instrument Pad Reinforcement, 193cm x 315cm with 8cm x 64cm Split



## EENT SPLIT DRAPE, 193cm x 315cm

- CONTROL\* Fabric Reinforcement
- 7cm x 71cm Split
- Tube Holders
- Instrument Pad


Code	Description	Packaging
<b>89352</b>	EENT Split Drape, Sterile	Items/Case = 18
79351	EENT Split Drape, Non-Sterile	Items/Case = 25



## HEAD DRAPE, 127cm x 91cm

- Turban, 83cm x 50cm

Code	Description	Packaging
<b>89431</b>	Head Drape, Sterile	Items/Case = 38
79430	Head Drape, Non-Sterile	Items/Case = 60



## HEAD DRAPE, 104cm x 81cm

- Detachable Turban, 91cm x 61cm

Code	Description	Packaging
<b>89438</b>	Head Drape, Sterile	Items/Case = 38

All codes in **bold** indicate a sterile product.

† Reinforcements include:


**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# ENT DRAPES & PACKS




## DENTAL DRAPE, 120cm x 150cm

- Teardrop Fenestration 10cm x 15cm; 40cm from the top
- Absorbent Reinforcement 30cm x 40cm

Code	Description	Packaging
------	-------------	-----------

<b>89543</b>	HALYARD BASICS* Dental Drape, Sterile	Items/Case = 15
--------------	---------------------------------------	-----------------



## EENT FENESTRATED DRAPE, 193cm x 229cm

- 6cm x 5cm Fenestration
- CONTROL Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
------	-------------	-----------

<b>89421</b>	EENT Fenestrated Drape, Sterile	Items/Case = 18
--------------	---------------------------------	-----------------

All codes in **bold** indicate a sterile product.

† Reinforcements include:


**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

# EQUIPMENT COVERS



## MAYO STAND COVERS

Code	Description	Packaging
<b>88665</b>	Mayo Stand Cover, Reinforced HALYARD BASICS* 58cm x 137cm, Sterile	Items/Case = 80
78656	Mayo Stand Cover, Reinforced HALYARD BASICS* 58cm x 137cm, Non-Sterile	Items/Case = 140
<b>89601</b>	Mayo Stand Cover, Reinforced, 58cm x 137cm, Sterile	Items/Case = 54
79607	Mayo Stand Cover, Reinforced, 58cm x 137cm, Non-Sterile	Items/Case = 1200
<b>89602</b>	Mayo Stand Cover, Reinforced, X-Large, 76cm x 145cm, Sterile	Items/Case = 44
79602	Mayo Stand Cover, Reinforced, X-Large, 76cm x 145cm, Non-Sterile	Items/Case = 120
<b>14007</b>	Mayo Stand Cover, Padded, X-Large, 76cm x 145cm, Sterile	Items/Case = 44
76358	Mayo Stand Cover, Padded, X-Large, 76cm x 145cm, Non-Sterile	Items/Case = 68
<b>89603</b>	Mayo Stand Cover, Reinforced, XX-Large, 80cm x 145cm, Sterile	Items/Case = 44
79603	Mayo Stand Cover, Reinforced, XX-Large, 80cm x 145cm, Non-Sterile	Items/Case = 120



## BACK TABLE COVERS

Code	Description	Packaging
<b>88666</b>	Back Table Cover, Zone Reinforced HALYARD BASICS* 112cm x 229cm, Sterile	Items/Case = 45
78664	Back Table Cover, Zone Reinforced HALYARD BASICS* 112cm x 229cm, Non-Sterile	Items/Case = 60
<b>89632</b>	Back Table Cover, HALYARD BASICS*, 100cm x 150cm, Sterile	Items/Case = 50
79482	Back Table Cover, 115cm x 152cm, Non-Sterile, Standard Fold	Items/Case = 90
<b>42217</b>	Back Table Cover, Standard, 112cm x 229cm, Sterile	Items/Case = 28
42218	Back Table Cover, Standard, 112cm x 229cm, Non-Sterile, Unfolded	Items/Case = 640
42231	Back Table Cover, Standard, 112cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 400
<b>42224</b>	Back Table Cover, Standard, 152cm x 229cm, Sterile	Items/Case = 36
42225	Back Table Cover, Standard, 152cm x 229cm, Non-Sterile, Unfolded	Items/Case = 480
42233	Back Table Cover, Standard, 152cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 400
<b>42306</b>	Back Table Cover, Heavy Duty, 152cm x 229cm, Sterile	Items/Case = 28
42307	Back Table Cover, Heavy Duty, 152cm x 229cm, Non-Sterile, Unfolded	Items/Case = 400
42366	Back Table Cover, Heavy Duty, 152cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 200
<b>42309</b>	Back Table Cover, Heavy Duty, 177cm x 279cm, Sterile	Items/Case = 15
42310	Back Table Cover, Heavy Duty, 177cm x 279cm, Non-Sterile, Unfolded	Items/Case = 280
42367	Back Table Cover, Heavy Duty, 177cm x 279cm, Non-Sterile, Fan Fold	Items/Case = 160
<b>29258</b>	Back Table Cover, Padded Heavy Duty, 152cm x 229cm, Sterile	Items/Case = 20
29259	Back Table Cover, Padded Heavy Duty, 152cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 22
29269	Back Table Cover, Padded Heavy Duty, 152cm x 229cm, Non-Sterile, Unfolded	Items/Case = 300
<b>37498</b>	Back Table Cover, Two Tier Heavy Duty, 152cm, Sterile	Items/Case = 10

## OVERHEAD / CV TABLE COVERS

Code	Description	Packaging
<b>89621</b>	Overhead Table Cover, 216cm x 193cm, 112cm x 152cm Skirt, Sterile	Items/Case = 14

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency


**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


# GENERAL DRAPES & PACKS

## ACCESSORIES




### ABSORBENT TOWELS, 38cm x 56cm

Code	Description	Packaging
<b>89701</b>	Absorbent Towels, White, Sterile, Pack of 2	Items/Case = 100
79720	Absorbent Towels, White, Non-Sterile	Items/Case = 1000



### TAPE STRIPS, 10cm x 55cm


Code	Description	Packaging
<b>89053</b>	Tape Strips, Sterile, Pack of 2	Items/Case = 300



### DIATHERMY POUCH, 38cm x 41cm

- Adhesive Backing Strip

Code	Description	Packaging
<b>89945</b>	Diathermy Pouch, Sterile	Items/Case = 150
79945	Diathermy Pouch, Non-Sterile	Items/Case = 220



### UTILITY DRAPE WITH TAPE, 38cm x 64cm

Code	Description	Packaging
<b>89731</b>	Utility Drape with Tape, Sterile, Pack of 2	Items/Case = 100
79731	Utility Drape with Tape, Non-Sterile	Items/Case = 550

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# GENERAL DRAPES & PACKS

## BASIC PACKS

### BASIC PACK VI

<b>Code</b>	<b>Description</b>	<b>Packaging</b>	
<b>88161</b>	Basic Pack VI	Items/Case = 7	
<b>Sterile Component Code</b>	Non-Sterile Component Code	<b>Contents</b>	
<b>89618</b>	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm	
<b>95111</b>	71311	1-ULTRA Surgical Gown, In Overwrap with Towel, Large	
<b>95121</b>	71321	1-ULTRA Surgical Gown, X-Large	
<b>89601</b>	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm	
<b>89731</b>	79731	4-Utility Drapes, 66cm x 38cm, with Tape 1-Suture Bag	
<b>89101</b>	79101	1-Half Drape, 152cm x 112cm	
<b>89121</b>	79121	1-Large Drape, 152cm x 193cm	

### BASIC PACK XV

<b>Code</b>	<b>Description</b>	<b>Packaging</b>	
<b>88178</b>	Basic Pack XV	Items/Case = 15	
<b>Sterile Component Code</b>	Non-Sterile Component Code	<b>Contents</b>	
<b>89618</b>	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm	
<b>89601</b>	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm	
<b>89121</b>	79121	1-Large Drape, 152cm x 193cm	

### SET-UP PACK I

<b>Code</b>	<b>Description</b>	<b>Packaging</b>	
<b>88171</b>	Set-Up Pack I	Items/Case = 14	
<b>Sterile Component Code</b>	Non-Sterile Component Code	<b>Contents</b>	
<b>89618</b>	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm	
<b>95111</b>	71311	1-Ultra Surgical Gown In Overwrap with Towel, Large	
<b>95121</b>	71321	1-ULTRA Surgical Gown, X-Large	
<b>89601</b>	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag	

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency


**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


# GENERAL DRAPES & PACKS

## BASIC DRAPES




### UTILITY DRAPE, 90cm x 90cm

Code	Description	Packaging
<b>89353</b>	Utility Drape, Sterile	Items/Case = 120
76323	Utility Drape, Non-Sterile	Items/Case = 200




### HALF DRAPE, 152cm x 101cm

Code	Description	Packaging
<b>47614</b>	Half Drape, Sterile	Items/Case = 54
47615	Half Drape, Non-Sterile	Items/Case = 1000




### MEDIUM DRAPE, 193cm x 112cm

Code	Description	Packaging
<b>89111</b>	Medium Drape, Sterile	Items/Case = 50



### MEDIUM DRAPE WITH TAPE, 193cm x 105cm

Code	Description	Packaging
<b>89023</b>	Medium Drape with Tape, Sterile	Items/Case = 50



### MEDIUM DRAPE, 185cm x 101cm

Code	Description	Packaging
47619	Medium Drape, Non-Sterile	Items/Case = 800

All codes in **bold** indicate a sterile product.

† Reinforcements include:


**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

# GENERAL DRAPES & PACKS



## LARGE DRAPE, 139cm x 193cm


Code	Description	Packaging
<b>47522</b>	Large Drape, Sterile	Items/Case = 36
47624	Large Drape, Non-Sterile	Items/Case = 600



## THREE-QUARTER DRAPE, 152cm x 193cm


- 39CONTROL\* Plus Fabric Reinforcement
- No Adhesive Tape

Code	Description	Packaging
<b>47628</b>	Three-Quarter Drape, Sterile	Items/Case = 20



## X-LARGE DRAPE, 193cm x 243cm


Code	Description	Packaging
<b>47630</b>	X-Large Drape, Sterile	Items/Case = 20
47631	X-Large Drape, Non-Sterile	Items/Case = 24



## SMALL ABSORBENT DRAPE WITH TAPE, 57cm x 60cm

- CONTROL\* Plus Absorbent Fabric

Code	Description	Packaging
<b>89828</b>	Small Absorbent Drape with Tape, Sterile	Items/Case = 72
79828	Small Absorbent Drape with Tape, Non-Sterile	Items/Case = 108



## STANDARD ABSORBENT DRAPE WITH TAPE, 91cm x 91cm

- CONTROL\* Plus Absorbent Fabric

Code	Description	Packaging
<b>89343</b>	Standard Absorbent Drape with Tape, Sterile	Items/Case = 50
79343	Standard Absorbent Drape with Tape, Non-Sterile	Items/Case = 70

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric


**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.



# GENERAL DRAPES & PACKS


## BREAST / CHEST DRAPES



### CHEST U-DRAPE, 193cm x 305cm

- CONTROL\* Plus Fabric Reinforcement
- 30cm x 76cm Split
- Tube Holders

Code	Description	Packaging
79335	Chest U-Drape, Non-Sterile	Items/Case = 22



### BREAST DRAPE, 254cm x 193cm x 315cm

- 15cm x 15cm Fenestration
- CONTROL\* Fabric Reinforcement
- Hook and Loop Tube and Line Control
- Armboard Covers

Code	Description	Packaging
<b>89236</b>	Breast Drape, Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency


**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# GENERAL DRAPES & PACKS


## FENESTRATED DRAPES



### EPIDURAL DRAPE WITH TAPE, 50cm x 75cm

- 8cm x 10cm Oval Fenestration, Centred
- Spunbond Film
- Positioning Tape


Code	Description	Packaging
<b>89636</b>	HALYARD BASICS* Epidural Drape with Tape, Sterile	Items/Case = 40



### VARIABLE FENESTRATED DRAPE, 75cm x 100cm

- Variable Fenestration
- Half Circle On Each Drape with 12cm Diameter
- Spunbond Film
- Positioning Tape


Code	Description	Packaging
<b>89680</b>	HALYARD BASICS* Variable Fenestrated Drape, Sterile	Items/Case = 30



### CV LINE DRAPE, 160cm x 240cm

- CONTROL\* Fabric Reinforcement
- 12cm x 12cm Fenestration

Code	Description	Packaging
<b>89835</b>	CV Line Drape, Sterile	Items/Case = 10
79835	CV Line Drape, Non-Sterile	Items/Case = 14



### CV FULL BODY DRAPE, 172cm x 279cm

- CONTROL\* Plus Fabric Reinforcement
- 10cm Circular Adhesive

Code	Description	Packaging
31445	CV Full Body Drape, Non-Sterile	Items/Case = 27

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency


**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# GENERAL DRAPES & PACKS


## LAPAROSCOPY / LAPAROTOMY DRAPES



### LAPAROSCOPIC / CHOLECYSTECTOMY DRAPE, 264cm x 193cm x 305cm

- 30cm x 34cm Adjustable Fenestration, Lateral Offset
- Hook and Loop Tube and Line Control
- CONTROL\* Instrument Pad Reinforcement
- Armboard Covers


Code	Description	Packaging
<b>89229</b>	Laparoscopic / Cholecystectomy Drape, Sterile	Items/Case = 12



### LAPAROSCOPIC ABDOMINAL DRAPE, 254cm x 193cm x 315cm

- 30cm x 33cm Abdominal Fenestration
- CONTROL\* Fabric Reinforcement
- Clear Instrument and Line Control Pockets
- Hook and Loop Tube and Line Control
- Armboard Covers


Code	Description	Packaging
<b>89233</b>	Laparoscopic Abdominal Drape, Sterile	Items/Case = 15
79233	Laparoscopic Abdominal Drape, Non-Sterile	Items/Case = 18



### LAPAROSCOPIC DRAPE, 264cm x 89cm x 315cm

- 31cm x 33cm Fenestration
- CONTROL\* Plus Fabric Reinforcement
- Clear Film Pouch
- Tube Holders
- Armboard Covers
- Under Leg Wings

Code	Description	Packaging
<b>29240</b>	Laparoscopic Drape, Sterile	Items/Case = 15
29241	Laparoscopic Drape, Non-Sterile	Items/Case = 19



### HALYARD BASICS\* LAPAROSCOPIC/CHOLECYSTECTOMY DRAPE, 264cm x 193cm x 305cm

- 31cm x 33cm Abdominal Fenestration
- Hook and Loop Tube and Line Control
- Fabric Reinforcement
- Armboard Covers

Code	Description	Packaging
<b>89538</b>	Laparoscopic/Cholecystectomy Drape, Sterile	Items/Case = 25
79538	Laparoscopic/Cholecystectomy Drape, Non-Sterile	Items/Case = 26

All codes in **bold** indicate a sterile product.

† Reinforcements include:


**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


# GENERAL DRAPES & PACKS



## PAEDIATRIC LAPAROTOMY DRAPE, 152cm x 274cm

- 7cm x 20cm Fenestration
- CONTROL\* Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
<b>89241</b>	Paediatric Laparotomy Drape, Sterile	Items/Case = 15
79240	Paediatric Laparotomy Drape, Non-Sterile	Items/Case = 25



## LAPAROTOMY DRAPE WITH INCISE, 254cm x 183cm x 305cm

- 20cm x 33cm Adhesive Incise Area
- CONTROL\* Plus Fabric Reinforcement
- Tube Holders
- Armboard Covers

Code	Description	Packaging
<b>89282</b>	Laparotomy Drape with Incise, Sterile	Items/Case = 15
79282	Laparotomy Drape with Incise, Non-Sterile	Items/Case = 15



## ROBOTIC ABDOMINAL / PELVIC DRAPE, 264cm x 406cm

- 30cm x 10cm x 48cm Trapezoidal Abdominal Fenestration
- Clear Film Poly Pouches
- Tube Holders and Line Control
- CONTROL\* Plus Fabric Reinforcement

Code	Description	Packaging
<b>29283</b>	Robotic Abdominal/Pelvic Drape, Sterile	Items/Case = 8
29284	Robotic Abdominal/Pelvic Drape, Non-Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency


**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# GENERAL DRAPES & PACKS


## MINOR PROCEDURE DRAPES



### MINOR PROCEDURE FENESTRATED DRAPE, 76cm x 76cm

- 4cm x 5cm Rectangular Fenestration
- CONTROL\* Fabric Reinforcement
- Positioning Tape


Code	Description	Packaging
<b>89206</b>	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 100
79206	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 250



### MINOR PROCEDURE FENESTRATED DRAPE, 122cm x 152cm

- 7cm Diameter Circular Fenestration
- CONTROL\* Fabric Reinforcement
- Positioning Tape


Code	Description	Packaging
<b>89207</b>	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 40
79207	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 60



### MINOR PROCEDURE FENESTRATED DRAPE, 122cm x 152cm

- 12cm x 15cm Rectangular Fenestration
- CONTROL\* Fabric Reinforcement
- Positioning Tape


Code	Description	Packaging
<b>89208</b>	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 40
79208	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 60



### MINOR PROCEDURE FENESTRATED DRAPE, 193cm x 315cm

- 15cm x 15cm Fenestration
- CONTROL\* Fabric Reinforcement
- Positioning Tape

Code	Description	Packaging
<b>89209</b>	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 19
79209	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 23



### CIRCUMCISION DRAPE, 31cm x 51cm

- 2.5cm Circular Fenestration
- CONTROL\* Fabric Reinforcement

Code	Description	Packaging
<b>44972</b>	Circumcision Drape, Sterile	Items/Case = 40

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency


**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# GENERAL DRAPES & PACKS

## THYROID DRAPES



### THYROID DRAPE, 193cm x 315cm

- 12cm x 12cm Diamond Fenestration
- CONTROL\* Plus Fabric Reinforcement
- Instrument Pad
- Tube Holders

Code	Description	Packaging
<b>89255</b>	Thyroid Drape, Sterile	Items/Case = 13
79255	Thyroid Drape, Non-Sterile	Items/Case = 18

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric


**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# GENERAL DRAPES & PACKS


## UNIVERSAL PACKS

### UNIVERSAL PACK – MAJOR WITH LONG SIDE DRAPES



<b>Code</b>	<b>Description</b>	<b>Packaging</b>
<b>88759</b>	Universal Pack – Major	Items/Case = 5
<b>Sterile Component Code</b>	<b>Non-Sterile Component Code</b>	<b>Contents</b>
<b>89601</b>	79601	1-Back Table Cover, Fully Reinforced, 137cm x 248cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Absorbent Towel, 38cm x 56cm 1-Suture Bag
<b>34902</b>	34903	1-Bottom Drape with CONTROL* Plus Reinforcement, 193cm x 195cm
<b>34906</b>	34901	1-Top Drape with CONTROL* Plus Reinforcement, 273cm x 127cm
<b>34904</b>	34905	2-Side Drapes with CONTROL* Plus Reinforcement, 195cm x 105cm

### UNIVERSAL PACK – MINOR WITH SHORT SIDE DRAPES



<b>Code</b>	<b>Description</b>	<b>Packaging</b>
<b>88769</b>	Universal Pack – Minor	Items/Case = 14
<b>Sterile Component Code</b>	<b>Non-Sterile Component Code</b>	<b>Contents</b>
<b>89611</b>	79613	1-Back Table Cover, Fully-Reinforced, 112cm x 224cm
<b>89601</b>	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
	79097	2-Side Drapes with CONTROL* Reinforcement, 114cm x 84cm
	79096	1-Bottom Drape with CONTROL* Reinforcement, 165cm x 193cm
	79095	1-Top Drape with CONTROL* Reinforcement, 276cm x 150cm

All codes in **bold** indicate a sterile product.

† Reinforcements include:


**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

# GENERAL DRAPES & PACKS



## UNIVERSAL PACK I

Code	Description	Packaging
<b>88878</b>	Universal Pack – I	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
<b>88666</b>	78658	1-HALYARD BASICS* Back Table Cover, Zone Reinforced, 112cm x 229cm
<b>88665</b>	78656	1-HALYARD BASICS* Mayo Stand Cover, Reinforced, 58cm x 137cm
	79756	2-HALYARD BASICS* Side Drape, 193 x 112cm
	79754	1-HALYARD BASICS* Top Drape, 274 x 127cm
	79755	1-HALYARD BASICS* Bottom Drape, 193 x 173cm



## UNIVERSAL PACK III WITH MEDIUM SIDE DRAPES

Code	Description	Packaging
<b>88774</b>	Universal Pack – III	Items/Case = 8
Sterile Component Code	Non-Sterile Component Code	Contents
<b>89601</b>	79601	1-Back Table Cover, Fully-Reinforced, 133cm x 248cm 1-Mayo Stand Cover, 58cm x 137cm 2-Absorbent Towels, 38cm x 56cm 1-Suture Bag 2-Side Drapes with CONTROL* Reinforcement, 145cm x 99cm 1-Bottom Drape with CONTROL* Reinforcement, 193cm x 193cm 1-Top Drape with CONTROL* Reinforcement, 276cm x 150cm

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency


**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.




# GENERAL DRAPES & PACKS



## U-BAR\* PACK

Code	Description	Packaging
<b>88334</b>	U-BAR* Pack	Items/Case = 16
Sterile Component Code	Non-Sterile Component Code	Contents
<b>89601</b>	79601	1-Back Table Cover, Zone Reinforced, 112cm x 224cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag 1-U-Drape with CONTROL* Reinforcement, 180cm x 315cm, with 10cm x 102cm Split 1-Bar Drape with CONTROL* Reinforcement, 105cm x 228cm



## U-BAR\* PACK – IV

Code	Description	Packaging
<b>88773</b>	U-BAR* Pack – IV	Items/Case = 10
Sterile Component Code	Non-Sterile Component Code	Contents
<b>89601</b>	79601	1-Back Table Cover, Reinforced, 133cm x 248cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-U-Drape with CONTROL* Reinforcement, 192cm x 315cm, with 11cm x 51cm Split 1-Bar Drape with CONTROL* Reinforcement, 248cm x 145cm

## HALYARD BASICS\* U-BAR PACK I

Code	Description	Packaging
<b>88875</b>	HALYARD BASICS* U-BAR Pack I	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
<b>89532</b>	79532	1-HALYARD BASICS* U-Drape, 203cm x 315cm, with 10cm x 112cm Split
<b>89533</b>	79533	1-HALYARD BASICS* Bar Drape, 244cm x 112cm
<b>88666</b>	78658	1-HALYARD BASICS* Back Table Cover, Zone-Reinforced, 112cm x 229cm
<b>88665</b>	78656	1-HALYARD BASICS* Mayo Stand Cover, Reinforced, 58cm x 137cm

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency


**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# GENERAL DRAPES & PACKS


## UNIVERSAL DRAPES



### TOP DRAPE, 273cm x 127cm

- CONTROL\* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88759


Code	Description	Packaging
<b>34906</b>	Top Drape, Sterile	Items/Case = 15
34901	Top Drape, Non-Sterile	Items/Case = 275



### BOTTOM DRAPE, 193cm x 195cm

- CONTROL\* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88759

Code	Description	Packaging
<b>34902</b>	Bottom Drape, Sterile	Items/Case = 12
34903	Bottom Drape, Non-Sterile	Items/Case = 242



### SIDE DRAPE, 195cm x 105cm

- CONTROL\* Plus Fabric Reinforcement
- Positioning Tape
- From Universal Pack 88759

Code	Description	Packaging
<b>34904</b>	Side Drape, Sterile	Items/Case = 20
34905	Side Drape, Non-Sterile	Items/Case = 440

All codes in **bold** indicate a sterile product.

† Reinforcements include:


**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


# GENERAL DRAPES & PACKS



## TOP DRAPE, 276cm x 150cm

- CONTROL\* Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88769


Code	Description	Packaging
79095	Top Drape, Non-Sterile	Items/Case = 60



## BOTTOM DRAPE, 165cm x 193cm

- CONTROL\* Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88769

Code	Description	Packaging
79096	Bottom Drape, Non-Sterile	Items/Case = 70



## SIDE DRAPE, 114cm x 84cm

- CONTROL\* Fabric Reinforcement
- Positioning Tape
- From Universal Pack 88769

Code	Description	Packaging
79097	Side Drape, Non-Sterile	Items/Case = 200

All codes in **bold** indicate a sterile product.

† Reinforcements include:


**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# NEUROSURGICAL DRAPE



## CRANIOTOMY DRAPE WITH POUCH, 310cm x 188cm x 340cm

- CONTROL\* Plus Fabric Reinforcement
- 20cm x 30cm Oval Fenestration
- Clear Gusseted Anesthesia Panels
- Screen and Suction Port

Code	Description	Packaging
<b>89062</b>	Craniotomy Drape with Pouch, Sterile	Items/Case = 9
79075	Craniotomy Drape with Pouch, Non-Sterile	Items/Case = 15

All codes in **bold** indicate a sterile product.

† Reinforcements include:


**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# OBSTETRICS & GYNAECOLOGY DRAPES & PACKS




## C-SECTION DRAPES

### C-SECTION DRAPE WITH CLEAR SCREEN, 264cm x 193cm x 305cm

- Triangular Fenestration
- 360° Fluid Collection Pouch
- Clear Screen
- Tube Holders
- Armboard Covers


Code	Description	Packaging
<b>44968</b>	C-Section Drape with Clear Screen, Sterile	Items/Case = 6
44967	C-Section Drape with Clear Screen, Non-Sterile	Items/Case = 100



### C-SECTION DRAPE WITH POUCH AND FENESTRATION, 264cm x 193cm x 305cm

- 36cm x 31cm Adhesive Area with Fenestration
- 360° Fluid Collection Pouch
- CONTROL\* Fabric Reinforcement
- Tube Holders
- Armboard Covers


Code	Description	Packaging
<b>44978</b>	C-Section Drape with Pouch and Fenestration, Sterile	Items/Case = 6
44977	C-Section Drape with Pouch and Fenestration, Non-Sterile	Items/Case = 100



### C-SECTION DRAPE WITH FENESTRATED INCISE, 254cm x 183cm x 305cm

- 38cm x 33cm Adhesive Incise Area with Triangular Fenestration
- 270° Fluid Collection Pouch, Malleable Band, with Suction Port
- Tube Holders
- Armboard Covers

Code	Description	Packaging
<b>89528</b>	C-Section Drape with Fenestrated Incise, Sterile	Items/Case = 8
79528	C-Section Drape with Fenestrated Incise, Non-Sterile	Items/Case = 150



### C-SECTION DRAPE WITH FULL INCISE, 254cm x 183cm x 305cm

- 38cm x 33cm Adhesive Incise Area
- 270° Fluid Collection Pouch, Malleable Band, with Suction Port
- Tube Holders
- Armboard Covers

Code	Description	Packaging
<b>89529</b>	C-Section Drape with Full Incise, Sterile	Items/Case = 8
79529	C-Section Drape with Full Incise, Non-Sterile	Items/Case = 150

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric


**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


# OBSTETRICS & GYNAECOLOGY DRAPES & PACKS

## LAPAROSCOPY PACKS

### LAPAROSCOPY PACK I



Code	Description	Packaging
<b>88088</b>	Laparoscopy Pack I	Items/Case = 10
Sterile Component Code	Non-Sterile Component Code	Contents
<b>89618</b>	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm 1-Absorbent Towel, 38cm x 56cm
<b>89411</b>		1-Under Buttocks Drape, Plastic, 75cm x 112cm
<b>89408</b>	79408	2-Leggings with 15cm Cuff, 78cm x 122cm
<b>89119</b>	79117	1-Laparoscopy "T" Drape with CONTROL* Plus Reinforcement, 249cm x 104cm x 251cm <ul style="list-style-type: none"> <li>• 15cm x 20cm Rectangular Abdominal Fenestration</li> <li>• 7cm x 15cm Perineal Fenestration with Flap</li> <li>• Tube Holders</li> </ul>



### LAPAROSCOPY / PELVISCOPY PACK I

Code	Description	Packaging
<b>88222</b>	Laparoscopy/Pelviscopy Pack	Items/Case = 6
Sterile Component Code	Non-Sterile Component Code	Contents
<b>89618</b>	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
<b>95111</b>	71311	1-ULTRA Surgical Gown, In Overwrap, with Towel, Large
<b>95111</b>	71311	1-ULTRA Surgical Gown, with Towel, Large
<b>89601</b>	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
<b>89413</b>		1-Under Buttocks Drape, Plastic-Reinforced, Large 101cm x 116cm
<b>89408</b>	79408	2-Leggings, with 15cm Cuff, Sterile Pair, 78cm x 122cm 1-Laparoscopy / Pelviscopy Drape with CONTROL* Plus Reinforcement (without Leggings), 264cm x 89cm x 264cm <ul style="list-style-type: none"> <li>• 18cm x 23cm Trapezoidal Abdominal Fenestration</li> <li>• 8cm x 15cm Perineal Fenestration with Flap</li> </ul>

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency


**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# OBSTETRICS & GYNAECOLOGY DRAPES & PACKS


## LAPAROSCOPY DRAPES



### LAPAROSCOPY / PELVISCOPY DRAPE WITH FLUID COLLECTION POUCH, 284cm x 76cm x 314cm

- 17cm x 23cm Trapezoidal Abdominal Fenestration
- 7cm x 20cm Perineal Fenestration with Flap
- Fluid Collection Pouch with Screen and Suction Port
- CONTROL\* Plus Fabric Reinforcement
- Attached Leggings
- Hook and Loop Line Controls
- Armboard Covers


Code	Description	Packaging
<b>89219</b>	Laparoscopy/Pelviscopy Drape Drape with Pouch, Sterile	Items/Case = 10



### LAPAROSCOPY / PELVISCOPY DRAPE, 254cm x 76cm x 312cm

- 17cm x 23cm Trapezoidal Abdominal Fenestration
- 7cm x 20cm Perineal Fenestration with Flap
- CONTROL\* Plus Fabric Reinforcement
- Attached Leggings
- Hook and Loop Tube and Line Controls

Code	Description	Packaging
<b>89224</b>	Laparoscopy/Pelviscopy Drape, Sterile	Items/Case = 10
79224	Laparoscopy/Pelviscopy Drape, Non-Sterile	Items/Case = 10



### LAPAROSCOPY "T" DRAPE, 249cm x 104cm x 251cm

- 15cm x 20cm Rectangular Abdominal Fenestration
- 7cm x 15cm Perineal Fenestration with Flap
- Tube Holders

Code	Description	Packaging
79117	1-Laparoscopy "T" Drape with CONTROL* Plus Reinforcement, Non-Sterile	Items/Case = 26

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency


**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


# OBSTETRICS & GYNAECOLOGY DRAPES & PACKS

## LITHOTOMY PACKS



### LITHOTOMY PACK II

Code	Description	Packaging
<b>88521</b>	Lithotomy Pack II	Items/Case = 10
Sterile Component Code	Non-Sterile Component Code	Contents
<b>89618</b>	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
		1-Absorbent Towel, 38cm x 56cm
<b>89413</b>	79413	1-Under Buttocks Drape, Plastic-Reinforced, 102cm x 117cm
<b>89408</b>	79408	2-Leggings with 15cm Cuff, 78cm x 122cm
<b>89101</b>	79101	1-Half Drape, 152cm x 111cm



### LITHOTOMY PACK V

Code	Description	Packaging
<b>88551</b>	Lithotomy Pack V	Items/Case = 12
Sterile Component Code	Non-Sterile Component Code	Contents
<b>89618</b>	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
<b>89701</b>	79720	1-Absorbent Towel, 38cm x 56cm
<b>89408</b>	79408	2-Leggings, with 15cm Cuffs
	79551	1-Lithotomy "T" Drape with CONTROL* Plus Reinforcement, 152cm x 78cm x 190cm with 8cm x 15cm Perineal Fenestration

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric


**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.



# OBSTETRICS & GYNAECOLOGY DRAPES & PACKS

## LITHOTOMY DRAPES




### LITHOTOMY DRAPE, 254cm x 76cm x 221cm

- 7cm x 15cm Perineal Fenestration
- CONTROL\* Plus Fabric Reinforcement
- Tube Holders
- Attached Leggings

Code	Description	Packaging
------	-------------	-----------

<b>89203</b>	Lithotomy Drape, Sterile	Items/Case = 16
--------------	--------------------------	-----------------




### GYN / UROLOGY DRAPE WITH POUCH, 160cm x 284cm x 332cm

- Pouch with Mesh Screen and Port
- Attached Leggings
- CONTROL\* Plus Fabric Reinforcement
- 7cm x 15cm Perineal Fenestration
- Tube Holders
- Armboard Covers

Code	Description	Packaging
------	-------------	-----------

<b>89204</b>	Gyn/Urology Drape with Pouch, Sterile	Items/Case = 10
--------------	---------------------------------------	-----------------

79204	Gyn/Urology Drape with Pouch, Non-Sterile	Items/Case = 10
-------	---	-----------------




### LITHOTOMY "T" DRAPE, 152cm x 78cm x 190cm

- 7cm x 15cm Perineal Fenestration
- CONTROL\* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
------	-------------	-----------

79551	Lithotomy "T" Drape, Non-Sterile	Items/Case = 55
-------	----------------------------------	-----------------



### LEGGINGS, 78cm x 122cm, WITH 15cm CUFF

Code	Description	Packaging
------	-------------	-----------

<b>89408</b>	Leggings, Sterile Pair	Items/Case = 26
--------------	------------------------	-----------------

79408	Leggings, Non-Sterile	Items/Case = 70
-------	-----------------------	-----------------

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency


**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# OBSTETRICS & GYNAECOLOGY DRAPES & PACKS


## UNDER BUTTOCKS DRAPES



### UNDER BUTTOCKS DRAPE, ZONE-REINFORCED, MEDIUM, 96cm x 112cm


- CONTROL\* Plus Fabric Reinforcement

Code	Description	Packaging
79411	Under Buttocks Drape, Zone-Reinforced, Medium, Non-Sterile	Items/Case = 115



### UNDER BUTTOCKS DRAPE, PLASTIC REINFORCED, LARGE, 101cm x 116cm


Code	Description	Packaging
<b>89413</b>	Under Buttocks Drape, Plastic Reinforced, Large, Sterile	Items/Case = 50



### UNDER BUTTOCKS DRAPE WITH FLUID COLLECTION POUCH, 101cm x 112cm

- CONTROL\* Plus Absorbent Pad
- 72cm x 66cm Triangular Pouch

Code	Description	Packaging
<b>89414</b>	Under Buttocks Drape, with Fluid Collection Pouch, Sterile	Items/Case = 30
79414	Under Buttocks Drape, with Fluid Collection Pouch, Non-Sterile	Items/Case = 80



### UNDER BUTTOCKS DRAPE WITH FLUID COLLECTION POUCH, 101cm x 112cm

- Screen and Suction Port
- Adhesive Pouch Closure
- 69cm x 61cm x 7cm Triangular Pouch

Code	Description	Packaging
<b>89415</b>	Under Buttocks Drape with Fluid Collection Pouch, Sterile	Items/Case = 40
79415	Under Buttocks Drape with Fluid Collection Pouch, Non-Sterile	Items/Case = 100

All codes in **bold** indicate a sterile product.

† Reinforcements include:


**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


# OPHTHALMIC DRAPES & PACKS



## BASIC OPHTHALMIC DRAPE, 132cm x 142cm

- Lightweight Fabric
- 12cm x 9cm Adhesive Incise Area
- Fluid Collection Pouch, with Eye Wick

Code	Description	Packaging
<b>89304</b>	Basic Ophthalmic Drape, Sterile	Items/Case = 40
79304	Basic Ophthalmic Drape, Non-Sterile	Items/Case = 100



## OCUARTS\* FULL BODY OPHTHALMIC DRAPE, 183cm x 292cm

- 11cm x 11cm Adhesive Incise Area
- 5cm x 3cm Oval Fenestration
- 12cm x 15cm Optional Incise Drape
- 270° Fluid Collection Pouch, Foam Band
- Breather Bar
- Tube Holders

Code	Description	Packaging
<b>89356</b>	OCUARTS* Full Body Ophthalmic Drape, Sterile	Items/Case = 10
79356	OCUARTS* Full Body Ophthalmic Drape, Non-Sterile	Items/Case = 20

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric


**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


# ORTHOPAEDIC DRAPES & PACKS

## ORTHOPAEDIC PACKS

### ORTHO PACK I – MAJOR



Code	Description	Packaging
<b>88411</b>	Ortho Pack I – Major	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
<b>89618</b>	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
<b>44682</b>	44664	2-AERO CHROME* Breathable High Performance Surgical Gown, with Towel, X-Large
<b>89601</b>	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
<b>89331</b>	79330	1-Impervious Split Drape, 152cm x 178cm with 10cm x 53cm Split 1-Ortho U-Drape with CONTROL* Plus Reinforcement, 223cm x 274cm with 10cm x 102cm Split
<b>89381</b>	79360	1-Ortho Bar Drape with CONTROL* Plus Reinforcement, 254cm x 152cm



### ORTHOARTS\* HIP PACK WITH POCKETS

Code	Description	Packaging
<b>88452</b>	ORTHOARTS* Hip Pack with Pockets	Items/Case = 3
Sterile Component Code	Non-Sterile Component Code	Contents
<b>89906</b>	79006	1-Back Table Cover, Heavy-Duty, 152cm x 228cm
<b>95111</b>	71311	1-ULTRA Surgical Gown, in Overwrap, with Towel, Large
<b>44682</b>	44664	2-AERO CHROME* Breathable High Performance Surgical Gown, with Towel, X-Large
<b>89601</b>	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
<b>89821</b>	79515	1-Stockinette, Large
<b>89121</b>	79121	1-Large Drape, 152cm x 193cm
<b>89333</b>	79333	1-Impervious Split Drape, Large, 152cm x 205cm with 10cm x 79cm Split
<b>89345</b>	79345	1-ORTHOARTS* Hip Drape with CONTROL* Plus Reinforcement, 289cm x 223cm x 320cm, with 21cm x 17cm Conformable Teardrop Fenestration 1-Yankauer Suction Tip 1-304cm Suction Tubing

All codes in **bold** indicate a sterile product.

† Reinforcements include:


**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# ORTHOPAEDIC DRAPES & PACKS



## ORTHOARTS\* LOWER EXTREMITY

Code	Description	Packaging
<b>88471</b>	ORTHOARTS* Lower Extremity Pack	Items/Case = 3
<b>Sterile Component Code</b>	Non-Sterile Component Code	Contents
<b>89906</b>	79006	1-Back Table Cover, Heavy-Duty, 152cm x 228cm
<b>95111</b>	71311	1-ULTRA Surgical Gown, In 60cm x 120cm Overwrap, with Towel, Large
		1-Extremity Prep Pad
<b>95221</b>	74921	1-ULTRA Fabric-Reinforced Surgical Gown, with Towel, X-Large
<b>44686</b>	44668	1-AERO CHROME* Breathable High Performance Specialty Surgical Gown, X-Large, X-Long
<b>89601</b>	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm
		1-Suture Bag
		1-304cm Suction Tubing
<b>89821</b>	79515	1-Stockinette, Large
<b>89331</b>	79330	1-Impervious Split Drape, 152cm x 178cm with 10cm x 53cm Split
<b>89278</b>	13233	1-ORTHOARTS* Lower Extremity Drape with CONTROL* Plus Reinforcement, with 6cm Conformable Circular Fenestration, 3-Drape Tape Strips
<b>89053</b>		3-Tape Strips, Sterile

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency


**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# ORTHOPAEDIC DRAPES & PACKS


## EXTREMITY DRAPES



### LOWER EXTREMITY DRAPE, 274cm x 203cm x 315cm

- Armboard Covers
- 6cm Conformable Circular Fenestration
- CONTROL\* Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
<b>89016</b>	Lower Extremity Drape, Sterile	Items/Case = 12



### LOWER EXTREMITY DRAPE, 223cm x 332cm

- CONTROL\* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders


Code	Description	Packaging
<b>89276</b>	Lower Extremity Drape, Sterile	Items/Case = 14



### ORTHOARTS\* LOWER EXTREMITY DRAPE, 289cm x 223cm x 320cm

- Armboard Covers
- 6cm Conformable Circular Fenestration
- Heavy-Duty Fabric
- CONTROL\* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
<b>89278</b>	ORTHOARTS* Lower Extremity Drape, Sterile	Items/Case = 7
13233	ORTHOARTS* Lower Extremity Drape, Non-Sterile	Items/Case = 84



### ORTHOARTS\* EXTRA LONG LOWER EXTREMITY DRAPE, 289cm x 223cm x 381cm

- Armboard Covers
- 6cm Conformable Circular Fenestration
- Heavy-Duty Fabric
- CONTROL\* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
<b>29219</b>	ORTHOARTS* Extra Long Lower Extremity Drape, Sterile	Items/Case = 8
29220	ORTHOARTS* Extra Long Lower Extremity Drape, Non-Sterile	Items/Case = 8

All codes in **bold** indicate a sterile product.

† Reinforcements include:


**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


# ORTHOPAEDIC DRAPES & PACKS



## UNIVERSAL EXTREMITY DRAPE, 229cm x 333cm

- CONTROL\* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration


Code	Description	Packaging
<b>89279</b>	Universal Extremity Drape, Sterile	Items/Case = 14
79279	Universal Extremity Drape, Non-Sterile	Items/Case = 150



## BILATERAL LIMB DRAPE, 193cm x 304cm

- CONTROL\* Plus Fabric Reinforcement
- Two 6cm Conformable Circular Fenestrations
- Tube Holders


Code	Description	Packaging
<b>89291</b>	Bilateral Limb Drape, Sterile	Items/Case = 14



## BILATERAL LIMB DRAPE WITH ARMBEARD COVERS, 295cm x 223cm x 366cm

- CONTROL\* Plus Fabric Reinforcement
- Two 6cm Conformable Circular Fenestrations
- Tube Holders
- Armboard Covers

Code	Description	Packaging
<b>89380</b>	Bilateral Limb Drape with Armboard Covers, Sterile	Items/Case = 12



## UPPER EXTREMITY DRAPE, 193cm x 305cm

- 3cm Circular Fenestration
- CONTROL\* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
<b>89271</b>	Upper Extremity Drape, Sterile	Items/Case = 14
79273	Upper Extremity Drape, Non-Sterile	Items/Case = 20

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency


**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# ORTHOPAEDIC DRAPES & PACKS

## HAND DRAPES




### HAND DRAPE, 289cm x 360cm

- 157cm x 147cm Hand Table Cover
- 3cm Conformable Circular Fenestration
- CONTROL\* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
------	-------------	-----------

<b>77450</b>	Hand Drape, Sterile	Items/Case = 11
--------------	---------------------	-----------------

13227	Hand Drape, Non-Sterile	Items/Case = 108
-------	-------------------------	------------------



### ORTHOARTS\* HAND DRAPE, 193CM X 381CM

- 157cm x 147cm Gusseted Hand Table Cover
- 3cm Conformable Fenestration
- Heavy-Duty Fabric
- CONTROL\* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
------	-------------	-----------

<b>89272</b>	ORTHOARTS* Hand Drape, Sterile	Items/Case = 6
--------------	--------------------------------	----------------

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric


**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.



# ORTHOPAEDIC DRAPES & PACKS


## HIP DRAPES



### ORTHOARTS\* HIP DRAPE WITH LEG POCKETS, 289cm x 223cm x 320cm

- Leg Pockets On Both Sides
- Armboard Covers
- 21cm x 17cm Conformable Teardrop Fenestration
- Heavy-Duty Fabric
- CONTROL\* Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
<b>89345</b>	ORTHOARTS* Hip Drape with Leg Pockets, Sterile	Items/Case = 6
79345	ORTHOARTS* Hip Drape with Leg Pockets, Non-Sterile	Items/Case = 6



### ORTHOARTS\* HIP DRAPE, 289cm x 223cm x 320cm

- Armboard Covers
- 21cm x 17cm Conformable Teardrop Fenestration
- Heavy-Duty Fabric
- CONTROL\* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
<b>89346</b>	ORTHOARTS* Hip Drape, Sterile	Items/Case = 7
79346	ORTHOARTS* Hip Drape, Non-Sterile	Items/Case = 7



### HIP DRAPE WITH CLEAR LEG POCKETS, 294cm x 223cm x 345cm

- Clear Leg Pockets
- Armboard Covers
- 21cm x 17cm Conformable Teardrop Fenestration
- CONTROL\* Plus Fabric Reinforcement
- Hook and Loop Tube and Line Control

Code	Description	Packaging
<b>89349</b>	Hip Drape with Clear Leg Pockets, Sterile	Items/Case = 9
79349	Hip Drape with Clear Leg Pockets, Non-Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric


**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# ORTHOPAEDIC DRAPES & PACKS

## KNEE ARTHROSCOPY DRAPES


### KNEE ARTHROSCOPY DRAPE, 228cm x 315cm



- Fluid Collection Pouch with Suction Port
- 76cm x 99cm x 8cm Triangular Pouch
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
<b>89265</b>	Knee Arthroscopy Drape, Sterile	Items/Case = 10
79265	Knee Arthroscopy Drape, Non-Sterile	Items/Case = 150


### KNEE ARTHROSCOPY DRAPE WITH ARMBUARD COVERS, 290cm x 223cm x 315cm



- Armboard Covers
- Fluid Collection Pouch with Suction Port
- 76cm x 99cm x 9cm Triangular Pouch
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
<b>89371</b>	Knee Arthroscopy Drape with Armboard Covers, Sterile	Items/Case = 10
79371	Knee Arthroscopy Drape with Armboard Covers, Non-Sterile	Items/Case = 13


### ORTHOARTS\* KNEE ARTHROSCOPY DRAPE, 289cm x 223cm x 294cm



- Armboard Covers
- Fluid Collection Pouch with Reflux Shields and Suction Port
- 76cm x 99cm x 9cm Triangular Pouch
- Heavy-Duty Fabric
- CONTROL\* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
<b>89267</b>	ORTHOARTS* Knee Arthroscopy Drape, Sterile	Items/Case = 7
79267	ORTHOARTS* Knee Arthroscopy Drape, Non-Sterile	Items/Case = 7

### HALYARD BASICS\* KNEE ARTHROSCOPY PACK



Code	Description	Packaging
<b>88886</b>	HALYARD BASICS* Knee Arthroscopy Pack	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
<b>88660</b>	78660	1-HALYARD BASICS* Film Split Drape, 152cm x 178cm, with 10cm x 53cm Split
<b>88661</b>	78661	1-HALYARD BASICS* Knee Arthroscopy Drape with Fluid Collection Pouch, 229cm x 315cm
<b>89564</b>	79764	1-HALYARD BASICS* Back Table Cover, Zone-Reinforced, 152cm x 229cm
88665	78656	1-HALYARD BASICS* Mayo Stand Cover, Reinforced, 58cm x 137cm

All codes in **bold** indicate a sterile product.

† Reinforcements include:


**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


# ORTHOPAEDIC DRAPES & PACKS



## ORTHOARTS\* EXTRA LONG KNEE ARTHROSCOPY DRAPE, 294cm x 233cm x 447cm

- Armboard Covers
- Fluid Collection Pouch with Reflux Shields and Suction Port
- 95cm x 112cm Pentagonal Pouch
- Heavy-Duty Fabric
- CONTROL\* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
<b>36930</b>	ORTHOARTS* Extra Long Knee Arthroscopy Drape, Sterile	Items/Case = 10
42429	ORTHOARTS* Extra Long Knee Arthroscopy Drape, Non-Sterile	Items/Case = 10



## HALYARD BASICS\* KNEE ARTHROSCOPY DRAPE, 229cm x 315cm

- Fluid Collection Pouch with Suction Port
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
<b>88661</b>	Knee Arthroscopy Drape, Sterile	Items/Case = 20
78661	Knee Arthroscopy Drape, Non-Sterile	Items/Case = 30

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency


**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# ORTHOPAEDIC DRAPES & PACKS


## SHOULDER DRAPES



### ORTHOARTS\* BEACH CHAIR SHOULDER ARTHROSCOPY DRAPE, 411cm x 261cm

- Attached Fluid Collection Pouch
- 95cm x 97cm Pentagonal Pouch
- Hook and Loop Tube and Line Control
- 6cm Conformable Tear Drop Fenestration
- Heavy Duty Fabric


Code	Description	Packaging
<b>36932</b>	ORTHOARTS* Beach Chair Shoulder Arthroscopy Drape, Sterile	Items/Case = 10
42426	ORTHOARTS* Beach Chair Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 10



### BEACH CHAIR SHOULDER DRAPE WITH ATTACHABLE FLUID COLLECTION POUCH, 437cm x 254cm

- Attachable Fluid Collection Pouch with Suction Port, 89cm x 109cm
- CONTROL\* Plus Fabric Reinforcement
- Hook and Loop Tube and Line Control
- Pouch Attachment Zones
- 90cm x 110cm Trapezoidal Pouch

Code	Description	Packaging
<b>89066</b>	Beach Chair Shoulder Drape with Attachable Fluid Collection Pouch, Sterile	Items/Case = 6
79049	Beach Chair Shoulder Drape without Attachable Fluid Collection Pouch, Non-Sterile	Items/Case = 11
79061	Attachable Fluid Collection Pouch, Non-Sterile	Items/Case = 44



### BEACH CHAIR SHOULDER ARTHROSCOPY DRAPE II, 406cm x 259cm

- Attached Fluid Collection Pouch
- 103m x 63cm x 2cm Pentagonal Pouch
- Hook and Loop Tube and Line Control
- 6cm Conformable Tear Drop Fenestration
- Instrument Pad Fabric Reinforcement

Code	Description	Packaging
<b>89070</b>	Beach Chair Shoulder Arthroscopy Drape, Sterile	Items/Case = 10
79070	Beach Chair Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 12

All codes in **bold** indicate a sterile product.

† Reinforcements include:


**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

# ORTHOPAEDIC DRAPES & PACKS




## 2-PIECE SHOULDER ARTHROSCOPY DRAPE PACK

- Body Split Drape, 223cm x 294cm
- 10cm x 68cm Split
- CONTROL\* Plus Fabric Reinforcement
- Tube Holders
- Shoulder Arthroscopy Drape, 254cm x 157cm
- Fluid Collection Pouch with Suction Ports
- 55cm x 112cm Trapezoidal Pouch
- Tube Holders




Code	Description	Packaging
<b>89490</b>	2-Piece Shoulder Arthroscopy Drape Pack, Sterile	Items/Case = 12
79490	2-Piece Shoulder Arthroscopy Drape Pack, Non-Sterile	Items/Case = 12



## ORTHOARTS\* SHOULDER ARTHROSCOPY DRAPE, 254cm x 152cm

- Fluid Collection Pouch
- 79cm x 149cm x 48cm Trapezoidal Pouch
- Heavy-Duty Fabric


Code	Description	Packaging
<b>89492</b>	ORTHOARTS* Shoulder Arthroscopy Drape, Sterile	Items/Case = 12
79492	ORTHOARTS* Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 15



## ORTHOARTS\* SHOULDER SPLIT DRAPE, 284cm x 152cm

- 15cm x 20cm Dual Splits
- Heavy-Duty Fabric
- CONTROL\* Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
<b>89274</b>	ORTHOARTS* Shoulder Split Drape, Sterile	Items/Case = 12



## HALYARD BASICS\* SHOULDER SPLIT DRAPE, 285cm x 152cm

- 13cm x 34cm Split
- Fabric Reinforcement

Code	Description	Packaging
<b>89559</b>	Shoulder Split Drape, Sterile	Items/Case = 30
79759	Shoulder Split Drape, Non-Sterile	Items/Case = 50



## HALYARD BASICS\* SHOULDER ARTHROSCOPY DRAPE, 254cm x 157cm

- Fluid Collection Pouch with Suction Port
- Tube Holders

Code	Description	Packaging
<b>89558</b>	Shoulder Arthroscopy Drape, Sterile	Items/Case = 30
79758	Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 50

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency


**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# ORTHOPAEDIC DRAPES & PACKS


## U & BAR DRAPES



### U-DRAPE, 193cm x 305cm

- 10cm x 101cm Split
- CONTROL\* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape


Code	Description	Packaging
<b>89301</b>	U-Drape, Sterile	Items/Case = 17
13231	U-Drape, Non-Sterile	Items/Case = 190



### LARGE U-DRAPE, 223cm x 317cm

- 15cm x 137cm Split
- CONTROL\* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape

Code	Description	Packaging
<b>89311</b>	Large Ortho U-Drape, Sterile	Items/Case = 11
79363	Large Ortho U-Drape, Non-Sterile	Items/Case = 15



### ORTHOARTS\* U-DRAPE, 223CM x 274CM

- 10cm x 101cm Split
- Heavy-Duty Fabric
- CONTROL\* Plus Fabric Reinforcement

Code	Description	Packaging
<b>89347</b>	ORTHOARTS* U-Drape, Sterile,	Items/Case = 10
13294	ORTHOARTS* U-Drape, Non-Sterile	Items/Case = 132

All codes in **bold** indicate a sterile product.

† Reinforcements include:


**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


# ORTHOPAEDIC DRAPES & PACKS



## ORTHOARTS\* BODY SPLIT DRAPE, 223CM x 294CM

- 10cm x 68cm Split
- Heavy-Duty Fabric
- CONTROL\* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape


Code	Description	Packaging
<b>89275</b>	ORTHOARTS* Body Split Drape, Sterile	Items/Case = 10
79275	ORTHOARTS* Body Split Drape, Non-Sterile	Items/Case = 12



## SPLIT DRAPE, 152cm x 178cm

- Impervious
- Polyethylene Material
- 10cm x 53cm Split


Code	Description	Packaging
<b>89331</b>	Split Drape, Sterile	Items/Case = 42
79330	Split Drape, Non-Sterile	Items/Case = 115



## LARGE SPLIT DRAPE, 152cm x 205cm

- Impervious
- Polyethylene Material
- 11cm x 78cm Split

Code	Description	Packaging
<b>89333</b>	Large Split Drape, Sterile	Items/Case = 34
79333	Large Split Drape, Non-Sterile	Items/Case = 100



## HALYARD BASICS\* BODY SPLIT DRAPE, 224cm x 295cm

- 10cm x 67cm Split
- Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
<b>89557</b>	Body Split Drape, Sterile	Items/Case = 25
79757	Body Split Drape, Non-Sterile	Items/Case = 30

All codes in **bold** indicate a sterile product.

† Reinforcements include:


**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# ORTHO PAEDIC DRAPES & PACKS




## ORTHO SPLIT "T" DRAPE, 297cm x 210cm x 342cm

- Armboard Covers
- 10cm x 53cm Split
- CONTROL\* Fabric Reinforcement
- Tube Holders
- Positioning Tape

Code	Description	Packaging
------	-------------	-----------

<b>36931</b>	Ortho Split "T" Drape, Sterile	Items/Case = 10
--------------	--------------------------------	-----------------




## ORTHO BAR DRAPE, 254cm x 152cm

- CONTROL\* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
------	-------------	-----------

<b>89381</b>	Ortho Bar Drape, Sterile	Items/Case = 18
--------------	--------------------------	-----------------

79360	Ortho Bar Drape, Non-Sterile	Items/Case = 28
-------	------------------------------	-----------------



## VERTICAL ISOLATION DRAPE, 228cm x 320cm

- Clear Film
- Adhesive Incise Area 40cm x 76cm
- Attachable Fluid Collection Pouch with Suction Port
- 42cm x 84cm x 17cm Trapezoidal Pouch
- Tube Holders
- Central Release Tape Strip For Horizontal Rod Attachment

Code	Description	Packaging
------	-------------	-----------

<b>77723</b>	Vertical Isolation Drape, Sterile	Items/Case = 14
--------------	-----------------------------------	-----------------

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency


**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.



# ORTHOPAEDIC DRAPES & PACKS



## ISOLATION SCREEN, 170cm x 290cm

- Clear Plastic Drape 170cm x 290cm
- Adhesive Backing Strips

Code	Description	Packaging
<b>89093</b>	Isolation Screen, Sterile	Items/Case = 24
79093	Isolation Screen, Non-Sterile	Items/Case = 30



## STOCKINETTES

Code	Description	Packaging
<b>89811</b>	Stockinette, Sterile, Medium, 23cm x 112cm	Items/Case = 36
79505	Stockinette, Non-Sterile, Medium, 24cm x 117cm	Items/Case = 75
<b>89821</b>	Stockinette, Sterile, Large, 30cm x 122cm	Items/Case = 25
79515	Stockinette, Non-Sterile, Large, 29cm x 117cm	Items/Case = 50

All codes in **bold** indicate a sterile product.

† Reinforcements include:


**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

# UROLOGY DRAPES & PACKS




## T.U.R. PACK

Code	Description	Packaging
<b>77885</b>	T.U.R. Pack	Items/Case = 12
	Non-Sterile Component Code	Contents
		1-Back Table Cover, Zone-Reinforced, 152cm x 191cm
		1-Absorbent Towel, 38cm x 50cm
76885	1-T.U.R. Drape, 289cm x 175cm x 170cm:	
		<ul style="list-style-type: none"> <li>• Attached Leggings</li> <li>• Fluid Collection Pouch with Screen, Suction Port and Tie</li> <li>• 54cm x 75cm Triangular Pouch</li> <li>• Finger Cot</li> <li>• 7cm Fenestration</li> <li>• 25cm x 20cm Fenestration</li> </ul>



## CYSTOSCOPY PACK I

Code	Description	Packaging
<b>88611</b>	Cystoscopy Pack I	Items/Case = 12
	<b>Sterile Component Code</b>	Non-Sterile Component Code
		Contents
		1-Back Table Cover, Zone-Reinforced, 112cm x 198cm
		1-Absorbent Towel, 38cm x 56cm
<b>89408</b>	79408	2-Leggings, with 15cm Cuff, 78cm x 122cm
		1-Cystoscopy "T" Drape with CONTROL* Plus Reinforcement, 152cm x 76cm x 213cm with 7cm x 15cm Perineal Fenestration and 23cm x 23cm Screen



## PERCUTANEOUS NEPHROLITHOTOMY DRAPE (PCNL), 299cm x 182cm

- CONTROL\* Fabric Reinforcement
- Tube Holders
- Fluid Collection Pouch
- 67cm x 61cm x 7cm Trapezoidal Pouch
- 18cm x 28cm Fenestration

Code	Description	Packaging
<b>42427</b>	Percutaneous Nephrolithotomy Drape, Sterile	Items/Case = 12

All codes in **bold** indicate a sterile product.

† Reinforcements include:

**CONTROL\* Plus** – Maximum Absorbency

**CONTROL\*** – Absorbent Reinforcement Fabric

**Instrument Pad** – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

# STANDARD PACKS

## ENT PACK

Code	Description	Packaging
<b>ENTEPHYH-1A</b>	ENT Pack Contents	Items/Case = 1
	1-Back Table Cover, Unfolded, 152cm x 229cm 1-Bag, Ziplock, 16cm x 25cm 2-Blade, Surgical, Carbon Steel, 15 1-Bowl, 1000ml 1-Bowl, 250ml, Clear 1-Cover, Light Handle, Flexible, Green 2-Dish, Kidney, 700ml, Clear 1-Drape, Head Turban, 127cm x 91cm 1-Drape, EENT, Split, 193cm x 315cm 3-Gallipot, 150ml 5-Gauze, 8 Ply, 10cm x 10cm, XRD 1-Handle, Yankaeur, Bulb Tip, Non-Vented 1-Mayo Stand Cover, 58cm x 137cm 2-Needle, 18G x 38mm 2-Needle, 23G x 25mm 1-Needle Counter, 20-Count, Foam Block, Magnet, Blade Remover, Yellow 2-Syringe, 10ml, L/L 2-Syringe, 20ml, L/L 1-Tray, Sharps Transfer, Yellow 1-Tubing, Suction, Soft, 4m	

## UNIVERSAL PACK

Code	Description	Packaging
<b>GENUPHYH-1A</b>	Universal Pack Contents	Items/Case = 1
	1-Back Table Cover, Unfolded, 152cm x 229cm 1-Bag, Ziplock, 16cm x 25cm 1-Blade, Surgical, Carbon Steel, 15 1-Bowl, 1000ml 1-Cover, Light Handle, Flexible, Green 2-Dish, Kidney, 700ml, Clear 1-Drape, Bottom, 193cm x 195cm 2-Drape, Side, 195cm x 105cm 1-Drape, Top, 273 x 127cm 2-Gallipot, 150ml 5-Gauze, 8 Ply, 10cm x 10cm, XRD 5-Label, 125mm x 5cm 1-Marker, Skin, Standard Tip 1-Mayo Stand Cover, 58cm x 137cm 2-Needle, 23G x 25mm 1-Needle Counter, 20-Count, Foam Block, Magnet, Blade Remover, Yellow 1-Prep Sponge, White, 8.4cm x 5cm x 2.4cm 1-Ruler, White, 15cm 5-Sponge, Abdominal, XRD, 30cm x 30cm 1-Syringe, 20ml, L/L 1-Tray, Sharps Transfer, Yellow	

# STANDARD PACKS

## CYSTO-GYNAE PACK

Code	Description	Packaging
<b>OBGGPHYH-1A</b>	Cysto-Gynae Pack Contents	Items/Case = 1
	1-Back Table Cover, Unfolded, 152cm x 229cm	
	1-Bowl, 1000ml	
	1-Dish, Kidney, 700ml, Yellow	
	1-Dish, Kidney, 700ml, Clear	
	2-Drape, Leggings, 78cm x 122cm	
	1-Drape, Lithotomy T	
	2-Gallipot, 150ml	
	5-Gauze, 8 Ply, 10cm x 10cm, XRD	
	1-Pad, Peri Winged 10cm x 27cm	

## CAESAREAN PACK

Code	Description	Packaging
<b>OBGCPHYH-1A</b>	Caesarean Pack Contents	Items/Case = 1
	1-Back Table Cover, Unfolded, 152cm x 229cm	
	1-Bag, Suture, with tape, 20cm x 25cm	
	1-Bag, Zip Lock, 16cm x 25cm	
	1-Blade, Surgical, Carbon Steel, 23	
	1-Blanket, Baby, 75cm x 100cm	
	1-Bowl, 1000ml	
	1-Catheter, Y-Suction, 10FG x 45cm	
	1-Cover, Light Handle, Flexible, Green	
	2-Dish, Kidney, 700ml, Clear	
	1-Drape, C-Section, Fenestrated Incise	
	2-Gallipot, 150ml	
	5-Gauze, 8 Ply, 10cm x 10cm, XRD	
	1-Handle, Yankauer, Bulb Tip, Non-Vented	
	1-Pad, Peri Winged 10cm x 27cm	
	1-Mayo Stand Cover, 58cm x 137cm	
	1-Needle Counter, 20-Count, Foam Block, Magnet, Blade Remover, Yellow	
	1-Prep Sponge, White, 8.4cm x 5cm x 2.4cm	
	5-Sponge, Abdominal, XRD, 30cm x 30cm	
	1-Syringe, 10ml, L/S	
	1-Tray, Sharps Transfer, Yellow	
	1-Tubing, Suction, 3.6m	

# STANDARD PACKS

## LOWER EXTREMITY PACK

<b>Code</b>	<b>Description</b>	<b>Packaging</b>
<b>ORTEPHYH-1A</b>	Lower Extremity Pack Contents	Items/Case = 1
	1-Back Table Cover, Unfolded, 152cm x 229cm	
	1-Bag, Ziplock, 16cm x 25cm	
	1-Bandage, 15cm x 2.7m, Medium Weight	
	1-Bowl, 1000ml	
	1-Cover, Light Handle, Flexible, Green	
	2-Dish, Kidney, 700ml, Clear	
	1-Drape, Extremity, Extra Long, 289cm x 321cm	
	1-Drape, Split, 152cm x 178cm	
	2-Gallipot, 150ml	
	5-Gauze, 8 Ply, 10cm x 10cm, XRD	
	5-Label, 125mm x 500mm, Blank	
	1-Marker, Skin, Standard Tip	
	1-Needle, 23G x 38mm	
	1-Needle Counter, 20-Count, Foam Block, Magnet, Blade Remover, Yellow	
	1-Prep Sponge, White, 8.4cm x 5cm x 2.4cm	
	1-Ruler, White, 15cm	
	5-Sponge, Abdominal, XRD, 30cm x 30cm	
	1-Stockinette, Impervious, Large, 30cm x 122cm	
	1-Syringe, 20ml, L/L	
	1-Tray, Sharps Transfer, Yellow	

# CUSTOM PROCEDURE PACKS (CPP)

## THE FUTURE OF CUSTOM PROCEDURE PACKS FROM THE NAME YOU TRUST

**For over 30 years HALYARD has been providing high quality surgical drapes and gowns, designed to deliver protection, comfort and optimal patient outcomes.**



Leveraging Halyard's global expertise in supply chain, component sourcing, utilisation reviews and cost efficiencies, we are now extending our offering to provide ANZ clinicians with HALYARD branded Custom Procedure Packs (CPPs).

**Based on years of working together with clinicians globally we know that some of the challenges you face in Operating Theatres include:**

- Operational inefficiencies
- Lack of visibility to costs and inventories
- Suboptimal utilisation of Operating Theatre supplies
- Procedural supply waste
- Need for improved throughput
- Staff satisfaction and retention

That's where HALYARD comes in. When you make the decision to switch to HALYARD CPPs, our clinical experts work with you to ensure EVERY component is used EVERY time. It's part of our PACT pledge and the HALYARD CPP Solution.

# INDEX

Code	Page	Code	Page	Code	Page	Code	Page	Code	Page	Code	Page
13227	40	42427	50	47619	15	79551	33	47622	16	89431	10
13231	46	42429	43	79117	31	79602	12	47630	16	89438	10
13233	38	42620	4	47624	16	79603	12	47628	16	89452	8
13294	46	42621	4	47631	16	79607	12	89203	33	89453	8
13296	7	43623	4	79147	16	79709	7	89204	33	89454	8
14007	12	43940	4	79204	33	79720	13	89206	21	89471	6
29219	38	43941	4	79206	21	79731	13	89207	21	89490	45
29220	38	44661	3	79207	21	79749	7	89208	21	89492	45
29240	19	44662	3	79208	21	79757	47	89209	21	89528	29
29241	19	44664	3	79209	21	79758	45	89219	31	89529	29
29258	12	44666	3	79224	31	79759	45	89224	31	89538	19
29259	12	44667	3	79233	19	79828	16	89229	19	89543	11
29269	12	44668	3	79240	20	79835	18	89233	19	89557	47
29283	20	44669	3	79255	22	79945	13	89236	17	89558	45
29284	20	44680	3	79265	42	88088	30	89241	20	89559	45
31445	18	44681	3	79267	42	88161	14	89255	22	89583	6
34901	26	44682	3	79273	39	88171	14	89265	42	89601	12
34902	26	44683	3	79275	47	88178	14	89267	42	89602	12
34903	26	44684	3	79279	39	88222	30	89271	39	89603	12
34904	26	44685	3	79282	20	88223	9	89272	40	89621	12
34905	26	44686	3	79304	35	88334	25	89274	45	89632	12
34906	26	44687	3	79330	47	88411	36	89275	47	89636	18
36930	43	44967	29	79333	47	88452	36	89276	38	89680	18
36931	48	44968	29	79335	17	88471	37	89278	38	89701	13
36932	44	44972	21	79343	16	88521	32	89279	39	89709	7
37498	12	44976	6	79345	41	88551	32	89282	20	89711	7
41725	4	44977	29	79346	41	88611	50	89291	39	89731	13
41727	4	44978	29	79349	41	88661	43	89301	46	89749	7
41739	4	70392	5	79351	10	88665	12	89304	35	89811	49
41740	4	72978	5	79356	35	88666	12	89311	46	89821	49
41741	4	72988	5	79360	48	88711	9	89331	47	89828	16
41870	4	76323	15	79363	46	88716	9	89333	47	89835	18
41871	4	76356	8	79368	6	88719	10	89343	16	89945	13
41872	4	76358	12	79371	42	88759	23	89345	41	90003	5
41873	4	77450	40	79408	33	88769	23	89346	41	90018	5
41874	4	77723	48	79411	34	88773	25	89347	46	90019	5
42217	12	77885	50	79414	34	88774	24	89349	41	90048	5
42218	12	78656	12	79415	34	88875	25	89352	10	90049	5
42224	12	78661	42	79430	10	88878	24	89353	15	90050	5
42225	12	78664	12	79452	8	88886	42	89356	35	90138	5
42231	12	79049	44	79453	8	89016	38	89368	6	ENTEPHYH-1A	51
42233	12	79061	44	79482	12	89023	15	89371	42	GENUPHYH-1A	51
42306	12	79070	44	79490	45	89053	13	89380	39	OBGPHYH-1A	52
42307	12	79075	28	79492	45	89062	28	89381	48	OBGCPHYH-1A	52
42309	12	79093	49	79505	49	89066	44	89408	33	ORTEPHYH-1A	53
42310	12	79095	27	79515	49	89070	44	89413	34	CPP	54
42366	12	79096	27	79528	29	89093	49	89414	34		
42367	12	79097	27	79529	29	47614	15	89415	34		
42426	44	47615	15	79538	19	89111	15	89421	11		



Infection Prevention Product Solutions  
Halyard\* Education Foundation  
Tools & Best Practice  
Commitment to Excellence and Sustainability  
Experienced Customer Support  
Knowledgeable Sales Force  
Innovation & Technology



**O&M Halyard Australia Pty Ltd.**  
**Customer Care: 1800 664 227**  
**or email [CustomerService\\_ANZ@hyh.com](mailto:CustomerService_ANZ@hyh.com)**

[www.halyardhealth.com.au](http://www.halyardhealth.com.au)



Halyard Health makes no representation that its products described in this catalogue will provide complete protection from exposure to contaminants, fluids or micro-organisms or will eliminate the risk of contracting infectious diseases.

\*Registered Trademark or Trademark of O&M Halyard or its affiliates. ©2022 O&M Halyard International. All Rights Reserved. 21:12-WFHY932v1