




Products supported by: 




OPTIVISC SINGLE

90MG/3ML

BIOTECH

Optivisc Single: One injection, lasting relief.





OPTIVISC 2ML

BIOTECH

Optivisc 2ml: A step towards pain-free living.





INSUPEN PEN NEEDLES G32X4MM PIC

INSUPEN PEN NEEDLES G32X4MM, 0.23mm/100pcs




SUPEN PEN NEEDLES G31X6MM

PIC

INSUPEN PEN NEEDLES G31X6MM, 0.25mm /100pcs





INSUPEN PEN NEEDLES G31X8MM

PIC

INSUPEN PEN NEEDLES G31X8MM, 0.25mm /100pcs

Pikdare Insupen Pen Needles: Your partner in diabetes management






FEMALE CATHETER

PENNINE

Pennine Female Catheters: Gentle and reliable.

FEMALE CATHETER, Intermittent 23cm
10ch PVC | FC-1410 & 12ch PVC | FC-1412






CPAP AUTO AEONMED AS100A

VIOBREATH - AEONMED

Your path to better sleep starts here.





Choose the perfect mask for your needs:

VIOBREATH - YUWELL



PILLOW MASK - YP-01



MASK YUWELL W/S - YF-01



MASK NASAL YUWELL - YN-02



Proudly Presents
The Newest and Most
Innovative Member
of the Nuvo Family
of Oxygen Concentrators



.....Mark 5
Nuvo
Lite

- “Lite” - a mere 30 lbs. (13.6kg)
- Quiet Operation
- Compact
- Lockable Flow Control Valve
- Liter Flow -- 0.125 - 5 LPM



Metal DISS Fitting



Lockable Flow Control Valve



Circuit Breaker Reset

Specifications*

Nuvo Lite Oxygen Concentrator

115 Volt - 60 Hz

520

230 Volt - 50 Hz

920

Nuvo Lite Oxygen Concentrator with Oxygen Monitoring System

525

925

Average Power Consumption
Oxygen Concentration (at 5 LPM)
Liter Flow
Sound Level
Outlet Pressure
Weight
Dimensions

310 Watts
90% (+6.5% / -3%)
0.125 to 5 LPM
40 dBA typical
7.0 PSI
30 lbs
14"W x 9"D x 23"H

290 Watts
90% (+6.5% / -3%)
0.125 to 5 LPM
40 dBA typical
7.0 PSI
13.6 kg
35 cm x 22 cm x 58 cm

Accessories and Re-order Items

Compressor Intake Filter
External Cabinet Filter
Humidifier
Cannula with 7' tubing
Nuvo Lite User's Guide

8400-1180
8400-1025
9012-8774
9012-8780
2010-8400

8400-1180
8400-1025
9012-8774
9012-8780
2010-8400CE

* NOM - Specifications are subject to change without notice



NidekMedical
PRODUCTS



CE
2862

..... 8 & 10 LPM
Nuvo

Specifications

.....Mark 5
Nuvo8



Model 585	Model 985
115 V / 60 Hz	230 V, 50/60 Hz
500 Watts	490 Watts



90% (+5.5% / -3%)
at 8 and 10 LPM



2 to 8 LPM
2 to 10 LPM



< 58 dBA




15 PSI (Nuvo 8)
20 PSI (Nuvo 10)



54-58 lbs
24.5-26 kg



Model 595	Model 1005
115 V / 60 Hz	230 V / 50 Hz
700 Watts	600 Watts



15.5" W x 15.6" D x 27.8" H
39 cm x 39 cm x 71 cm

.....Mark 5
Nuvo10



www.cimbiomedical.com.cy

Products supported by: **ΓΕΟΥ**